

730 Digitizer Family

16/8 Ch. 14-bit 500 MS/s Digitizer



APPLICATIONS

- Nuclear and Particle Physics
- Dark Matter and Astroparticle Physics
- Fast Neutron spectroscopy
- Homeland Security

FORM FACTOR



VME64 VME64X NIM DESKTOP

FUNCTIONS

ICH WV TS TDC QDC CFD PHA PSD ZS DAW

The CAEN top level digitizer

APPLICATION NOTES

- AN3250
- AN5157
- AN5995
- GD2827
- AN3251
- AN5830
- GD2783

Overview

The 730 is a family of CAEN Waveform Digitizers able to perform basic waveform recording and run online advanced algorithms (DPP) for digital pulse processing: charge integration and pulse shape discrimination with constant fraction timing, pulse height analysis, zero-length encoding, and dynamic acquisition window. The wide selection of DPP algorithms supported by this family makes the 730 a “must-have” for any kind of nuclear physics application.













Data is read by a Flash ADC, 14-bit resolution and 500 MS/s sampling rate, which is well suited for mid-fast signals as the ones coming from liquid or inorganic scintillators coupled to PMTs or Silicon Photomultipliers, but also for high precision detectors as Silicon or HPGe coupled with charged sensitive preamplifiers. In this case the sampling rate can be reduced using the decimation feature. The acquisition can be channel independent and it is possible to make coincidence/anti-coincidence logic among different channels and external veto/gating. Multiple boards can be synchronized to build up complex systems. In case of DPP mode, data can be saved in time-stamped list mode to support higher input rates and improving the throughput performances. Piled-up events can be rejected or saved for offline analysis. The acquisition in DPP mode is fully controlled by the CoMPASS and MC²Analyzer software, which manage the algorithm parameters, build the plots and saves the relevant energy, time, and PSD spectra. In case of waveform recording mode, the user can take advantage

of the CAENScope and WaveDump software to access and save the waveforms. Libraries and demo software in C and LabView are available for integration and customization of specific acquisition systems.

730 family comes in three form factors: VME (16 or 8 input channels), NIM (8 input channels) and Desktop (8 input channels). The communication to and from the board is provided through the following interfaces: USB (Desktop and NIM form factors), VMEbus (VME form factor), and Optical Link (all form factors).

Features

- 14-bit @ 500 MS/s
- Analog inputs on MCX coax. connectors
- VME64/VME64X (16/8 ch.), NIM (8 ch.) and Desktop (8 ch.) modules
- 0.5 and 2 Vpp input dynamic range with programmable DC offset adjustment
- Algorithms for Digital Pulse Processing (DPP)
- VME, USB and Optical Link communication interfaces
- Multi-board synchronization features
- Daisy chain capability
- Demo software tools, Control Software for waveform recording and DPP firmware, C and LabVIEW libraries

Firmware	Software	ICH	WV	TS	TDC	QDC	CFD	PHA	PSD	ZS	DAW
	 		•	•							
		•	•	•	•	•	•		•		
	 	•	•	•	•			•			
			•	•						•	
		•	•	•							•

Principle of Operation


CAEN Waveform Digitizers are devices able to continuously acquire analog input signals, which are sampled by fast ADCs and stored into digital memories where they are available for readout through different communication interfaces (USB, VMEBus, Optical Link). Depending on the FPGA firmware, the digitizer can work like an oscilloscope to acquire raw waveforms, or perform online processing to calculate parameters such as pulse height, charge, time stamp, pulse shape discrimination. In the latter case, the output data is a time-stamped list of parameters. Data reductions and zero suppression algorithms are also available. Digitization in CAEN digitizers is based on two main technics: **Flash ADC** and **Switched Capacitor Arrays**. Flash ADC are the fastest A/D converters, where the sampling and the analog-to-digital conversion are made practically at the same time. Flash ADC are so not affected by dead-time due to conversion. In the Switched Capacitor Arrays, the sampling and the A/D conversion take place at different times, thus introducing a dead-time. Despite of the dead-time, the Switched Capacitor Array Digitizers are able to sample the input pulse at very high frequency, up to 5 GS/s, with high channel density, while the highest Flash ADC frequency is 4 GS/s with a quite low number of channels. If compared to a commercial digital oscilloscope, the waveform digitizer presents a list of differences which make it an advanced instrument for many applications:

- waveform recording with no dead-time due to conversion (Flash ADC digitizers)
- Independent channel self-triggering and event acquisition
- On-line digital pulse processing (DPP) algorithms
- Data reduction
- Multi-board synchronization for system scalability
- Communication interfaces with high bandwidth readout

The benefits of the digital approach are great stability and reproducibility, ability to reprogram and adjust the algorithms to the application, ability to preserve the information of the signal along the entire acquisition chain, flexibility, better correction of baseline fluctuation, pile-up, ballistic deficit, etc.. All in one board.

CAEN Digitizer block diagram:

- A motherboard which contains one FPGA for the readout interfaces and the services, and defines the form factor.
- One or more daughterboards which define the type of digitizer and contain the signal conditioning input stage, the ADCs, the FPGA for the data processing and the memories.



CAEN has made a big effort in dev

The user can install a DPP algorithm that goes beyond the simple

Firmware Upgrade

Firmware



Waveform Recording Firmware

Time windowed waveform recording
Dead-timeless acquisition
Multi-buffer memories
Multi-board synchronization

Software



CAEN Digitizer readout application



Digitizer Software for Signal Inspection and waveform recording

Features

ICH

WV

TS

TDC

QDC

CFD

PHA

PSD

ZS

D

Feature Acronyms Legend

ICH

(Independent channels): the trigger and acquire data independently can be correlated using a user (anti) coincidence logic criterion

TDC

(Time to Digital Converter): stamp information from the fast digital discriminator filter resolution (sub-ns resolution)

PHA

(Pulse Height Analysis): pul exponential signals or comin sensitive preamplifier


DAW

(Dynamic Acquisition Window): acquisition of a dynamically a number of ADC samples

Developing algorithms for the Digital Pulse Processing (DPP).

firmware on the FPGA of the digitizer (firmware upgrade), run it on-line and implement new acquisition and waveform recording.

725 / 730 Digitizer Family



Pulse Height Analysis

Digitizer running DPP-PHA firmware becomes a digital replacement of Shaping Amplifier and Peak Sensing ADC (Multi-Channel Analyzer) by means of digital trapezoidal filters.



Pulse Shape Discrimination

Digitizer running DPP-PSD firmware becomes a digital replacement of dual gate QDC, discriminator and gate generator.



Zero Length Encoding

New

Digitizer running DPP-ZLE firmware becomes able to transfer digitized waveforms in compressed mode thanks to an enhanced Zero Suppression algorithm.



Dynamic Acquisition Window

New

Digitizer running DPP-DAW firmware allows the user to acquire data from each channel independently and to adapt the acquisition window to the real pulse duration (over-threshold).



Multiparametric DAQ Software for Physics Applications



Digital MCA Data Acquisition and Analysis Software



Multiparametric DAQ Software for Physics Applications



Software interface for DPP-ZLEplus management



Software interface for DPP-DAW management

- ICH
- WV
- TS
- TDC
- QDC
- CFD
- PHA
- PSD
- ZS
- DAW

- ICH
- WV
- TS
- TDC
- QDC
- CFD
- PHA
- PSD
- ZS
- DAW

- ICH
- WV
- TS
- TDC
- QDC
- CFD
- PHA
- PSD
- ZS
- DAW

- ICH
- WV
- TS
- TDC
- QDC
- CFD
- PHA
- PSD
- ZS
- DAW

da

channels can independently and user-defined

fine time interpolation or with high resolution

height of any signal from charge

low): adjusted



(Waveforms): acquisition of a programmable number of ADC samples (raw waveforms)



(Charge to Digital Converter): gated integrator (charge)



(Pulse Shape Discrimination): particle discrimination via double gated integration



(Timestamp): trigger coarse time stamp with low resolution (10-20 ns)



(Constant Fraction Discriminator): digital Constant Fraction Discriminator with programmable delay and fraction. Works in combination with TDC



(Zero Suppression): acquisition of a programmable number of ADC samples with zero suppression

Technical Specifications

GENERAL

Form Factor

1-unit wide, 6U VME64/VME64X

1-unit wide NIM

154x50x164 mm³ (WxHxD) Desktop

ANALOG INPUT

Channels

16/8 channels, single ended (VME)

8 channels, single ended (NIM, Desktop)

Impedance

50 Ω

Connector

MCX

Full Scale Range (FSR)

2 Vpp (0.5 Vpp software selectable)

Bandwidth

250 MHz

Offset

Programmable DAC for DC offset adjustment. Range: ± 1 V @ 2 Vpp, ± 0.25 V @ 0.5 Vpp

DIGITAL CONVERSION

Resolution

14 bits

Sampling rate

500 MS/s simultaneously on each channel

ADC CLOCK GENERATION

Clock source: internal/external

On-Board PLL provides ADC sampling clock generation from an internal (50 MHz loc. oscillator) or external reference (50 MHz or 62.5 MHz; other options on request) on front panel CLK-IN connector

MEMORY

640 kS/ch or 5.12 MS/ch Multi-Event Buffer divisible into $1 \div 1024$ buffers with independent read and write access.

Programmable event size and pre-/post-trigger

TRIGGER

Trigger source

Self-trigger: channel over/under threshold for either Common or Individual (DPP firmware only) trigger generation

External-trigger: Common by TRG-IN or Individual by LVDS connectors (DPP firmware only)

Software-trigger: Common by software command

Trigger propagation

TRG-OUT (VME) / GPO (NIM, Desktop) digital output

Trigger Time Stamp

Waveform recording firmware/DPP-DAW/DPP-ZLE: 31-bit counter, 16 ns resolution, 17 ns range(*); 48-bit extension by firmware

DPP-PHA/PSD: 31-bit counter, 2 ns resolution, 4 s range; 47-bit extension by firmware; 10-bit and 2 ps fine time stamp by digital CFD; 64-bit extension by software

SYNCHRONIZATION

Clock propagation

Daisy chain (VME only) through CLK-IN/CLK-OUT connectors

One-to-many clock distribution from an external clock source

Clock Cable delay compensation

Acquisition Synchronization

Sync Start/Stop through digital I/O (S-IN, TRG-IN or GPI input, TRG-OUT or GPO output)

External Trigger Time Stamp reset

LVDS I/O (VME only)

16 general purpose LVDS I/Os controlled by FPGA

Busy, Data Ready, Memory full, Individual Trig-Out and other functions can be programmed

An Input Pattern from the LVDS I/Os can be associated to each trigger as an event marker

ANALOG MONITOR (VME only)

12-bit/125 MHz DAC FPGA controlled output with four operating modes: Trigger Majority / Test Pulses / Memory Occupancy / Voltage Level

COMMUNICATION INTERFACE

Optical Link

CAEN CONET proprietary protocol, up to 80 MB/s transfer rate

Daisy chainable: it is possible to connect up to 8/32 ADC modules to a single Optical Link Controller (Mod. A2818/A3818)

USB (NIM and Desktop direct, VME via V1718 bridge)

USB 2.0 compliant

Transfer rate up to 30 MB/s

VME

VME 64X compliant

Data transfer mode: BLT32, MBLT64 (70 MB/s using CAEN Bridge), CBLT32/64, 2eVME, 2eSST (up to 200 MB/s)

POWER CONSUMPTIONS

Desktop: 2.8 A @ 12 V (Typ.)

NIM: N6730

4.9 A @ +6 V

250 mA @ -6 V

VME: V1730

8.2 A @ +5 V

840 mA @ +12 V

-12 V not used

N6730B

6.5 A @ +6 V

250 mA @ -6 V

V1730B

10.2 A @ +5 V

840 mA @ +12 V

-12 V not used

V1730C

5 A @ +5 V

400 mA @ +12 V

400 mA @ -12 V

V1730D

6.5 A @ +5 V

400 mA @ +12 V

400 mA @ -12 V

Ordering Options

Code	Description	Form Factor
WDT5730SAAAA	DT5730S - 8 Ch. 14 bit 500 MS/s Digitizer: 640kS/ch, CE30, SE	Desktop
WDT5730SBXAA	DT5730SB - 8 Ch. 14 bit 500 MS/s Digitizer: 5.12MS/ch, CE30, SE	Desktop
WN6730SAAAA	N6730S - 8 Ch. 14 bit 500 MS/s Digitizer: 640kS/ch, CE30, SE	NIM
WN6730SBXAA	N6730SB - 8 Ch. 14 bit 500 MS/s Digitizer: 5.12MS/ch, CE30, SE	NIM
WV1730SAAAA	V1730S - 16 Ch. 14 bit 500 MS/s Digitizer: 640kS/ch, CE30, SE	6U-VME64
WV1730SBXAA	V1730SB - 16 Ch. 14 bit 500 MS/s Digitizer: 5.12MS/ch, CE30, SE	6U-VME64
WV1730SCXAA	V1730SC - 8 Ch. 14 bit 500 MS/s Digitizer: 640kS/ch, CE30, SE	6U-VME64
WV1730SDXAA	V1730SD - 8 Ch. 14 bit 500 MS/s Digitizer: 5.12MS/ch, CE30, SE	6U-VME64
WVX1730SAXAA	VX1730S - 16 Ch. 14 bit 500 MS/s Digitizer: 640kS/ch, CE30, SE	6U-VME64X
WVX1730SBXAA	VX1730SB - 16 Ch. 14 bit 500 MS/s Digitizer: 5.12MS/ch, CE30, SE	6U-VME64X
WVX1730SCXAA	VX1730SC - 8 Ch. 14 bit 500 MS/s Digitizer: 640kS/ch, CE30, SE	6U-VME64X
WVX1730SDXAA	VX1730SD - 8 Ch. 14 bit 500 MS/s Digitizer: 5.12MS/ch, CE30, SE	6U-VME64X
WFWDPPPTFAA30	DPP-PHA - Digital Pulse Processing for Pulse Height Analysis (16 ch x730)	6U-VME64 VME64X
WFWDPPNGAA30	DPP-PSD - Digital Pulse Processing for Pulse Shape Discrimination (16ch x730)	6U-VME64 VME64X
WFWDPPPTFAD30	DPP-PHA - Digital Pulse Processing for Pulse Height Analysis for (8ch x 730)	ALL
WFWDPPNGAD30	DPP-PSD - Digital Pulse Processing for Pulse Shape Discrimination (8ch x730)	ALL
WFWDPPDWAA30	DPP-DAW - Digital Pulse Processing with Dynamic Acquisition Windows (16ch x730)	6U-VME64 VME64X
WFWDPPDWAD30	DPP-DAW - Digital Pulse Processing with Dynamic Acquisition Windows (8ch x730)	ALL
WFWDPPZLAA30	DPP-ZLEplus - Digital Pulse Processing Zero Length Encoding for (16ch x730)	6U-VME64 VME64X
WFWDPPZLAD30	DPP-ZLEplus - Digital Pulse Processing Zero Length Encoding for (8ch x730)	ALL
WFWDPPS3001A	DPP-SUP - Super Licence for 16ch x 730 Digital Pulse Processing	6U-VME64 VME64X
WFWDPPS3001D	DPP-SUP - Super Licence for 8ch x 730 Digital Pulse Processing	ALL

Accessories

A2818

PCI CONET Controller



A3818

PCI Express CONET2 Controller



A654

MCX to LEMO Cable Adapter



A659

MCX to BNC Cable Adapter



A317

Clock Distribution Cable



A318

SE to Differential Clock Cable Adapter



AI2700

Optical Fiber Series



DT4700

Clock Generator and FAN-OUT



Cables for CONET Optical Link Networks

(*) Trigger Logic and Trigger Time Stamp counter operate at 125 MHz (i.e. 8 ns or 4 ADC clock cycles), while the counter value is read at a frequency of 62.5 MHz (i.e. 16 ns).