



Rev. 1 - 7th December 2022

SP5600C

Educational Gamma Kit



1.1 Purpose of this Guide



This QuickStart Guide contains basic information and examples that will let you use Educational Gamma kit in few steps.

1.2 Change Document Record

Date	Revision	Changes
September 2016	00	Initial release.
December 2022	01	Updated §Getting started, §Basic Measurements, §Educational Experiments and §Technical Support. Added New §PID (Product Identifier), §Hardware Description, §Software Description, §Appendix, §Instructions for Cleaning, §Device decommissioning and §Disposal.

1.3 Symbols, abbreviated terms and notation

AMC FPGA	Acquisition & Memory Controller FPGA
DPP	Digital Pulse Processing
FPGA	Field Programmable Gate Array
OS	Operating System
PSAU	Power Supply & Amplification Unit
ROC FPGA	ReadOut Controller FGPA
SiPM	Silicon Photo-Multiplier
GUI	Graphical User Interface
PSAU	Power Supply and Amplification Unit

1.4 Reference Documents

[RD1]	DT5720 User Manual
[RD2]	UM1935 – CAENDigitizer User & Reference Manual
[RD3]	GD2783 - First Installation Guide to Desktop Digitizers & MCA
[RD4]	GD7873 - Digital Pulse Processing for SiPM kit
[RD5]	DS2626 – SP5600 Power Supply and Amplification

<https://www.caen.it/support-services/documentation-area/>

1.5 Manufacturer Contacts



CAEN S.p.A.
Via Vetràia, 11 55049 Viareggio (LU) - ITALY
Tel. +39.0584.388.398 Fax +39.0584.388.959
www.caen.it | info@caen.it
© CAEN SpA – 2022

1.6 Limitation of Responsibility

If the warnings contained in this manual are not followed, CAEN will not be responsible for damage caused by improper use of the device. The manufacturer declines all responsibility for damage resulting from failure to comply with the instructions for use of the product. The equipment must be used as described in the user manual, with particular regard to the intended use, using only accessories as specified by the manufacturer. No modification or repair can be performed.

1.7 Disclaimer

No part of this manual may be reproduced in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of CAEN spa.


The information contained herein has been carefully checked and is believed to be accurate; however, no responsibility is assumed for inaccuracies. CAEN spa reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it.

1.8 Made in Italy

We remark that all our boards have been designed and assembled in Italy. In a challenging environment where a competitive edge is often obtained at the cost of lower wages and declining working conditions, we proudly acknowledge that all those who participated in the production and distribution process of our devices were reasonably paid and worked in a safe environment (this is true for the boards marked "MADE IN ITALY", while we cannot guarantee for third-party manufactures).



Table of Contents

1.1	Purpose of this Guide	2
1.2	Change Document Record.....	2
1.3	Symbols, abbreviated terms and notation	2
1.4	Reference Documents	2
1.5	Manufacturer Contacts	2
1.6	Limitation of Responsibility	3
1.7	Disclaimer	3
1.8	Made in Italy	3
Table of Contents		4
List of Figures		5
List of Tables		5
1	Safety Notices	6
	Carefulness with Radioactive Sources	8
2	Introduction	9
3	PID (Product Identifier).....	12
4	Hardware Description.....	14
	SP5600 - Power Supply and Amplification Unit	14
	DT5720A - Desktop Digitizer	14
	SP5606 - Mini Spectrometer	15
	A315 - Splitter	17
	SP5607 - Absorption tool	17
5	Getting started	19
	5.1 Software Installation Requirements.....	19
	5.2 Software Installation	24
6	Software Description	25
7	Basic Measurements	47
	7.1 Can you see γ spectrum?	47
	7.1.1 First Kit Configuration	47
	7.1.2 Second Kit Configuration.....	51
	7.1.3 Correctly bias the SiPM and set the gain	54
	7.1.4 External and internal trigger	55
8	Educational Experiments	59
	8.1 Detecting γ -radiation (SG6111A).....	60
	8.2 Poisson and Gaussian Distribution (SG6112A)	61
	8.3 Energy Resolution (SG6113A)	62
	8.4 System Calibration: Linearity and Resolution (SG6114A)	63
	8.5 A comparison of different scintillating crystals: Light Yield , Decay Time and resolution (SG6115A)	64
	8.6 γ -Radiation Absorption (SG6116A)	65
	8.7 Photonuclear cross-section/Compton Scattering cross-section (SG6117A)	66
	8.8 Study of the ^{137}Cs spectrum: the backscatter peak and X rays (SG6118A)	67
9	Appendix.....	68
10	Instructions for Cleaning.....	85
	10.1 Cleaning the Touchscreen	85
	10.2 Cleaning the air vents.....	85
	10.3 General cleaning safety precautions	85
11	Device decommissioning	86
12	Disposal	87
13	Technical Support.....	88

List of Figures

Fig. 3.1:PID location taking a CAEN Educational kit as an example (the number in the picture and the device model are purely indicative).	12
Fig. 3.2: The PID position is located on the plastic black case for SP5606 - Mini Spectrometer, SP5607 - Absorption Tool and SP5608 - Scintillating Tile.	12
Fig. 3.3: SP5600 - Power Supply and Amplification Unit: PID position is located on the back panel of the module hosting the power input. Same location is adopted for DT5720A - Desktop Digitizer also.	13
Fig. 3.4: PID position: on the metal shield for the A315 . Splitter or around the holder for the SP5650C - Sensor Holder with SiPM.	13
Fig. 5.1: Tracking the PSAU port assignment on a PC running Windows 10.	24
Fig. 6.1: Main GUI of the HERA software.	25
Fig. 6.2: Main GUI Description.	26
Fig. 6.3: About window.	28
Fig. 6.4: Bias Setting Window.	30
Fig. 6.5: Selection by Experiment.	31
Fig. 6.6: Example of experimental activity.	31
Fig. 6.7: HERA: Selection by kit.	32
Fig. 6.8: Hardware Management.	33
Fig. 6.9: PSAU General Tab.	42
Fig. 6.10: PSAU Channels Setting Tab.	43
Fig. 6.11: PSAU Commons Setting Tab.	43
Fig. 7.1: Power Supply & Amplification Unit Interactive Panel – Gamma spectroscopy settings.	50
Fig. 7.2: Scope trace of the spectrometer.	50
Fig. 7.3: Dark Counts Rate frequency versus Discriminator threshold.	51
Fig. 7.4: Scope trace of the spectrometer with a ²² Na source.	54
Fig. 7.5: The DIGITIZER control panel.	55
Fig. 7.6: Wave for a right PSAU trigger threshold.	55
Fig. 7.7: Frequency vs Threshold for SiPM 6x6 mm ² of spectrometer, in a run without crystal, with crystal and with Crystal and radioactive source.	56
Fig. 7.8: The Trigger digitizer parameters.	56
Fig. 7.9: The Waveform tab of the HERA Main GUI.	57
Fig. 7.10: ²² Na spectrum.	57







List of Tables

Tab. 2.1: Building blocks of the kits.	10
Tab. 4.1: BGO Scintillator Properties.	15
Tab. 4.2: LYSO(Ce) Scintillator Properties.	16
Tab. 4.3: CsI(Tl) Scintillator Properties.	16
Tab. 4.4: S13360-6050CS Features.	16
Tab. 4.5: Typical properties of the optical grease.	17
Tab. 5.1: Host PC requirements.	19
Tab. 6.1: PSAU library return codes.	42
Tab. 6.2: Digitizer library return codes.	44
Tab. 8.1: Physics Experiments performed via the Educational Kit – Premium version.	59
Tab. 9.1: Waveforms saving scheme.	69
Tab. 9.2: Histograms saving scheme.	70
Tab. 9.3: Two channels saving scheme.	72
Tab. 9.4: Charge vs Time saving scheme.	73
Tab. 9.5: Counting saving scheme.	74
Tab. 9.6: Physical structure of the .dat files. Each coloured box indicates a single byte (8 bits).	77
Tab. 9.8: Data saving scheme of the Experiments Section B1.	78


1 Safety Notices

N.B. Read carefully the “Precautions for Handling, Storage and Installation” document provided with the product before starting any operation.

The following HAZARD SYMBOLS may be reported on the unit:

	Caution, refer to product manual
	Caution, risk of electrical shock
	Protective conductor terminal
	Earth (Ground) Terminal
	Alternating Current
	Three-Phase Alternating Current

The following symbol may be reported in the present manual:

	General warning statement
---	---------------------------

The symbol could be followed by the following terms:

- **DANGER:** indicates a hazardous situation which, if not avoided, will result in serious injury or death.
- **WARNING:** indicates a hazardous situation which, if not avoided, could result in death or serious injury.
- **CAUTION:** indicates a situation or condition that, if not avoided, could cause physical injury, or damage the product and / or its environment.

CAUTION: To avoid potential hazards



**USE THE PRODUCT ONLY AS SPECIFIED.
ONLY QUALIFIED PERSONNEL SHOULD PERFORM SERVICE PROCEDURES**

CAUTION: Avoid Electric Overload



**TO AVOID ELECTRIC SHOCK OR FIRE HAZARD, DO NOT POWER A LOAD
OUTSIDE OF ITS SPECIFIED RANGE**

CAUTION: Avoid Electric Shock



**TO AVOID INJURY OR LOSS OF LIFE, DO NOT CONNECT OR DISCONNECT
CABLES WHILE THEY ARE CONNECTED TO A VOLTAGE SOURCE**

CAUTION: Do Not Operate without Covers



**TO AVOID ELECTRIC SHOCK OR FIRE HAZARD, DO NOT OPERATE THIS
PRODUCT WITH COVERS OR PANELS REMOVED**

CAUTION: Do Not Operate in Wet/Damp Conditions



**TO AVOID ELECTRIC SHOCK, DO NOT OPERATE THIS PRODUCT IN WET
OR DAMP CONDITIONS**

CAUTION: Do Not Operate in an Explosive Atmosphere



**TO AVOID INJURY OR FIRE HAZARD, DO NOT OPERATE THIS PRODUCT
IN AN EXPLOSIVE ATMOSPHERE**

Do Not Operate with Suspected Failures. If you suspect this product to be damaged, please contact Technical Support.



**THIS DEVICE SHOULD BE INSTALLED AND USED BY SKILLED TECHNICIAN
ONLY OR UNDER HIS SUPERVISION**



**DO NOT OPERATE WITH SUSPECTED FAILURES.
IF YOU SUSPECT THIS PRODUCT TO BE DAMAGED, PLEASE CONTACT
THE TECHNICAL SUPPORT**



**THE SAFETY OF ANY SYSTEM THAT INCORPORATES THE DEVICE IS UNDER
THE RESPONSIBILITY OF THE ASSEMBLER OF THE SYSTEM**

See Chap. 13 for the Technical Support contacts.

Carefulness with Radioactive Sources

The Physics experiments related to Beta spectroscopy proposed in this manual needed radioactive sources.

There are two radioactive source types for educational purpose: sealed and unsealed sources. In the following experiments, sealed sources have been used. This source type is typically easier to use because the radioactive material is deposited in a plastic disk and sealed inside with a durable epoxy. Problems related to possible spills or decontamination are negligible.

Sealed gamma or beta sources of low activity, such as 0,1 μCi or a little bit more, can be handled directly without significant risk, although it is good practice to utilize tongs. Otherwise, sealed gamma sources with high activity, such as 10 μCi or more, should only be handled with tongs.

Nevertheless, when working with radioactive sources, mitigation of radiation exposures is very important. The basic principles of ALARA can give instructions. ALARA (As Low As Reasonably Achievable) is a radiation safety principle for minimizing radiation doses and releases of radioactive materials by employing all reasonable methods. ALARA is not only a sound safety principle for all radiation safety programs but is a regulatory requirement.

The three main principles are related to:

- Time: minimizing the time of exposure is the simplest way to directly reduce radiation exposure;
- Distance: doubling the distance between the radiation source and human body means to reduce radiation exposure by a factor of 4;
- Shielding: using absorber materials, such as lead for X-rays and gamma rays and Plexiglas for beta particles, to reduce the radiation reaching the body from a radioactive source is an effective way to reduce radiation exposures.

The radioactive sources for educational purpose have a low level of activity and their storing is a relatively simple matter. Solid sealed sources can be safely stored in their own plastic containers of shipment and then they can be put together in a locked cabinet, possibly with an additional shielding of lead sheets or bricks.



Important Note: Beta and Gamma Radioactive Sources are not included in the SP5600AN Educational Kit – Premium Version.

2 Introduction

CAEN brings the experience acquired in more than 40 years of collaboration with the High Energy & Nuclear Physics community into the University educational laboratories. Thanks to the most advanced instrumentation developed by CAEN for the major experiments Worldwide, together with the University teaching experience at the University of Insubria, a series of experiments covering several applications has been carried out.

CAEN realized different modular Educational Kits. The set-ups are all based on Silicon Photomultipliers (SiPM) state-of-the-art sensor of light with single photon sensitivity and unprecedented photon number capability.

The **Educational Gamma Kit, SP5600C**, is the system solution to perform several gamma spectroscopy experiments based on Silicon Photomultipliers.









The Educational Gamma Kit comprises:

- Nr. 1 Power Supply & Amplification Unit (PSAU, ID code SP5600). The PSAU supplies the bias for the sensors, features a variable amplification factor up to 50 dB and integrates a feedback circuit to stabilize the sensor gain against temperature variations. Moreover, the PSAU includes one leading edge discriminator/channel and a coincidence circuit for flexible event trigger logic. Sensors housed in dedicated mechanical holders can be directly connected to the PSAU. The PSAU technical specifications are reported in the relevant data sheet, together with the front and rear panel description.
- Nr. 1 Desktop Waveform digitizer (ID code DT5720A), with 2 input channels sampled at 250 MS/s by a 12-bit ADC. The DT5720A runs the Digital Pulse Processing for enhanced triggering and integration capabilities. The Digitizer technical specifications are reported in the relevant User's manual, together with the front and rear panel description.
- Nr. 1 Mini-Spectrometer for gamma ray detection (ID code SP5606) composed of a mechanical structure that houses a scintillating crystal, coupled to a dedicated SiPM. Three different crystals are available: CsI, LYSO and BGO. The spectrometer is equipped with a bottom support to allow an easy connection to the SP5600 via the splitter A315, to avoid saturation effects. The spectrometer head was designed to make easy a full range of basic gamma spectrometry educational experiments. The main features of the spectrometer are reported in the data sheet.
- Nr.1 The Gamma absorption tool (ID code SP5607) allows to perform gamma attenuation measurements. It is a modular tool, and its design allows an easy connection to the SP5606 bottom support. The SP5607 is composed of spacers, Aluminium and PMMA absorbers (one 4 mm thick, five 10 mm thick).
- Nr.1 The Splitter (ID code A315) divides one input on two output signals. The splitter is used to connect the SP5606 to the SP5600 by avoiding saturation effects. The A315 technical specifications are reported in the data sheet.
- External AC/DC stabilized 12 V power supplies (Meanwell GS40A12-P1J 40 W, 12 V DC Output, 3.34 A).
- Nr.1 Kit cables (ID code FKITSP56) composed of: n.1 LEMO-LEMO cable, n.2 MCX-LEMO cables, n.1 MCX-MCX cables, n.1 Power Cord Adapter (1IN / 3 OUT).
- Nr.1 Pen-Vac Vacuum Pickup Tool (VPV) is an ideal tool for manually offloading absorbers from SP5607. The pickup tool is a self-contained unit and can lift up to 500 grams.
- Nr.1 Optical grease.
- USB cables.
- A LabView™ based software: HERA (Handy Educational Radiation Application).

The different building blocks of the kit can be assembled in a customized configuration, according to the specific application and the user's requirements.

The purpose of this guide is to provide a hands-on primer on the use of the essential functionalities of the kit.



Item description	Code	Image	SP5600AN Educational Premium kit	SP5600C Educational Gamma kit	SP5600D Educational Beta kit	SP5600E Educational Photon kit
SP5600 Power Supply and Amplification Unit	WSP5600XAAAA		yes	yes	yes	yes
DT5720A Desktop Digitizer	WDT5720XAAAA		yes	yes	yes	yes
SP5601 Led Driver	WSP5601XAAAA		yes	no	no	yes
SP5650C Sensor Holder with SiPM	WSP5650XCAAA		yes	no	no	yes
SP5606 - Mini Spectrometer	WSP5606XAAAA		yes	yes	no	no
SP5607 Absorption Tool	WSP5607XAAAA		yes	yes	no	no
SP5608 Scintillating Tile	WSP5608XAAAA		yes	no	yes	no
A315 - Splitter	WA315XAAAAAA		yes	yes	no	no

Tab. 2.1: Building blocks of the kits.

CAUTION: to manage the product, consult the operating instructions provided.

When receiving the unit, the user is strictly recommended to:

- Inspect containers for damage during shipment. Report any damage to the freight carrier for possible insurance claims.
- Check that all the components received match those listed on the enclosed packing list as in Tab. 2.1. (CAEN cannot accept responsibility for missing items unless we are notified promptly of any discrepancies.)
- Open shipping containers; be careful not to damage contents.
- Inspect contents and report any damage. The inspection should confirm that there is no exterior damage to the unit such as broken knobs or connectors and that the front panel and display face are not scratched or cracked. Keep all packing material until the inspection has been completed.
- If damage is detected, file a claim with carrier immediately and notify CAEN service.
- If equipment must be returned for any reason, carefully repack equipment in the original shipping container with original packing materials if possible. Please, contact CAEN service.
- If equipment is to be installed later, place equipment in original shipping container and store in a safe place until ready to install



DO NOT SUBJECT THE ITEM TO UNDUE SHOCK OF VIBRATIONS



DO NOT BUMP, DROP OR SLIDE SHIPPING CONTAINERS



DO NOT LEAVE ITEMS OR SHIPPING CONTAINERS UNSUPERVISED IN AREAS WHERE UNTRAINED PERSONNEL MAY MISHANDLE THE ITEMS



USE ONLY ACCESSORIES WHICH MEET THE MANUFACTURER SPECIFICATIONS

3 PID (Product Identifier)

PID is the CAEN product identifier, an incremental number greater than 10000 that is unique for each product¹. The PID is on a label affixed to the kit suitcase by the opening mechanism (Fig. 3.1).



Fig. 3.1: PID location taking a CAEN Educational kit as an example (the number in the picture and the device model are purely indicative).

The PIS is even stored on each educational kit subparts as shown in Fig. 3.2, Fig. 3.3 and Fig. 3.4.



Fig. 3.2: The PID position is located on the plastic black case for SP5606 - Mini Spectrometer, SP5607 - Absorption Tool and SP5608 - Scintillating Tile.

¹ The PID substitutes the serial number previously identifying the boards.



Fig. 3.3: SP5600 - Power Supply and Amplification Unit: PID position is located on the back panel of the module hosting the power input. Same location is adopted for DT5720A - Desktop Digitizer also.



Fig. 3.4: PID position: on the metal shield for the A315 . Splitter or around the holder for the SP5650C - Sensor Holder with SiPM.



Note: The serial number is still valid to identify older boards, where the PID label is not present.

4 Hardware Description

As previously mentioned, the Educational Gamma Kit is composed of a lot of hardware devices. To better understand their use, this section provides a brief description of the main kit units.

SP5600 - Power Supply and Amplification Unit



- Variable amplification gain (up to 50 dB)
- Low noise, to guarantee high performances of the sensor even with small signals
- Wideband, to comply with the fast sensor response
- Fast leading-edge discriminator and time coincidence
- Provides the bias for the sensors with gain stabilization
- USB 2.0 interface
- Dimension: 150 x 50 x 70 mm³ (WxHxD)

The SP5600 is a general-purpose Power Supply and Amplification Unit, integrating up to two SiPMs in a mother & daughter architecture allowing easy mounting and replacement of the sensors. The basic configuration features two channels with independent gain control up to 50 dB and provides the bias voltage (up to 130 V) to the sensors with gain stabilization. Each channel can provide a digital output generated by the fast-leading edge discriminators. A timing coincidence of the two channels is also available [RD5].

DT5720A - Desktop Digitizer



- 2 Channel 12 bit 250 MS/s Digitizer
- Digital Pulse Processing for Charge Integration DPP-CI for SiPM
- Best suited for PMT and SiPM/MPPC readout at low and high rates
- Mid-High speed signals (Typ: output of PMT/SiPM)
- Good timing resolution with fast signals (rise time < 100 ns)
- Optical Link and USB 2.0 interfaces
- Dimension: 154 x 50 x 164 mm³ (WxHxD)

The DT5720A is a 2 Channel CAEN Waveform Digitizers able to perform waveform recording and run online advanced algorithms of charge integration (DPP-CI), i.e. the digital version of the traditional QDC (Charge-to-Digital Converter) [RD4].

Data is read by a Flash ADC, 12-bit resolution and 250 MS/s sampling rate, which is well suited for mid-fast signals as the ones coming from liquid or inorganic scintillators coupled to PMTs or Silicon Photomultipliers. The acquisition can be either channel independent or common through an external signal and the acquired data can be saved for offline analysis.

The acquisition in DPP-CI mode for SiPM is fully controlled by the Hera software, which manages the algorithm parameters, builds plots and saves the relevant information through the USB 2.0 interface of the digitizer (data transfer up to 30 MB/s).

The digitizer runs on real time:

- Self-Trigger using CR-RC digital Time filter algorithm
- Input signal baseline (pedestal) calculation
- Charge Integration (with programmable gate parameters) with pedestal subtraction for energy calculation.

SP5606 - Mini Spectrometer



- Mechanical structure for optimal SiPM to crystal coupling
- Included crystals: LYSO, BGO, CsI
- Crystals dimension: 6 x 6 x 15 mm³
- One SiPM 6 x 6 mm² embedded

SP5606 is a mini spectrometer for gamma ray detection. The spectrometer is composed of a mechanical structure that houses a scintillating crystal, coupled to a dedicated 6 x 6 mm² SiPM. Three different crystals are available: CsI, LYSO and BGO. The spectrometer is equipped with a bottom support for an easy connection to the SP5600 via the splitter A315, to avoid saturation effects.

More details about the mini spectrometer are listed here below.

Scintillating Crystals

The SP5606 is provided with three different crystals: BGO (Bismuth Germanate [Bi₄Ge₃O₁₂]), LYSO(Ce) (Cerium-doped Lutetium Yttrium Orthosilicate), CsI(Tl) (Thallium-doped Cesium Iodide). The crystals are polished on all sides and coated with white epoxy on 5 faces. One 6 x 6 mm² face is open in order to be coupled with the Silicon Photomultiplier.

○ BGO information:

A relatively hard, high density, non-hygroscopic crystal with good gamma ray absorption. Often used for PET imaging and high energy physics applications as Compton shields. The main characteristics of the crystal are summarized in the Tab. 4.1.

Properties	Value
Cleavage Planes	None
Decay Constant (ns)	300
Density (g cm ⁻³)	7.13
Emission Spectral Range (nm)	350-650
Melting Point (K)	1323
Peak Scintillation Wavelength (nm)	480
Photons/MeV	8500
Radiation Length (cm)	1.13
Refractive Index at Peak Emission	2.15
Solubility (g/100g H ₂ O @ 300K)	Insoluble
Stability	Good
Structure	Cubic
Transmission Range (nm)	470-7500

Tab. 4.1: BGO Scintillator Properties.

○ LYSO(Ce) information:

Non-hygroscopic scintillator that is both bright and fast. It is often employed in applications where fast timing is needed such as PET and TOF PET. The main characteristics of the crystal are summarized in the Tab. 4.2.

Properties	Value
Cleavage Planes	None
Decay Constant (ns)	40
Density (g cm ⁻³)	7.1
Emission Spectral Range (nm)	380-480
Melting Point (K)	2323
Peak Scintillation Wavelength (nm)	420
Photons/MeV	32000
Radiation Length (cm)	1.15
Refractive Index at Peak Emission	1.81
Solubility (g/100g H ₂ O @ 300K)	Insoluble
Stability	Good
Structure	Cubic

Tab. 4.2: LYSO(Ce) Scintillator Properties.

○ CsI(Tl) information:

This scintillator offers a high light yield and emits at a wavelength very suitable for silicon photomultipliers (SiPMs). Typical applications include arrays of this material used in security imaging systems, such as baggage scanners. The main characteristics of the crystal are summarized in the Tab. 4.3.

Properties	Value
Cleavage Planes	None
Decay Constant (ns)	1000
Density (g cm ⁻³)	4.51
Emission Spectral Range (nm)	350-725
Gamma and X-ray absorption coefficients (cm ⁻¹)	0.48 at 660keV 10.00 at 100KeV
Melting Point (K)	894
Peak Scintillation Wavelength (nm)	550
Photons/MeV	52000
Radiation Length (cm)	1.86
Refractive Index at Peak Emission	1.78
Solubility (g/100 g H ₂ O @ 300 °K)	44.0
Stability	Slightly Hygroscopic
Structure	BCC
Thermal Conductivity (W·m ⁻¹ ·K ⁻¹) @ 300K	1.13
Transmission Range (nm)	240-70000

Tab. 4.3: CsI(Tl) Scintillator Properties.

Silicon Photomultiplier

The embedded detector is the Hamamatsu Silicon Photomultiplier (SiPM), S13360-6050CS. A small printed circuit board inside the mini-spectrometer hosts the 6 x 6 mm² detector.

The main features of the S13360-6050CS SiPM are listed in Tab. 4.4.

Properties	Value
Package type	Ceramic
Active Area	6 × 6 mm ²
Number of pixels	14400
Pixel Pitch	50 μm
Spectral response range	270 to 900 nm
Peak sensitivity wavelength (typ.)	450 nm
Dark count/ch (typ.)	2000 kcps
Terminal capacitance/ch (typ.)	1280 pF
Gain (typ.)	1.7×10 ⁶
Measurement condition	Ta=25 °C

Tab. 4.4: S13360-6050CS Features.


Optical Grease

The optical coupling grease is a non-curing colourless coupling gel, clear and colourless having moderate viscosity and providing excellent transmission properties well into the near-ultraviolet region. It should be stored at temperatures below 26 °C, preferably above 5 °C, but it retains clarity and fluid property down to -60 °C.

Typical Properties	Value	
Colour	Clear	
Refractive Index @25 °C	1.466	
Specific Gravity	1.06	
Penetration	300	
Light Transmittance @ 1 cm	300 nm	99.45%
	400 nm	99.99%
	450 nm	99.99%
	500 nm	99.99%
	633 nm	99.99%
	850 nm	99.99%
	1310 nm	99.65%
	1550 nm	99.38%

Tab. 4.5: Typical properties of the optical grease.

A315 - Splitter



The Mod. A315 splits one input on two output signals. All the connectors are LEMO female type. The splitter is adapted for 50 Ohm lines to avoid reflections of the signal. The device is completely passive (no power supply is required); the amplitude on each output is one half of that on the input.

SP5607 - Absorption tool



- Spacers: one 4mm thick, five 10 mm thick
- Aluminium Absorbers: one 4mm thick, five 10 mm thick
- PMMA Absorbers: one 4mm thick, five 10 mm thick

The Gamma absorption tool allows performing gamma attenuation measurements. It is a modular tool, and thanks to its design, it can be easily connected to the SP5606 bottom support. This tool includes several spacers and two different absorber materials: Aluminium and Polymethyl methacrylate (PMMA).

As well known, the Aluminium is lightweight, durable, malleable, and corrosion-resistant metal and it is widely used in the aerospace, transportation, and construction industries.

Out of scientific community, the PMMA is more known as Plexiglas and acrylic. The first use of PMMA as a dental device was for the fabrication of complete denture bases. Its qualities of biocompatibility, reliability, relative ease of manipulation, and low toxicity were soon seized upon and incorporated by many different medical specialties. PMMA has been used for bone cement, contact and intraocular lens, screw fixation in bone, filler for bone cavities and skull defects. Moreover, it is widely used in medical dosimetry as a water equivalent solid-state organic material. Water and PMMA are found to be the closest soft-tissue and water substitutes, respectively. Soft-tissue and water equivalence of dosimetry materials need to be evaluated for a range of photon energies and field sizes before their application in complex radiation beams. The behaviour of the PMMA under the action of gamma radiation is generally studied to measure the absorbed dose rate to water using a PMMA phantom.


The following table summarizes additional information on Aluminium and PMMA.

Material	Mass Density (g cm^{-3})	$(Z/A)_{\text{eff}}$	Mole fraction (%)
PMMA	1.185	0.539	H 53.333, C 33.333, O 13.333,
Aluminium	2.698	0.482	Al 100

5 Getting started

This chapter will guide you through the drivers installation of PSAU and Digitizer, as well as the installation of HERA (Handy Educational Radiation Application) software and the first measurements.


5.1 Software Installation Requirements

OS	Hardware	CAEN drivers required
 Microsoft Windows 10 (64-bit)	2 available USB2.0 ports	DT5720 USB driver (32/64-bit) SP5600 USB driver (32/64-bit)


Tab. 5.1: Host PC requirements.

- Download the USB drivers for both DT5720 and SP5600 compliant to the Windows version 64-bit on CAEN website: Educational kit webpage > “Download” > “Software” tab > Driver section (login is required before the download).
- Install the DT5720 drivers** following the instruction of the setup wizard. The OS will automatically recognize the DT5720 when it is connected to the PC. If the automatic installation fails, perform it manually from the Device Manager by selecting the driver update and pointing to the driver folder you downloaded from CAEN website.


For example, once connected and powered on the digitizer, you can do it going to Control Panel -> System & Security -> System -> Device Manager. In the Device Manager window, find the unknown **CAEN DT5xxx USB 1.0** in the list **Other Devices**:




Right click on **CAEN DT5xxx USB 1.0** and select **Update Driver Software** option in the scroll menu.




Select **Browse My Computer** for driver software.



Click **[Browse]** to point to the Windows drivers' folder you have previously unpacked, click **[OK]** to include the path in the search and click **[Next]** to continue.



When the driver installation will be completed, click **Close** to close the window.




Refer to [RD3] for detailed installation OS-dependent.


- Connect to the PC and power ON the **SP5600**; the PC will recognize as a new peripheral by the OS. Perform the driver installation manually from the Device Manager by selecting the driver update and pointing to the driver folder you downloaded from CAEN website. Finally, a COMM port will be associated to SP5600.

For example (Windows 10 – 64bit), once connected and powered on the SP5600, you can follow the previous instructions going to Control Panel -> System & Security -> System -> Device Manager -> Controller USB [Ports (COM)] Manager.


Right click on **USB Serial Device** and select **Update Driver Software** option in the scroll menu.




Select **Browse My Computer** for driver software.



Click **[Browse]** to point to the Windows drivers' folder you have previously unpacked, click **[OK]** to include the path in the search and click **[Next]** to continue.



When the driver installation will be completed, click **Close** to close the window.



Finally, a COM port will be associated to SP5600; please check the port number as shown in Fig. 5.1.




Fig. 5.1: Tracking the PSAU port assignment on a PC running Windows 10.




Important Note: HERA rel. 1.0.0 Build: 1.5.21.0103 or higher:

- does not require LabVIEW™ Run-Time Engine. or LabVIEW™ version 2018 (or higher). The installation of LabVIEW™ Run-Time Engine 2018 is already implemented in the HERA.
- does not work with a digitizer USB Driver release < 3.4.7, if running in a 32-bit Windows environment.

5.2 Software Installation

Download the standalone HERA Software full installation package on CAEN website: Educational kit webpage > “Download” > “Software” tab > Application SW section (login is required before the download).

Unpack the installation package, login as administrator, launch the setup file, and complete the Installation wizard.



The setup automatically creates a link on the PC Desktop.

6 Software Description

When the installation procedure has been completed, the user can run the program by clicking the correspondent icon.

HERA (Handy Educational Radiation Application) is a user-friendly software platform allowing the user to manage the following CAEN Educational kits: SP5600E – Photon kit, SP5600D – Beta kit, SP5600C – Gamma kit, SP5600AN – Premium Version kit.

The simple graphical interfaces help the user to perform its own experimental activity. As shown in the opening window in Fig. 6.1, several ways of operative openings are available. This initial access multiplicity makes the software very flexible and suitable both for expert users as well as for beginner ones. With a simple selection, the user can decide how to execute the activity by choosing the direct access to the suggested experiments or access to devices management.

Via this main GUI, it is possible to visualize the devices status, server messages and, to access to data (log file, data stored, configuration files, etc.).





Fig. 6.1: Main GUI of the HERA software.

Before running the software, the user should wait the hardware connection. The software recognises the hardware automatically and start the connection. Two connection indicators, "Online Hardware", are present on the opening window:

- **Green light** means that the connection is ok.
- **Red light** means that there is no connection.
If the PSAU is power on, but the light colour is red, the software can be forced to search for a new connection via the rescanning procedure from the Verbose Menu: Hardware-> PSAU Server Scanning -> Restart Scanning.



- **Yellow light** means that either the DGTZ is not a DT5720A/C, or its firmware is not compliant with Hera software, and another firmware type is probably running on the board.

The special firmware compatible with HERA Software is the Digital Pulse processing for Charge Integration for SiPM Kit (DPP-CI for SiPM) for DT5720A and the Digital Pulse Processing for Charge Integration and Pulse Shape Discrimination (DPP-PSD) for DT5720C [COMING SOON]. The firmware can be download from CAEN Website. Without any licenses, it will run in a 30-minute-per-power-cycle fully functional trial version.

To upload the firmware on the digitizer, use the CAENUpgrader Software (free download on CAEN Website): <https://www.caen.it/products/caenupgrader/>





Fig. 6.2: Main GUI Description.

The main features of the GUI are:

- **Timer/Watch:** Time and Date.
- **Online Hardware:** Indicators of the Digitizer (DGTZ) and Power Supply and Amplification Unit (PSAU) status.
- **User Interface Messages:** Operation messages related to user activity.
- **Server Messages:** Messages related to the hardware server activity.
- **Sub-panel shell:** Sub-panel of the initial menu choice.
- **Online Help:** QuickStart guide available for each software window and experimental activity.
- **Verbose Menu:** It is organized into several items (File, View, User & Administration, Activities, Hardware, About menus), each one allowing the user to perform several actions.

– File Menu


The user can assess further functionalities by pressing the “File” label on the top left of the GUI. As shown in the picture below, the File Menu is composed by four sections.



The first section gives access to the logfiles and to a special tool for the data file conversion in .txt format. The second section allows the user to save and retrieve a configuration file containing the parameters settings for the DGTZ and the PSAU and, moreover, the default configurations recovery. The third section is focused on data storage management. The user can open the data directory located in ProgramData folder or copy it in another PC location. The fourth section allows closing the software.


– View Menu

Through this menu, the user can enable/disable the display of the server messages in the bottom part of the window.



– User & Administration Menu

This menu section gives high flexibility to the user in managing the accounts and deciding what each user can access of the software.




For this purpose, two interesting functionalities are available:

User Account

The software can create several accounts via the "User Login" selection. This procedure is advantageous when several people or groups work on the same computer. Once created the account, the software automatically produces the related folder in which data, configuration files, and images can be saved and stored.

Via “Logout” selection, the account goes back to the default "Guest" user.



Software Access Mode

Three access levels to the software are implemented and are available through the “Limit Access...” selection in the drop-down menu. The first one, “Level 1”, just gives access to the Hardware Management. The second one, “Level 2”, allows the user to access the Hardware Management and the guided procedures to perform the experiments listed in the CAEN Educational Handbook. This access level does not include analysis tools. The third one, “Level 3”, gives full access to all software functionalities and all the analysis tools are included.

The initial option, “Selection by Kit”, is accessible by all the three access levels.

The user needs the Master Password to change the access mode. The Password is unique, not changeable, and not declared in the embedded Help to give this type of modification power to the tutors only.

This functionality allows to the tutors deciding what each user can access and therefore, structuring the courses depending on the course attendee levels.

The Master Password is the build of the HERA release in use. The build is displayed in the "About HERA" window via the about label in the GUI verbose menu.





Fig. 6.3: About window.


– Activities Menu

The “Activities” drop-down menu allows the user to direct access to the experimental activities and to the hardware control panels. Moreover, through it, it is possible the use the analysis tools without connecting to the devices.



– Hardware Menu

The “Hardware” drop-down menu is made of three sections useful for the hardware management. The first section gives access to the PSAU Control Panel and to launch again the rescanning procedure when its connection is lost.



Through the second section it is possible to select the digitizer model (DT5720A or DT5720C) and to access the DGTZ Control Panel.

The last section is very important for its preventing action related to the possible detectors damaging:

Detector Safety

The Bias Voltage Limits can be modified from the main window of the GUI only, before selecting the experiment or hardware. The user can set the detector safety condition via the “Hardware” drop-down menu.

This functionality is very important for preventing action related to the possible detectors damaging. The SP5600 module houses two detectors (SiPM). The module provides independent bias voltages (up to 130 V) to the sensors with gain stabilization. The user can apply a safety measure to prevent detector damage due to a wrong and too high bias voltage. Via “Safety and Serial Number Setting...” selection, the user can set the recommended operating voltage for each channel and, discretionary, the serial number to identify the detector itself. The software stabilizes the maximum value of bias voltage that can be applied to the sensor as a percentage (2,5%) of the operating one. To change the voltage limit is requested to modify the value in the “Bias Setting window” (see Fig. 6.4).

– About

The “About” leads to a new window including all information related to the software (release, build, etc).




Fig. 6.4: Bias Setting Window.

The Main GUI clearly shows that several ways of operative openings are available:

- **“Selection by Experiment”**: access to experiments frame covering Nuclear and Particle Physics fields.
- **“Selection by Kit”**: access to operative options allowed by the educational kit in use.
- **“Hardware Management”**: direct access to the management of the device parameters and data readout.

The chosen option can be run by double clicking on the relative box or by selecting it and then by the press on the “Select” button.

The user can easily access to the GUI description via the “Help” button. Each window of the software is equipped with a dedicated “Help” button that must be closed before starting any activity.

Selection by Experiment

This option allows the user to access the experiment menu listed in the CAEN Educational Handbook. By selecting the Physics topic of interest, a series of experiments can be performed. The software programs a predefined settings of the devices and gives a detailed guide into the “Help” button. The option “Selection by Experiment” can be run by double click on the relative icon or by selecting it and then by pressing the “Select” button.




Fig. 6.5: Selection by Experiment.

The “Help” button is present in all the windows in use and provides guides and advice about the experimental procedures.




Fig. 6.6: Example of experimental activity.

Selection by Kit

This option allows the user to access the experiment menu listed in the CAEN Educational Handbook.

By selecting the Physics topic of interest, a series of experiments can be performed. The software programs a predefined settings of the devices and gives a detailed guide into the “Help” button.

The option “Selection by Experiment” can be run by double click on the relative icon or by selecting it and then by pressing the “Select” button.

The “Help” button is present in all the windows in use and provides guides and advice about the experimental procedures.




Fig. 6.7: HERA: Selection by kit.

Hardware Management

The main units of the Educational Gamma kit, which are common among all the systems, are:

- Power Supply and Amplification Unit (PSAU) - SP5600
- Desktop Digitizer (DGTZ) - DT5720A

The “Hardware Management” section allows the user to manage all the parameters of both PSAU and DGTZ giving the highest flexibility in the operating modes.

With few easy steps, the setting of bias voltage, gain, thresholds, and digital outputs are possible. The digitized signals can be monitored for a real-time fine-tuning of the set-up. Energy spectra, trends of the charge as a function of the time, signal frequency versus threshold, and frequency counting are also displayed in the visualization tabs of the main GUI.




Fig. 6.8: Hardware Management.


All tabs allow the user to save plots and data on file for the offline analysis processes.

- Main GUI Description


The Main GUI of the HERA software consists of several visualization tabs. These tabs allow the user to visualize and manage the signals of the detector. The “Waveform”, “Histogram”, “Two Channels” and “Charge vs Time” tabs refer to the digitizer (DTGZ). The other two, “Frequency vs Threshold” and “Counting”, refer to the Power Supply and Amplification Unit (PSAU).

Control keyboard


The keyboard on the left side of the GUI allows the user to control and manage the acquisition tabs and to monitor the system status.



“Queue Backlog” indicates the number of acquired data (elements) that are waiting for displaying or saving to file, for both the "Waveforms" tab and "Histogram" tab. This element number should normally be equal to zero unless some extra time-consuming operation occurs.



Two *light indicators* provide the system status related to Data Acquisition and Storage.




“Histo Settings...” button makes available the selection of the different types of histograms for the “Histogram” tab.

“Settings...” button leads to an additional window providing the options for the channels enabling/disabling for both PSAU and DTGZ, and for setting the run preset according to a fixed time or number of events.



“Start” button must be used to launch the acquisition and to visualize the results on the related tabs. To stop the acquisition and/or change to another visualization tab, the user can press the *“Stop”* button. The single shot mode of the waveform can be activated via *“Pause”* button, then the *“Start”* button will change its name to *“Single shot”*. In that case, the plot will be updated and frozen with a single trigger. The continuous data stream can be activated again by pressing the *“Pause”* button.



The last three buttons of the keyboard are related to data storage. The *“Print”* button sends the result visualized in the tab to the selectable printer. The *“Export”* button opens an additional window to export data in two formats. If the "Clipboard" box is selected, a bitmap image to the Clipboard is exported. If the "Excel" box is selected, just numerical data are exported. *“Save on Time”* (or *“Save on #Events”*) button saves data according to the setting previously defined via *“Settings...”* button.

The display of the tabs is equipped with a Graph Palette that allows the user to interact with a graph.



This palette appears always with the following buttons, in order from the left to the right:


- Cursor Movement Tool moves the cursor on the display. If the cursor is not present, it does not work.
- Zoom acts by zooming in and out the display.
- Panning Tool picks up the plot and moves it around on the display.



Note: If one of the controls of the Graph Palette does not work, that interaction type is not permitted for the specified graph. The graph palette is always composed of the three controls, even if not all of them are used.

Waveform tab

The “Waveform” tab shows the traces of the analog and digital signals read out from the digitizer. The signals visualization can be enabled/disabled by selecting the related box on the legend on the right side.



The analog signals are the traces of the input channels (Ch.1 and Ch.0) and the virtual probe, i.e. the baseline signal.

The digital signals are the Gate (red), the Over Threshold (violet), the HoldOff (brown) and the Flat (yellow):

- “Gate” represents the width of the signal integration.
- “Holdoff” means the veto width for the generation of other gates.
- “Over Threshold” is generated when the signal is over the set threshold.
- “Flat” stands for the veto width for the baseline calculation.

The “Digital Signal View” section on the window right side includes the graphical controls for the digital traces. These traces can be amplified via the “Expansion” cursor and moved in a vertical direction via the “Offset” cursor. The “Record length” control allows the user to change the time scale of the acquisition window, from 4 μ s to 12 μ s.


The “Swap Trigger” button enables/disables the rising edge stabilization of the gate on the time scale.

The acquisition conditions for the data saving, previously set via the “Setting...” button on the control keyboard, can be applied simply by pushing the “Save on...” button during the acquisition run.


If the system is not in acquisition mode, the “Save on...” button only allows the data storage of the displayed waveform without any constraints in time or in events. The data will be saved in .TDMS format. The waveforms data format is described in detail in the Appendix.

Histogram tab

The "Histogram" tab shows the histogram of the active channel according to the PSAU and DGTZ settings.




Via "Histo Settings..." button, the user can enable/disable the general statistics on the display right side, choose the Y scale as linear or logarithmic, and the histogram types. HERA software supports two different histograms:



The "*Running Histogram*" accumulates data until the number of entries defined in the "*Number of points*" parameter is reached. The user can set this value in the Histogram Setting Window via "Histo Settings..." button of the Control keyboard. Data is overwritten by the new events. The number of bins, bin size and starting X point are automatically set by the software.

This kind of histogram processing is useful when the hardware conditions are changed during the measurement, to check how the system response evolves. For example, the user can try to change the L.E.D. intensity during the acquisition and observe how the histogram changes.



The "*Unlimited Histogram*" accumulates data with no limits in the number of entries of the Y-axis. Differently from the "Running Histogram", the user must provide the properties of the X scale in order to determine the histogram range:

- the *origin* of the histogram means the minimum plotted charge value;
- the *number of bins* determines the end of the plotting window;
- the *bin size*, i.e. bin width.

The Unlimited Histogram can be used to make comparisons between measurements taken in the same setup conditions. Note: the hardware setup must not be changed during the measurement.

The acquisition conditions for data saving, previously set via the "Setting..." button, can be applied through the "Save on..." button during the acquisition run.

If the system is not in acquisition mode, the "Save on..." button allows the data storage without any constraints in time or in events. The data will be saved in .TDMS or ASCII formats. The histograms data format is described in detail in the Appendix.



Note: Because of the automatic and variable setting of the bin size, the Running Histogram is not suitable for comparison purposes. This histogram type does not guarantee the same acquisition conditions. Conversely, the Unlimited Histogram is suitable for comparison among spectra due to its setting properties in terms of the number of bins, bin size and starting X point.


Two Channels tab

The “Two Channels” tab allows the user to manage the histogram plots from the two channels of the digitizer simultaneously. It is possible to plot and reset the two histograms, the histogram sum, and their correlation.

The *graph on the left side* contains the histogram plots of both channels and the sum plot. Each plot can be enabled or disabled through the relative box in the legend. The sum of the histograms is defined as the histogram resulting from adding channel0's histogram to channel1's histogram multiplied by an *alpha factor*. Common x-axis origin, number of bins and bin size can be set via “Histo Settings...” button for all the spectra. All graphs can be reset at the same time via the “Reset” button in the lower part of the window.

The *graph on the right side* shows a scatter plot of the signals from the two sensors, after being integrated in the specified time window (look at the displayed “Gate” in the *Waveforms tab*). This tab might help for specific applications relying on simultaneous use of the two detectors, e.g. when using the scintillator tiles for cosmic ray experiments or two spectrometry heads for ^{22}Na positron annihilation detection.


The histograms saving is described in detail in the Appendix.



Charge vs Time tab


The “*Charge vs time*” tab plots the signal charge versus time. The user can change the number of charges for the plotted mean. The plot can be stored pushing the “*Save on...*” button.

During the acquisition, the conditions for data saving, previously set via the “*Setting...*” button, can be applied simply by pushing the “*Save on...*” button. If the system is not in acquisition mode, the “*Save on...*” button allows the user to storage the data of the displayed plot without any condition related to the number of entries or acquisition time and the data will be saved in .TDMS format. The Charge vs Time data format is described in detail in the Appendix.




Frequency vs threshold tab

The “*Frequency vs Threshold*” tab allows the user to interact with the PSAU in order to produce the so-called “SiPM staircase”: the plot shows the frequency of the signals which are over the threshold, during a threshold scan from the minimum up to maximum threshold value.



After pressing the “*Start*” button on the control keyboard, the user can change the *limits of the scan*, the *step*, the number of read point which produces the *mean* plotted value and the *gate width* for the counting via an additional window, “*Frequency Scan Setting*”.




Once the acquisition is completed, the user can choose two threshold values and respectively write them in the two boxes, in the lower part of the window. These two threshold values will be transmitted to the “*Counting*” tab for further analysis.

To store the information included in this tab, the user can employ the “*Export*” button on the control keyboard.

Counting tab

The “Counting” tab shows three plots for which the user can change the number of points for the plotted mean value and the gate width for the counting.

The graph in the lower part plots the frequency trend of the signals over the threshold set in the PSAU Interactive Panel for the active channel or over the threshold value corresponding to the Low/High button selection, set in the previous tab. The two histograms show the distribution of the mean number of counts referred to the two threshold values (Low and High). The mean values of these distributions are displayed in the blue box, together with their uncertainty. The number of histogram entries and bins can be set by the controls “# of Entries” and “Histo # bins”.



The “Save on..” button allows the user to save the plots content at the end of the acquisition. The data format is described in detail in the Appendix.

- Power Supply & Amplification Unit (PSAU) Interactive Panel

The PSAU Interactive Panel is fully dedicated to the management of the Power Supply and Amplification Unit (PSAU). It is composed of three tabs: “General”, “Channels” and “Commons”.

General tab

The “General” tab contains “Board ID and Global Status” frame with information about the PSAU firmware release, Serial Number and COM Port. The “Last PSAU Hardware Error” frame shows the last Error Code of the library which the PSAU stands on. Moreover, the Temperature History plot shows the temperature of the board and of both two detectors.




Fig. 6.9: PSAU General Tab.

The PSAU library return codes displayed in the General tab are summarized in the following table.

Error code	Value	Meaning
PSAU_Success	0	Operation completed successfully
PSAU_InvalidComPortError	-1	Error related to the COM port
PSAU_TooManyClientsError	-2	Max. nr. of PSAUs simultaneously manageable exceeded
PSAU_CommunicationError	-3	Communication error
PSAU_InvalidHandleError	-4	Invalid device handler
PSAU_InvalidCommandError	-5	Unspecified error
PSAU_InvalidParameterError	-6	Invalid command error
PSAU_DeviceNotFound	-7	Invalid parameter error
PSAU_DeviceNotFound	-8	Device error (i.e., hardware or firmware issue)

Tab. 6.1: PSAU library return codes.

Channels tab

The "Channels" tab is composed of two sections: Setting and Monitor. The first one, on the left side, provides the switchers for the two channels enabling the settings of the bias voltage, the gain, the discriminators threshold, and the temperature compensation. The temperature compensation requires the setting of the coefficient “dV/dT” for both the channels. The compensation acts on the bias of the sensor to keep its gain constant, according to the voltage linear dependence as a function of the temperature. For both channels, the SiPM serial number is visualized according to the initial setting via “Safety and Serial Number Setting...” selection in the “Hardware” drop-down menu.

Two different graphical visualizations are provided to monitor the set parameters and verify channels status. An important note is shown to underline that the setting and monitoring of Bias Voltage have different resolutions due to the hardware.




Fig. 6.10: PSAU Channels Setting Tab.

Commons tab

The "Common" tab allows to user to set the width of both signals produced as digital outputs. The output level can be set as NIM or TTL standard and the polarity of the discriminator edge can be selected. The coincidence can be activated when both PSAU channels are switched on. The coincidence signal is provided on digital output of the selected channel and its width can be set in the Coincidence section of this tab.




Fig. 6.11: PSAU Commons Setting Tab.



Important Note: The "Update" button must be selected for all settings change to apply them correctly. The lack of this operation leaves the default settings unchanged. The "Auto Update" button executes the updating process automatically.

- Desktop Digitizer Interactive Panel

The Digitizer Interactive Panel allows the user to:

- check the digitizer connection and status (online/offline)
- check the model and serial number, revision, and firmware of the device
- overwrite default values with new ones for both input channels
- set Coincidence, Trigger mode, Gate and Baseline parameters.

The Digitizer Window is composed of two tabs: “General” and “Channels”.

General tab

The “General” tab contains the “Unit ID and Characteristics” and “Last Digitizer Error received” frames:

- *Handle Number*: once the device is opened, the function returns a handle that becomes the unique identifier of that device; any access operation to the device will take place according to its handle.
- *ROC & AMC Firmware release*: these fields contain the current firmware release running on the mainboard (i.e. on the ROC FPGA) and on the mezzanine (i.e. on the AMC PFGA). Moreover, a message box related to firmware compatibility.
- *Serial and Model Number, PCB revision*
- *Last DGTZ Error received*: any error given back by the CAEN Digitizer library which the program stands on, is reported in the field code.

The DGTZ library return codes are summarized in the Tab. 6.2.


Error code	Value	Meaning
CAEN_DGTZ_Succes	0	Operation completed successfully
CAEN_DGTZ_CommError	-1	Communication error
CAEN_DGTZ_GenericError	-2	Unspecified error
CAEN_DGTZ_InvalidParam	-3	Invalid parameter
CAEN_DGTZ_InvalidLinkType	-4	Invalid Link Type
CAEN_DGTZ_InvalidHandler	-5	Invalid device handler
CAEN_DGTZ_MaxDevicesError	-6	Maximum number of devices exceeded
CAEN_DGTZ_BadBoardType	-7	Operation not allowed on this type of board
CAEN_DGTZ_BadInterruptLev	-8	The interrupt level is not allowed
CAEN_DGTZ_BadEventNumber	-9	The event number is bad
CAEN_DGTZ_ReadDeviceRegisterFail	-10	Unable to read the registry
CAEN_DGTZ_WriteDeviceRegisterFail	-11	Unable to write into the registry
CAEN_DGTZ_InvalidChannelNumber	-13	The Channel is busy
CAEN_DGTZ_ChannelBusy	-14	The channel number is invalid
CAEN_DGTZ_FPIOModelInvalid	-15	Invalid FPIO Mode
CAEN_DGTZ_WrongAcqMode	-16	Wrong acquisition mode
CAEN_DGTZ_FunctionNotAllowed	-17	This function is not allowed for this module
CAEN_DGTZ_Timeout	-18	Communication Timeout
CAEN_DGTZ_InvalidBuffer	-19	The buffer is invalid
CAEN_DGTZ_EventNotFound	-20	The event is not found
CAEN_DGTZ_InvalidEvent	-21	The event is invalid
CAEN_DGTZ_OutOfMemory	-22	Out of memory
CAEN_DGTZ_CalibrationError	-23	Unable to calibrate the board
CAEN_DGTZ_DigitizerNotFound	-24	Unable to open the digitizer
CAEN_DGTZ_DigitizerAlreadyOpen	-25	The Digitizer is already open
CAEN_DGTZ_DigitizerNotReady	-26	The Digitizer is not ready to operate
CAEN_DGTZ_InterruptNotConfigured	-27	The Digitizer has not the IRQ configured
CAEN_DGTZ_DigitizerMemoryCorrupted	-28	The digitizer flash memory is corrupted
CAEN_DGTZ_DPPFirmwareNotSupported	-29	The digitizer DPP firmware is not supported in this lib version
CAEN_DGTZ_InvalidLicense	-30	Invalid Firmware License
CAEN_DGTZ_InvalidDigitizerStatus	-31	The digitizer is found in a corrupted status
CAEN_DGTZ_UnsupportedTrace	-32	The given trace is not supported by the digitizer
CAEN_DGTZ_InvalidProbe	-33	The given probe is not supported for the given digitizer's trace
CAEN_DGTZ_UnsupportedBaseAddress	-34	The Base Address is not supported, as in case of DT and NIM devices
CAEN_DGTZ_NotYetImplemented	-99	The function is not yet implemented

Tab. 6.2: Digitizer library return codes.

Channels tab

The “Channels” tab consists of five sections: *Channel Setting*, *Coincidence*, *Trigger*, *Gate*, and *Baseline*.

◆ Channel Setting Section



The *Channel Setting* section contains:


the switches to enable/disable the two channels of the digitizer.

“*Input DC Offset*” is a percentage shift of the input range scale ($=2 V_{pp}$), allowing the dynamic range to be shifted from -2.0/0 V up to 0/2.0 V. -50% is its minimum value and it corresponds to -2.0/0 V dynamic range. 0% corresponds to a -1.0/+1.0 V dynamic range, and +50% corresponds to 0/2.0 V dynamic range.

“*Trigger Threshold*” is related to the settings of the *Trigger* section of this software panel, and it is available only when the internal trigger mode is selected.

The internal trigger mode uses a CR-RC digital filtering algorithm. After digitalization, the DPP applies the digital filter to the raw input pulse to create a shaped bipolar pulse (called *DELTA*). The trigger and internal gate are generated as soon as the *DELTA* signal is greater than a programmable digital threshold, which is the “trigger threshold”.

◆ Coincidence Section




The *Coincidence* section allows the user to select the coincidence mode if both the channels are switched on.

“*Coincidence Time*” represents the width of the discriminator signal of each channel. Two signals are in coincidence if all of them exceed their own threshold during this time width.

“*GP Output*” allows the user to choose the signal output on the “GPO” of the digitizer front panel between: *Coincidence*, *Gate* and *Discrimination*.

◆ Trigger section




The *Trigger* section allows the user to select external or internal trigger mode. If “*External Mode*” is selected, the digitizer waits for a trigger signal on the “TRG IN” front panel connector. If the “*Internal Mode*” is selected, the digitizer is able to self-detect the signals, according to the trigger parameters.

The purpose of the digital filter is to improve the signal-to-noise ratio by attenuating the low frequencies, (using a numerical differentiator filter) and to smooth out the high frequency noise (using a smoothing function). This filter averages a certain number of samples within a moving window.

“*Mean*” represents the number of sampling used by the average window; the selectable values are 1, 2, 4, 8, 16 and 32.


Rise Time is the rise time of the input signal, used in the calculation of the signal *DELTA*.

◆ General Purpose I/O section



General Purpose I/O section allows the user to set input and output levels as NIM or TTL.

◆ Baseline Section



The *Baseline* section contains the controls for the baseline evaluation.

The “*Mean*” parameter is the number of samples for the average calculation of the baseline. The value 0 disables the baseline restoration. The “*Threshold*” represents the value on *DELTA*, over that the baseline calculation is frozen, and “*no flat*” is the veto for the calculation of baseline.

◆ *The Gate Section*

The *Gate* section consists of three parameters.

- “*Gate*” represents the width of the gate signals.
- “*Pre-Gate*” is the time between the gate generation and the trigger leading edge.
- “*Holdoff*” is a veto for the generation of other gates.

◆ *Baseline Section*

The *Baseline* section contains the controls for the baseline evaluation.

The “*Mean*” parameter is the number of samples for the average calculation of the baseline. The value 0 disables the baseline restoration. The “*Threshold*” represents the value on DELTA, over that the baseline calculation is frozen, and “*no flat*” is the veto for the calculation of baseline.



Important Note: The “*Update*” button must be selected every time to communicate and apply the selected parameters to the DGTZ. The lack of this operation leaves the default settings unchanged. The “*Auto Update*” button executes the updating process automatically.

7 Basic Measurements

This manual section is dedicated to the simple and practical use to perform the first basic measurements by using the Educational Gamma Kit.


7.1 Can you see γ spectrum?

7.1.1 First Kit Configuration

- Required elements: PSAU + Digitizer + Spectrometer + [Oscilloscope].
- Cabling instructions:
 - The main units of the kit must be power on.




SP5600 – Power Supply and Amplification Unit




DT5720A – Desktop Digitizer


- SP5600 (PSAU) and DT5720A (DGTZ) shall be connected to the PC via USB cables.




- In order to avoid saturation during the spectroscopy measurements, the output signal of the spectrometer is divided using the A315 splitter.



- Power Cable of the Spectrometer shall be connected to the PSAU channel (for example channel 0).




- One splitter branch shall be connected to PSAU channel 0 in order to be amplified by the module. The other splitter branch shall be connected to channel 0 input on the front panel of the digitizer and will be digitized.



- The Digital Output of the used PSAU channel shall be connected to the TRG IN on the front panel of the Digitizer.




Important Note: How can you connect the *Power Cable* of the Spectrometer board to SP5600?



Normally, the board power cable is already connected to the bottom part of the spectrometer. If the connection doesn't exist, take careful to attach the marked extremity of the cable connector to the "J2" side of the board connector as shown in picture.

To connect the board power cable to the SP5600 channel, put the marked cable side close to the channel label.

- Getting the system alive:
 - Power on the kit elements
 - Run the program by clicking the HERA icon and wait for the hardware connection. The software recognises the hardware automatically and starts the connection. Keep attention to the two connection indicators status.
 - Select the Hardware Management access.



- Activate the connected PSAU channel and Digitizer channel

Once the system is running, the first action to take is properly biasing the detector and setting the gain to the PSAU amplifier.

As far as the optimal sensor bias, it is suggested to stick to the value reported on the sensor ID card, which may be set in the Bias & Gain tab of the PSAU panel (Fig. 7.1).




Fig. 7.1: Power Supply & Amplification Unit Interactive Panel – Gamma spectroscopy settings.

At the same time, the amplification factor can be set at high values, e.g. 40 dB, due to the PSAU digital output works only as trigger.

Moreover, for the sake of clarity, the feedback system for the SiPM gain stabilization against temperature variations can be disabled.

Putting on the oscilloscope the spectrometer signal output before (Fig. 7.2) or after the splitting, it's hard to distinguish the different traces corresponding to avalanches in the cells triggered by the thermal generation of the charge carriers or by the photons associated to the avalanche development (optical cross-talk). This fact is due to a combination of factors: no signal amplification and big cells number.




Fig. 7.2: Scope trace of the spectrometer.



Important Note: In order to see multiple traces, it is better to put the spectrometer signal into the PSAU, set the gain and use the analog output to see the bands on the oscilloscope.

The intensity of the light measured by the detecting system (SiPM + scintillating crystal) is proportional to the γ source activity. However, this information is affected and biased by stochastic effects characteristic of the sensor and occurring within the time window: spurious avalanches due to thermally generated carriers (a.k.a. Dark Count Rate (DCR)).

Before to acquire a γ source spectrum, take care to measure this entry-level parameter of the SiPM. It is a standard procedure to quantify the DCR as the counting frequency with a threshold corresponding to 0.5 x single photo-electron (p.e.) peak ($DCR_{0.5}$). The DCR of the mini spectrometer SiPM is 2÷6 MHz.

In order to avoid the system be blind to the radioactive source, due to this high DCR, a proper cut-off threshold has to be selected. The DCR vs discriminator threshold can be precisely measured by using the “PSAU Staircase” tab in the SiPM Kit Control Software.

The Fig. 7.3 shows the DCR staircase: once the SiPM is biased, and a convenient gain is chosen, through the “Frequency vs Threshold” tab in the right part of the GUI, the threshold scan of the signals frequency can be performed.

Through two staircases it is possible to choose a right cut-off threshold: one acquired with the source and one acquired without it.






Fig. 7.3: Dark Counts Rate frequency versus Discriminator threshold.

7.1.2 Second Kit Configuration


- Required elements: PSAU + Digitizer + Spectrometer + Splitter + Radioactive Source + [Oscilloscope]
- Mechanical instructions:
 - In the educational Gamma kit, there are one mini spectrometer and three scintillating crystals.
 - Before starting any application with radioactive sources, remove the screws and open the spectrometer.



- Choose the scintillating crystal for the application (It doesn't matter if you use the hands to manage the crystal, but a small pliers could help you).
- Take care to spread homogeneously optical grease on the open face of the scintillating crystal (the only one not polished).
- Insert the crystal greased side in the spectrometer, put on top the radioactive source, close it with the screws and start the application.




- Cabling instructions:
 - The PSAU and the Digitizer shall be connected to the PC via the USB.



- The Power Cable of the Spectrometer shall be connected to the PSAU in ch0 position (for example).
- In order to avoid saturation, the output signal of the spectrometer is divided using the A315 splitter: one branch shall be connected to channel 0 input on the front panel of the digitizer and will be digitized. The other branch will be amplified by the SP5600 module (ch0) and


the digital output shall be used as DT5720A “trigger IN” for the integration signal by the on-board leading edge discriminator.



- Getting the system alive:
 - Power on the kit elements




SP5600 – Power Supply and Amplification Unit



DT5720A – Desktop Digitizer


- Run the program by clicking the HERA icon and wait for the hardware connection. The software recognises the hardware automatically and starts the connection. Keep attention to the two connection indicators status.
- Select the Hardware Management access.



- Activate the connected PSAU channel and Digitizer channel.

7.1.3 Correctly bias the SiPM and set the gain

- The first action to take is biasing the SiPM and setting the gain: set the bias according to the detector datasheet (in this example bias = 55.4 V, as before). For the gain, since the SiPM signal in the spectrometer application will be split, a high value can be used (e.g. 40 dB, as before).



- The discriminator threshold shall be defined looking at the spectrum and evaluating the dark count rate. Once fixed threshold if the radioactive source is properly positioned the spectrum can be recorded.

In Fig. 7.4 a typical scope trace of the spectrometer with ^{22}Na source.




Fig. 7.4: Scope trace of the spectrometer with a ^{22}Na source.

- Software: in the Digitizer Interactive Panel (Fig. 7.5).
- Select EXTERNAL trigger mode.
- Select the active channel (for example channel 0).
- Accept default values for the Gate and Baseline sub-panels.




Fig. 7.5: The DIGITIZER control panel.

7.1.4 External and internal trigger

The kit allows the user to trigger the digitizer for the acquisition of a source spectrum in two different ways: internally to the digitizer, selecting an appropriate threshold value (on delta signal) or externally to the digitizer, setting the discriminator threshold in the PSAU, through the PSAU discriminator tab, and using the correspondent NIM/TTL signal produced by the PSAU as a trigger signal for the digitizer itself.

If the digitizer external trigger mode is selected, a threshold value on the PSAU discriminator tab shall be selected, as suggested in this application. With a wrong choice of this value, the PSAU triggers the signal on its noise. Increasing the threshold value, the particle signal shall be clear, becoming similar to the wave of Fig. 7.6.




Fig. 7.6: Wave for a right PSAU trigger threshold.

How select the “right” external threshold value?

Running a staircase (in the Frequency vs Threshold tab) from -5 mV to -110 mV, the Dark Count Rate decreases drastically from ~ 0.6 MHz to a value around 1 kHz (as shown in Fig. 7.7). This example is performed with a ^{22}Na source of an activity of 16 kBq, so, in order to remove all the dark count rate due to the nature of the sensor, a threshold bigger than 60 ÷ 70 mV can be selected. In fact, the frequency measured with source remains constant over this threshold value, while the Dark Count Rate drops to zero.




Fig. 7.7: Frequency vs Threshold for SiPM 6x6 mm² of spectrometer, in a run without crystal, with crystal and with Crystal and radioactive source.

If digitizer internal trigger mode is selected, set the digitizer trigger parameters and the channel threshold in the Digitizer Interactive Panel: before the calculation of delta (see the “Digital Pulse Processing for SiPM kit” document [RD4]), the input signal is filtered in order to reduce the high frequency noise, using a low pass filter that averages a certain number of samples within a moving window. The Fig. 7.8 shows the trigger parameters the user can set:

- *mean[#]* represents the number of double sampling periods used by the average window;
- *rise time[ns]* is the rise time of the input signal, used in the calculation on delta.




Fig. 7.8: The Trigger digitizer parameters.

In this example the external trigger mode is suggested. Please don't forget the connection of the LEMO cable between the digital output of the PSAU and the trigger input of the digitizer.

By now, the system is ready for digitizing the signal but, rather than doing it in a blind way, it is worth taking a guided tour of the system features, going to the Visualization panels of the HERA Main GUI and switching ON the Waveform Tab.



Important Note: Putting 0 as number of samples used to calculate the mean, the baseline restoration is disabled, therefore the histogram doesn't start from zero [ADC channels].




Fig. 7.9: The Waveform tab of the HERA Main GUI.

The Gate actually defines the integration time. For the sake of clarity in the display, every signal can be offset and magnified, enabled or disabled.

Fig. 7.9 is showing the WAVE panel for the SiPM coupled to LYSO crystal and ^{22}Na source, for optimal tuning of the pre-gate and the gate width, depending on the time development of the signal. The signal time development is dominated not only by the sensor response, but also by the decay time of the scintillating crystal. Moreover, the waveform is positive because the digitized signal comes directly from the spectrometer (after splitting). This configuration allows to avoid saturation effects and overshoot due to AC coupling of the SP5600.

As long as the Gate is properly defined, the system is ready to record the spectrum, displayed in the Histogram tab. As exemplary illustration, Fig. 7.10 shows a spectrum of ^{22}Na . It was obtained with the mini-spectrometer with the Hamamatsu S13360- 6050CS SiPM coupled with a 6 x 6 x 15 mm² LYSO crystal.




Fig. 7.10: ^{22}Na spectrum.

The threshold is fixed in order to cut the system noise peak, so the spectrum shows the following peaks:

- The backscatter at the low energy and the Compton distribution for the first peak;
- First γ peak @ 511KeV: Electron-positron annihilation;
- Compton distribution for the second peak;
- Second γ peak @ 1275KeV: Nuclear Energy Transition.

The energy resolution of the peaks represents the spectroscopic capability of the detector at different energy values.

8 Educational Experiments

The Educational Gamma kit allows to perform experiments that have to do with radioactive γ decays.

In 1895 the radioactivity was discovered by H. Becquerel and in 1903 the Nobel Prize in Physics was assigned to Curie spouses for their studies on these phenomena. Radioactivity is around us and getting to know it experimentally is essential for physics students.

When an unstable nucleus decays in a cascade leading to a stable nuclide, it emits α or β or γ quanta or a combination of them. Gamma rays are high energy photons and the spectroscopy of the emitted γ rays is instrumental for understanding the mechanism of the interaction with matter, the fundamentals about detection and the underlying nuclear physics. Moreover, it is relevant in basic and applied fields of science and technology, from nuclear to medical physics, from archaeometry to homeland security.

This section represents an overview of the experiments proposed by CAEN using the Educational kit of your choice. Each experiment has its own identification code (reference ID). For each ID, a step by step guide that includes a detailed description to perform the data analysis of the physical process is available on the CAEN Educational web page. The experiments address the essence of the phenomenon as well as exemplary illustrations of their use in medical imaging and industry, complemented by basic and advanced statistical exercises.

The experiments proposed by CAEN in Gamma Spectroscopy field are listed in Tab. 8.1.

Section	Subsection	Reference ID	Experiment
Nuclear Physics and Radioactivity	Gamma Spectroscopy	6111	Detecting γ -Radiation
		6112	Poisson and Gaussian Distributions
		6113	γ Spectrum and Energy Resolution
		6114	System Calibration: Linearity and Resolution
		6115	Scintillator Crystals Comparison: Light Yield and Decay Time
		6116	γ -Radiation Absorption
		6117	Photonuclear cross-section/Compton Scattering cross-section
		6118	Study of the ^{137}Cs spectrum: the backscatter peak and X rays

Tab. 8.1: Physics Experiments performed via the Educational Kit – Premium version.

8.1 Detecting γ -radiation (SG6111A)

Purpose of the experiment:

Gamma radioactivity detection by using a system composed of a scintillating crystal coupled to a photon detector.

Fundamentals:

Gamma rays interact with matter by three processes: Compton Scattering, Photoelectric Effect and Pair Production (whenever the energy exceeds the 1.022 MeV threshold corresponding to the e^+e^- rest mass). The cross section of each process depends on the energy of the gamma ray.

The Compton Effect is the inelastic scattering between the incoming photon and an atomic electron. In the Photoelectric Effect, the incident gamma ray transfers all of its energy to a bound electron which acquires a kinetic energy equal to the incoming gamma energy decreased by the binding energy.


These processes convert, totally or partially, the gamma ray energy into kinetic energy of electrons (or positrons, in case of pair production). The interaction of the charged particles with the atomic and molecular systems of the medium results in excited states whose decay, possibly mediated, leads to light in the visible or UV region, eventually detected by the light sensor. A wide range of scintillator products is available today, differing for the light yield, the material properties, the time characteristics of the scintillation light and, last but not least, cost. The choice of the scintillator is essentially dependent on the specific targeted application.

Requirements:

Gamma Radioactive Source.

Carrying out the experiment:


The selected scintillator crystal shall be coupled to the SiPM in the SP5607, through a thin layer of index matching grease to maximize the light collection. In order to avoid saturation, the output of the SiPM is divided using the A315 splitter: one branch is connected to the DT5720A and will be digitized. The other branch will be amplified by the SP5600 module, generating the trigger for the integration signal by the on-board leading edge discriminator or simply counting the pulses induced by the detected gamma ray.



Experimental setup block diagram.

Results:

The student may get acquainted with the presence of radioactivity with a simple preliminary measurement, namely comparing the counting frequency as a function of the discriminator threshold with/without the source. Presuming the source, essentially in contact to the crystal, to be point like with respect to the crystal surface, and assuming its activity is known, the student may estimate for every threshold value the detection efficiency and the signal over noise ratio, building up an efficiency-purity plot. Exemplary results obtained with a ^{137}Cs Source are shown. Moving away the source from the crystal, the law governing the variation of the flux can also be investigated.



Sensor output frequency as a function of the threshold in mV, with and without ^{137}Cs source.

8.2 Poisson and Gaussian Distribution (SG6112A)

Purpose of the experiment:

Study the statistical distribution of the counting rates of a gamma radioactive source. Comparison of the data to the Poisson distribution, turning into a Gaussian as the mean number of counts grows. The study can be performed both experimentally, with the SiPM kit or simulating it with the emulation kit.

Fundamentals:

The number of radioactive particles detected over a time Δt is expected to follow a Poisson distribution with mean value μ . It means that for a given radioactive source, the probability that n decays will occur over a given time period Δt is given by:

$$P_{\mu}(n) = \frac{\mu^n}{n!} e^{-\mu}$$

Where μ is proportional to the sample size and to the time Δt and inversely proportional to the half-life $T_{1/2}$ of the unstable nucleus. As long as μ grows, the probability $P_{\mu}(n)$ is well approximated by a Gaussian distribution:

$$P(n) = \frac{1}{\sqrt{2\pi}\sigma} e^{-\frac{(n-\mu)^2}{2\sigma^2}}$$


Where $\sigma = \sqrt{\mu}$ is the standard deviation.

Requirements:

Gamma Radioactive Source

Carrying out the experiment:


The selected scintillator crystal shall be coupled to the SiPM in the SP5607, through a thin layer of index matching grease to maximize the light collection. In order to avoid saturation, the output of the SiPM is divided using the A315 splitter: one branch is connected to the DT5720A and will be digitized. The other branch will be amplified by the SP5600 module, generating the trigger for the integration signal by the on-board leading edge discriminator or simply counting the pulses induced by the detected gamma ray. The discriminator threshold shall be defined looking at the spectrum and evaluating the dark count rate. Once this is properly set, the counting experiment shall be performed.




Experimental setup block diagram.

Results:

Changing the counting window and/or the activity of the source or the threshold, the number of counts changes, with a probability density function moving from a Poissonian to a Gaussian shape. The student may play with the data, fitting them and comparing the expectations to the measurement.



Poissonian distribution.



Gaussian distribution.

8.3 Energy Resolution (SG6113A)

Purpose of the experiment:

The analysis of the spectrum of the deposited energy by a γ ray in a detector discloses the essence of the interaction of high energy photons with matter and allows to learn by doing the detector related effects.

Fundamentals:


For γ -energy less than 2 MeV, the interaction with matter is dominated by Compton scattering and Photo-absorption. The analysis of the Compton continuum of the deposited energy and of the photo-peak conveys information on the characteristics of the decaying isotope as well as the effects due to the system noise, the detected photon statistics, the stochastic terms in the detector and the intrinsic resolution of the scintillator. The experiment presumes to use ^{137}Cs with its decays detected by a CsI crystal coupled to a Silicon Photomultiplier. The ^{137}Cs source is particularly interesting due to its low energy X ray line at 33 keV and the high energy gamma emission at 662 keV. The former is relevant to optimize the lower detection limit of the system; the latter is a standard to evaluate the energy resolution. The use of the 2 lines and the analysis of the Compton spectrum characteristics allow to perform a rough measurement of the linearity with a single isotope.

Requirements:

Gamma Radioactive Source.

Carrying out the experiment:


The CsI scintillator crystal shall be coupled to the SiPM in the SP5607, through a thin layer of index matching grease to maximize the light collection. In order to avoid saturation, the output of the SiPM is divided using the A315 splitter: one branch is connected to the DT5720A and will be digitized. The other branch will be amplified by the SP5600 module, generating the trigger for the integration signal by the on-board leading edge discriminator. The discriminator threshold shall be defined looking at the spectrum and evaluating the dark count rate. Once this is properly set and the radioactive source is properly positioned, the spectrum can be recorded.



Experimental setup block diagram.

Results:

The figure shows a typical gamma spectrum, recorded with a very low energy threshold. The left over from the system noise is clearly visible, as well as the low energy line at 33 keV and the photo-peak. For this specific spectrum, the energy resolution on the 662 keV peak corresponds to



The ^{137}Cs spectrum; activity of the source: 180 kBq.

$$\text{Energy Resolution} = \frac{\text{FWHM}_{\text{peak}}}{\mu_{\text{peak}}} * 100 \sim 10\%$$

$\text{FWHM}_{\text{peak}}$ = full width at half maximum of the peak

μ_{peak} = channel number of the peak centroid.

8.4 System Calibration: Linearity and Resolution (SG6114A)

Purpose of the experiment:

Recording and comparing the γ energy spectra of several radioactive sources is the main goal of the experiment. The photo-peaks are used to calibrate the response of the system and to measure the energy resolution.

Fundamentals:

Linearity and energy resolution are the main figures of merit of a spectrometric system. In the proposed experiment, based on a scintillating crystal coupled to a Silicon Photomultipliers, deviations in the linearity may be due to the sensor or the front-end electronics saturation. The student is guided through the analysis of the response curve using a series of isotopes up to the MeV energy by a ^{60}Co source and to disentangle the different effects.

At the same time, the energy resolution of the system is measured by the width of the photo-peaks and the results compared to what is expected by the fluctuations in the number of detected scintillation photons, the system noise, the sensor stochastic effects, the intrinsic resolution of the scintillator.


This is following an initial activity on the optimization of the operating parameters by an analysis of the photo-peak position and the resolution for a single isotope.

Requirements:

Gamma Radioactive Sources.


Carrying out the experiment:

The scintillator crystal shall be coupled to the SiPM in the SP5607, through a thin layer of index matching grease to maximize the light collection. In order to avoid saturation, the output of the SiPM is divided using the A315 splitter: one branch is connected to the DT5720A and will be digitized. The other branch will be amplified by the SP5600 module, generating the trigger for the integration signal by the on-board leading edge discriminator. The discriminator threshold shall be defined looking at the spectrum and evaluating the dark count rate. Once this is set and the radioactive source is properly positioned, the spectrum can be recorded.




Results:

By fitting the photo-peaks with a Gaussian curve, the system linearity as a function of energy is verified. The peak widths are determining the energy resolutions. At more advanced level, the interpretation of the results accounting for the system properties may be performed.



Energy calibration.



Energy dependence of the system resolution.

8.5 A comparison of different scintillating crystals: Light Yield , Decay Time and resolution (SG6115A)

Purpose of the experiment:

Compare the basic characteristics of different scintillating crystals, namely the light yield and the decay time of the scintillation light. Verify the effect on the energy resolution.

Fundamentals:

Scintillating materials have different characteristics related to the light yield and the characteristics time of the emission. The CAEN spectrometer is provided with three different crystals: BGO (Bismuth Germanate), LYSO(Ce) (Cerium-doped Lutetium Yttrium Orthosilicate), CsI(Tl) (Thallium-doped Cesium Iodide). All of them have the same volume ($6 \times 6 \times 15 \text{ mm}^3$), are polished on all sides and coated with a white epoxy on 5 faces. One $6 \times 6 \text{ mm}^2$ face is open in order to be coupled with the Silicon Photomultiplier. The main characteristics of the crystals are summarized in the following table:

	BGO	LYSO(Ce)	CsI(Tl)
Density (g/cm^3)	7.13	7.4	4.51
Decay Time (ns)	300	40	1000
Light Yield (ph./MeV)	8200	27000	52000
Peak emission (nm)	480	420	560
Radiation length (cm)	1.13	1.14	1.85
Reflective index	2.15	1.82	1.78


The light yield is having an impact on the energy resolution. This is also affected by the decay time, constraining the integration time and implying a different effect of the sensor stochastic effects (dark counts and afterpulses).

Requirements:

Gamma Radioactive Source.

Carrying out the experiment:


The scintillator crystal shall be coupled to the SiPM in the SP5607, through a thin layer of index matching grease to maximize the light collection. In order to avoid saturation, the output of the SiPM is divided using the A315 splitter: one branch is connected to the DT5720A and will be digitized. The other branch will be amplified by the SP5600 module, generating the trigger for the integration signal by the on-board leading edge discriminator. The discriminator threshold shall be defined looking at the spectrum and evaluating the dark count rate. Once this is set and the radioactive source is properly positioned, the spectrum can be recorded. The procedure shall be repeated for every crystal.



Experimental setup block diagram.

Results:

The crystal characteristics are investigated recording a source spectrum (for example ^{137}Cs) with the three different crystals, optimizing the integration time as a function of the scintillation decay time. According to table, the Light Yield of the three crystal is very different. LYSO(Ce) has a light yield three times greater than the BGO, and CsI(Tl) light yield is twice than LYSO(Ce). The analysis of the signal waveform or the trend of the charge vs integration time leads to the measurement of the time characteristics of the scintillator.



^{137}Cs energy spectra. Blue spectrum corresponds to the acquisition through LYSO crystal, the red and black ones respectively with BGO and CsI crystals.

8.6 γ -Radiation Absorption (SG6116A)

Purpose of the experiment:

The main goal of the experiment is the measurement of the γ radiation attenuation coefficient for different materials and different energies.

Fundamentals:

The attenuation of a γ radiation flux passing through matter is described by the exponential law:

$$I(x) = I_0 * e^{-\mu x}$$


where I_0 is the incident photon flux and $I(x)$ measures the flux of γ rays emerging from a layer x of material without having interacted. The coefficient μ depends on the material properties (atomic number, density) and on the energy of the impinging photon. The student is guided towards the development of complementary measurement techniques based on counting and on the analysis of the spectrum, performing the experiment for different materials (including PMMA, a water equivalent solid state organic material used in medical dosimetry).

Requirements:

Gamma Radioactive Source.

Carrying out the experiment:


The scintillator crystal shall be coupled to the SiPM in the SP5607, through a thin layer of index matching grease to maximize the light collection. In order to avoid saturation, the output of the SiPM is divided using the A315 splitter: one branch is connected to the DT5720A and will be digitized. The other branch will be amplified by the SP5600 module, generating the trigger for the integration signal by the on-board leading edge discriminator. The discriminator threshold shall be defined looking at the spectrum and evaluating the dark count rate. Once this is set the SP5609 absorption tool shall be mounted. The experiment can be performed for every absorber thickness in counting mode and analysing the spectrum, measuring the events in the photo-peak for a constant pre-defined time interval.




Experimental setup block diagram.

Results:


Exemplary results are shown below, reporting the variation of the events in the photopeak for different absorber thickness, a plot verifying the exponential absorption law and the dependence of the absorption coefficient on the energy.



Gamma spectra acquired with different absorber thicknesses.



Linear dependence of logarithmic intensity of gamma rays as a function of penetration thickness.



Gamma attenuation coefficient as a function of energy.

8.7 Photonuclear cross-section/Compton Scattering cross-section (SG6117A)


Purpose of the experiment:

Determination of the ratio of the effective cross-sections due to Compton and Photoelectric effects as a function of photons energy.


Fundamentals:

In the energy range up to 2 MeV, gamma rays interact with matter by two processes:

- Photoelectric Effect, dominant at energy less than 100 KeV. In this process the photon energy is completely transferred to atomic electron bounded



- Compton Scattering, linked to the elastic collision between electrons and photons and relevant at 1 MeV energy level




The predominant mode of interaction depends on the energy of the incident photons and the atomic number of the material with which they are interacting. From the acquired γ -spectrum, it is possible to estimate the fraction of events due to Compton scattering and those caused by the photoelectric. The ratio of the event fractions is used to determine the ratio of the two effective cross-sections that depends on the detector size.

Requirements:

Gamma Radioactive Source.

Carrying out the experiment:


Spread the optical grease on the open face of the scintillating crystal, insert this crystal side in the SP5607 spectrometer. Connect the power cable to the SP5600 module and connect the other cable of the spectrometer to the splitter A315. Connect the two split outputs to SP5600 channel 0 and DT5720A channel 0. Use the SP5600 digital output as DT5720A "trigger IN". Use the default software values or optimize the parameters to choose the discriminator cut-off threshold in mV. Switch off the power supply, open the spectrometer and insert the radioactive gamma source to acquire the first spectrum. After that, switch off the power supply, open the spectrometer, change the radioactive gamma source and repeat the measurement.




Experimental setup block diagram.

Results:

By using several radioactive sources, the energy dependence of the ratio between the cross-sections of two phenomena can be examined, by verifying that the Photoelectric Effect cross section decreases with increasing energy compared to the Compton Scattering cross section for the used detector size.



Spectra of radioactive sources used to estimate the ratio of Photonuclear and scattering Compton cross sections.



Behaviour of the ratio between Photo-Peak and Compton contributions as a function of energy.

8.8 Study of the ^{137}Cs spectrum: the backscatter peak and X rays (SG6118A)

Purpose of the experiment:

Study the characteristics of the ^{137}Cs spectrum, with special relevance given to the low energy spectrum. The student can learn effects related to the experimental observation of a gamma decay and have basic information about the experimental setup used in gamma spectroscopy. Estimate the energy of the backscatter peak and of the K_α line.

Fundamentals:

The Compton effect is linked with experimental issues, since it is caused by the interaction of photons with the electrons instrument that measure the gamma radiation. In a real detector setup, some photons can and will undergo one or potentially more Compton scattering processes (e.g. in the housing material of the radioactive source, in shielding material or material otherwise surrounding the experiment) before entering the detector material. This leads to a peak structure, the so-called backscatter peak.


The basic principle for the backscatter peak formation is the following: gamma-ray sources emit photons isotropically, some photons will undergo a Compton scattering process with a scattering angle close to 180° and some of these photons will subsequently be detected by the detector. The result is an excess of counts in the Compton part of the spectrum, the so-called backscatter peak. This peak has an energy approximately equal to the photopeak energy minus the Compton edge energy.

Requirements:

Gamma Radioactive Source: ^{137}Cs

Carrying out the experiment:


The scintillating crystal shall be coupled to the SiPM in the SP5607, through a thin layer of index matching grease to maximize the light collection. In order to avoid saturation, the output of the SiPM is divided using the A315 splitter: one branch is connected to the DT5720A and will be digitized. The other branch will be amplified by the SP5600 module, generating the trigger for the integration signal by the on-board leading edge discriminator. The discriminator threshold shall be defined looking at the spectrum and evaluating the dark count rate. Once this is properly set and the radioactive source is properly positioned, the spectrum can be recorded.



Experimental setup block diagram.

Results:

The user can calibrate the system by using the spectrum itself. The backscatter peak and the K_α line can be identified. After calibrating the spectrum, it is possible to estimate the energy of the two peaks and compare them with theoretical predictions.



Plot of the ^{137}Cs spectrum acquired by HERA software. The backscatter peak and the K_α line are indicated with the red arrows.

9 Appendix

Data Storage

The HERA system allows the user to save data in several ways:

- Save data during the run (streaming mode).
- Save the data for an offline run (snapshot mode).
- Save an image.
- Export data in Excel.

The generation of the files due to the experimental activity deserves a dedicated discussion and it is referred to the step-by-step guides of each experiment.

All directories and files are generated under the following path: *C:\ProgramData\HERA\UserName* .

Where *UserName* is the name of the logged user. If no specific username is chosen, the default name used is "Guest".

The list of the directories created by the HERA system during Hardware Management usage is the following:

- PSAU Temperatures
- Waveform
- Histograms
- Charge-Time
- TDMS

HERA generates several file formats: ASCII (.txt), binary (.dat), and another special file format of Labview (.TDMS) format.

National Instruments defined a new flexible technical data management (TDM) data model, which is accessible through LabVIEW, LabWindows™/CVI™, Measurement Studio, SignalExpress, and DIAdem.

The TDMS file format saves both the raw data and the metadata in the binary format in one file with the .TDMS extension.

When creating or opening a .TDMS file, HERA automatically creates a .TDMS_index file, used to speed up random access to the .TDMS file.

The .TDMS files can be open via a simple Add-In for Microsoft Excel (<https://www.ni.com/example/27944/en/>) or by using the "Convert data files (.TDMS)" in the File Menu.

PSAU Temperature

Regardless of the type of acquisition and the tab in use, a new file (ASCII format) is stored at every change of date (Log file type).

.txt Structure			
Typical Filename	PSAU_Temperature_date(mmm-dd-yy).txt	Ex.: PSAU Temperature_ Jan-15-21.txt	
General	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).		
File Properties recorder	Dataset Conditions: absent		
Organisation:		Name	Type
Structure	5 columns		
Separation character	TAB		
Channels	Channel(s): 5	Date	Date (O.S. format)
		Time	24 h format
		Board	Decimal
		Temperature	float
		Ch.0 Temperature	Decimal float
		Ch.1 Temperature	Decimal float
Data Format	Decimal separator: point (.)		
Data Type Recorded	Single points of measures		
Length	Depends on the running time of the Main Program		

Waveform tab

In addition to the waveforms export in a Bitmap Image to the Clipboard or "Excel" numerical data via the "Export" button, the Waveforms can be saved in both modes, streaming, and snapshot. The generated files are in .TDMS format (see Tab. 9.1).

Folder	Streaming Mode	Snapshot Mode
Waveform	.TDMS <i>streaming</i>	--
TDMS	--	.TDMS

Tab. 9.1: Waveforms saving scheme.

- Streaming Mode

TDMS Structure		
Typical Filename	<i>Activity Acronym_Wave(Time or events xx)_date(mm-dd-yy)-T-time(hhmm).TDMS</i>	Ex.: <i>HRDW_Wave(#Evn 1000)_01-08-21-T-1155.TDMS</i>
General	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name
Existing Groups	Group(s): 1	Analog Waveforms
Channels	Channel(s): 2	Trace Ch. 1
		Trace Ch. 0
Channel Range	0..4095	
Data Format	DT_Float (floating point double precision, 64 bits)	
Data Type	Array	
Length	Depends on the acquisition time or # of triggers	

- Snapshot Mode

TDMS Structure		
Typical Filename	<i>Activity Acronym_Wave.TDMS</i>	Ex.: <i>SiPM_Wave.TDMS</i>
General	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name
Existing Groups	Group(s): 2	Analog Waveforms
		Digital Waveforms
Channels	"Analog Waveforms" Channels: 2	Trace Ch. 1
		Trace Ch. 0
	"Digital Waveforms" Channels: 5	Virtual Probe
		Gate
		Hold Off
		Over Threshold
		Flat
Channel Range	0..4095	
Data Format	DT_Long (Long 32 bits integer)	
Data Type	Arrays	
Length	Depends on the x scale extension of the Waveform plots originating the file	

Histogram tab

The histograms can be export in a Bitmap Image to the Clipboard and in "Excel" numerical data via the "Export" button. Moreover, as the waveforms saving, the histograms can be saved in streaming and snapshot mode. The generated file formats are summarized in Tab. 9.2.

Folder	Streaming Mode	Snapshot Mode
Histogram	.txt (ASCII)	.txt (ASCII) [Under request]
TDMS	--	.TDMS [Under request]

Tab. 9.2: Histograms saving scheme.

- Streaming Mode

.txt Structure		
Typical Filename	<i>Activity Acronym_Charge_Histo(Time or events xx)_date(mm-dd-yy)—Time (hhmm).txt</i>	Ex.: <i>HRDW_Charge_Histo(Time 10)_02-10-21 Time 1610.txt</i>
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name Type
Structure	2 columns	
Channels	Channel(s): 2	ADC Channel Decimal float
		Counts Integer
Data Format	Decimal separator: point (.)	
Data Type	Array(s)	
Recorded Length	Depends on the number of bins present in the Histogram	

- Snapshot Mode

.txt Structure		
Typical Filename	<i>Activity Acronym_Charge_Histo_date(mm-dd-yy)-Time-time(hhmm).txt</i>	Ex.: <i>HRDW_Charge_Histo_01-25-21 Time 1048.txt</i>
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name Type
Structure	2 columns	
Channels	Channel(s): 2	ADC Channel Decimal float
		Counts Integer
Data Format	Decimal separator: point (.)	
Data Type	Array(s)	
Recorded Length	Depends on the number of bins present in the Histogram	

<i>TDMS Structure</i>		
Typical Filename	Activity Acronym_Charge_Histo.TDMS	Ex.: HRDW_Charge_Histo.TDMS
General	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
File Properties recorder	Dataset Conditions: empty	
Organisation:		Name
Existing Groups	Group(s): 1	Charge Histogram
Channels	Channel(s): 2	X coord Histo(X)
Data Format	DT_Float (floating point double precision, 64 bits) DT_Long (long 32 bits integer)	
Data Type Recorded	Array	
Length	Depends on the number of bins present in the Histograms	

Two Channels tab

The storage of the histograms in the *Two Channels tab* can be occurred in snapshot mode only, as showed in Tab. 9.3.

Folder	Streaming Mode	Snapshot Mode
Histogram	--	.txt (ASCII), 2 separate files
TDMS	--	.TDMS (2 files) [Under request]

Tab. 9.3: Two channels saving scheme.

- Snapshot Mode

<i>.txt Structure</i>		
Typical Filename	Hardware Management (generic)Chx(channel number)_date(mm-dd-yy)-Time-time(hhmm).txt A file is generated for each channel.	Ex.: Hardware Management (generic)Ch1_02-09-21 Time 1238.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name Type
Structure	2 columns	
Channels	Channel(s): 2	Rate bin Decimal float Counts/bin Integer
Data Format	Decimal separator: point (.)	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histograms	

TDMS Structure		
Typical Filename	<i>Activity Acronym_2Ch_Charge_Histo.TDMS</i>	Ex.: <i>HRDW_2Ch_Charge_Histo.TDMS</i>
General	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
File Properties recorder	Dataset Conditions: empty	
Organisation:		Name
Existing Groups	Group(s): 2	Channel 0 Histogram Channel 1 Histogram
Channels	Channel(s): total 4, (2 per group)	X coord Histo(X)
Data Format	DT_Float (floating point double precision, 64 bits) DT_Long (long 32 bits integer)	
Data Type Recorded	Array	
Length	Depends on the number of bins present in the Histograms	

Charge vs Time tab

The Charge vs Time data can be saved in streaming and snapshot mode. The generated file formats are summarized in Tab. 9.4.

Folder	Streaming Mode	Snapshot Mode
Charge-Time	.TDMS streaming	--
TDMS	--	.TDMS (2 files) [Under request]

Tab. 9.4: Charge vs Time saving scheme.

- Streaming Mode

TDMS Structure		
Typical Filename	<i>Activity Acronym_ChargeVSTime(Time or events xx)_date(mm-dd-yy)-T-time(hhmm).TDMS</i>	Ex.: <i>HRDW_ChargeVSTime(Time 5)_02-09-21-T-1245.TDMS</i>
General	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name
Existing Groups	Group(s): 1	ChargeVSTime
Channels	Channel(s): 2	Charge DGTZ-Ch.0 Charge DGTZ-Ch.1
Data Format	DT_Float (floating point double precision, 64 bits)	
Data Type Recorded	Array	
Length	Depends on the acquisition time or # of triggers	

- Snapshot Mode

TDMS Structure		
Typical Filename	<i>Activity Acronym_ChargeVSTime.TDMS</i>	Ex.: <i>HRDW_ChargeVSTime.TDMS</i>
General	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
File Properties recorder	Dataset Conditions: empty	
Organisation:		
Existing Groups		Name
Channels	Group(s): 1	ChargeVSTime
	Channel(s): 2	Trace Ch. 0
		Trace Ch. 1
Channel Range	Full range	
Data Format	DT_Float (floating point double precision, 64 bits)	
Data Type	Array	
Length	Depends on the acquisition time or # of triggers	

Counting tab

The Counting tab data can be saved in snapshot mode only. The generated file formats are summarized in Tab. 9.5.

Folder	Streaming Mode	Snapshot Mode
Histogram	--	3 files .txt (ASCII)
TDMS	--	.TDMS

Tab. 9.5: Counting saving scheme.

- Snapshot Mode

.txt Structure			
Typical Filename	<ul style="list-style-type: none"> • <i>ActivitymAcronym_Counts_HistoFrequency_date(mm-dd-yy) Time time(hhmm).txt</i> • <i>ActivityAcronym_Counts_HistoLeft_date(mm-dd-yy) Time time(hhmm).txt</i> • <i>ActivityAcronym_Counts_HistoRight_date(mm-dd-yy) Time time(hhmm).txt</i> 	Ex.:	<ul style="list-style-type: none"> • <i>HRDW_Counts_HistoFrequency_01-23-21 Time 1219.txt</i> • <i>HRDW_Counts_HistoLeft_01-23-21 Time 1219.txt</i> • <i>HRDW_Counts_HistoRight_01-23-21 Time 1219.txt</i>
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).		
File Properties recorder	Dataset Conditions: string record (*)		
Organisation:		Name	Type
Structure	2 columns (all)		
Channels	Channels: 2 (Histo Left and Right)	X coord	Decimal float
		Histo(X)	Integer
Channels	Channels: (HistoFrequency)	Freq. or #of pulses	Decimal float
		Event Index	Integer
Data Format	Decimal separator: point (.)		
Data Type	Array(s)		
Recorded			
Length	Depends on the number of events present in the Plot		

TDMS Structure		
Typical Filename	ActivityAcronym_Counts_Histo.TDMS	Ex.: HRDW_Counts_Histo.TDMS
General	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
File Properties recorder	Dataset Conditions: empty	
Organisation:		Name
Existing Groups	Group(s): 3	Count Histogram Left Count Histogram Right Frequency
Channels	Channel(s): total 6, (2 per group)	X coord Histo(X)
Data Format	DT_Float (floating point double precision, 64 bits) DT_Long (long 32 bits integer)	
Data Type Recorded	Array	
Length	Depends on the number of bins present in the Histograms	

(*) Dataset Conditions description

STRING Structure	
type	ASCII STRING, readable by any text editor. Variable. Terminated by \n (new line == 0x13 character). The string can be read abruptly until footer, and reproduced in any desirable context (comment, reports, screen fields,...)
HEADER	*** Start of Header ***
FOOTER	*** End of Header ***
Number of Records	56
Structure of records	Description <tab> value <nl> .OR. Description <nl> Value <nl>
Data Type Recorded	Single points of measures
Length	Fixed: 56 + 5 lines.
Contents	Example of a typical Dataset Condition header follows (with sample values) *** Start of Header *** Signature: HERA Writer V 1.0 Separator TAB Decimal Separator . Date: 05/01/2021 Time: 10.51.45 Username: Description: Dataset Channels: 2 Dataset Samples: 200 Data taking Conditions: Dataset metadata Date / Time creation: 20210105-T105144 User:

```

Activity code: 25
Sub-Activity code: 0
Sub-Activity: D2 After-Pulses studies

PSAU settings:
Channel in use: 0

Channels setting follow:
Ch.: 0
SiPM Serial: *****
Bias Voltage [V]: 55,00
Gain [dB]: 32,00
Threshold [mV]: -15,00
Channel Temperature [°C]: 25,50
T Compensation: Off
dV/dT 50,00
Ch.: 1
SiPM Serial: *****
Bias Voltage [V]: 55,00
Gain [dB]: 32,00
Threshold [mV]: -15,00
Channel Temperature [°C]: 21,00
T Compensation: Off
dV/dT 20,00

DGTZ settings:
Model: DT5720A
Model #: 9
Serial #: 812
Channel 0 Status: On
Channel 1 Status: Off
DC Offset 0: 0
DC Offset 1: 0
Trigger Mode: FALSE
Trigger Rise Time: 8
Trigger Mean: 8
Trigger 0 Thresh.: 0
Trigger 1 Thresh.: 10
Gate Mode: FALSE
Gate Width|Pre|Hold: 340 | 154 | 304
Baseline Mean|Thresh.|NoFlatTime: 1024 | 8 | 4008
Coincidence Status: FALSE
Coincidence on GPO|Time: 0 | 0
*** End of Header ***

```

Notes: <nl> stays for “new line character”. <tab> stays for Tab character.

Important Note:

Legend of the *Activity Acronyms*

- *HRDW*: Hardware Management
- *SiPM*: SiPM Experiments
- *ADV-AfterP*: Advanced Statistics Experiment (After-Pulses)
- *BETA*: Beta Spectroscopy Experiments
- *GAMMA*: Gamma Spectroscopy Experiments
- *PHOTONS*: Photons Experiments
- *COSMICS*: Cosmic Rays Experiments



Files generated during experimental activities.

The files generated during experimental activities are saved in different data formats. In addition to the previous data saving, each experimental activity generates a directory every time the activity is undergoing or has been completed. The file and directory names, the structure, and contents of these directories must not be changed. Moreover, no files can be added to those directories because it would affect data analysis procedures with different issues, including wrong results or the inability in performing the analysis.

The file formats are ASCII (.txt) and binary (.dat).

The .dat files contain direct binary copy of data in memory. Data represent a single histogram and is composed of two arrays of the same number of elements. This number depends on the number of bins included in the Histogram saved.

Arrays represent respectively the bin values sequence (float) and the counts per bin (integer).

No header or footer is present, so no data length information is present.

Arrays are aligned one after the other and data representation is, in the order: Double Precision Float (64 bits) and Long integer (32 bits). Therefore, the physical structure of the file is the following:

First array	Element 0	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	Element 1	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	Element 2	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	...	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	Element n	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
Second array	Element 0	I32 3	I32 2	I32 1	I32 0				
Second array	Element 1	I32 3	I32 2	I32 1	I32 0				
Second array	Element 2	I32 3	I32 2	I32 1	I32 0				
Second array	...	I32 3	I32 2	I32 1	I32 0				
Second array	Element n	I32 3	I32 2	I32 1	I32 0				

Tab. 9.6: Physical structure of the .dat files. Each coloured box indicates a single byte (8 bits).

Since no array length is prepend, the only way to locate and separate the two blocks is to consider that the first one must occupy the 2/3 of the total number of bytes and the second one the remaining 1/3.

Data are in Little Endian coding (Windows). Double Float is in 64-bit IEEE double-precision format.

For example, a .dat file of 72,000 bytes, contains:

1. The DB_Float array in the first 48,000 bytes
2. The I32 array in the following 24,000 bytes

And this principle must be used to locate and separate them.

- **Section n B1: Gamma Spectroscopy**

The following table reports the organization of the data files generated by HERA during the experimental activities of Section B1: Gamma Spectroscopy.

File generating experiment	Folder	Generation during the run	Description
B1.1 <u>Counting TAB</u>	21-1-DateTime-Count_aaaa..	<ul style="list-style-type: none"> • .txt (ASCII) • Dataset Condition .txt 	<ul style="list-style-type: none"> • Summary of measured counts • Dataset Condition.txt logfile
B1.2 <u>Counting TAB</u>	21-2-DateTime-Count_aaaa..	<i>Not yet implemented!</i>	
B1.3 <u>Histogram TAB</u>	21-3-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> • .txt (ASCII) • .dat (histogram binary) • Dataset Condition .txt 	<ul style="list-style-type: none"> • Histogram of charge in binary • ASCII translation of the histogram • Dataset Conditions.txt logfile
B1.4 <u>Histogram TAB</u>	21-4-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> • <i>n</i> .dat (histogram binary) • <i>n</i> .txt (ASCII) • Dataset Conditions .txt 	<ul style="list-style-type: none"> • Histogram of charge in binary • ASCII translation of histogram • Dataset Condition.txt logfile
B1.5 <u>Histogram TAB</u>	21-5-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> • <i>n</i> .dat (histogram binary) • <i>n</i> .txt (ASCII) • Dataset Conditions .txt 	<ul style="list-style-type: none"> • Histogram of charge in binary • ASCII translation of histogram • Dataset Condition.txt logfile
B1.6 <u>Histogram TAB</u>	21-6-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> • <i>n</i> .dat (histogram binary) • <i>n</i> .txt (ASCII) • Dataset Conditions .txt 	<ul style="list-style-type: none"> • Histogram of charge in binary • ASCII translation of histogram • Dataset Condition.txt logfile
B1.7 <u>Histogram TAB</u>	21-7-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> • <i>n</i> .dat (histogram binary) • <i>n</i> .txt (ASCII) • Dataset Conditions .txt 	<ul style="list-style-type: none"> • Histogram of charge in binary • ASCII translation of histogram • Dataset Condition.txt logfile
in background for all experiments	PSAU Temperatures	.txt (ASCII) New file at every change of date (Log file type)	

Tab. 9.7: Data saving scheme of the Experiments Section B1.

Activity 21.1: B1.1- Detecting γ -Radiation

.txt Structure	<i>Summary of measured counts</i>		
Typical Filename	20-2 Dark_Gamma Spectroscopy.txt <i>Never change filename and contents of this file.</i>		
General	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).		
File Properties recorder	Dataset Conditions: absent		
Organisation:	Name and value per line	Name	Type
Structure	2 columns, 7 rows (fixed)		
Separation character	TAB		
Channels	Channel(s):	Rate Src	Decimal float
		Rate Src error	Decimal float
		Rate NO Src	Decimal float
		Rate NO Src error	Decimal float
Data Format	Decimal separator: O.S. dependent (" , " or " . ")		
Data Type Recorded	Single points of measures during execution of the experiment		
Length	Fixed		

.txt Structure (Dataset Conditions)	
Filename	Dataset Conditions.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
File Properties recorder	Dataset Conditions: string record (*) without Header and Footer
Organisation:	No further structures are present
Data Format	Decimal separator: point (.)

Activity 21.3: B1.3- γ Spectrum and Energy Resolution
.txt Structure

Typical Filename	21-3-Raws_En.Spectrum -(Time or Events xxx).txt	Ex.: 21-3-Raws_En.Spectrum -(Time 300).txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name Type
Structure	2 columns	
Channels	Channel(s): 2	ADC Channel Decimal float
		Counts Integer
Data Format	Decimal separator: point (.)	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.dat Structure (histogram Binary)

Typical Filename	21-3-Raws_En.Spectrum -(Time or Events xxx).dat	Ex.: 21-3-Raws_En.Spectrum -(Time 300).dat
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
File Properties recorder	Not applicable	
Organisation:		Name Type
Structure	Cluster of 2 elements	
Elements	Arrays	Not applicable Decimal float
		Not applicable Long Integer
Data Format	Not applicable	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.txt Structure (Dataset Conditions)

Filename	Dataset Conditions.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
File Properties recorder	Dataset Conditions: string record (*) without Header and Footer
Organisation:	No further structures are present
Data Format	Decimal separator: point (.)

Activity 21.4: B1.4- System Calibration: Linearity and Resolution
.txt Structure

Typical Filename	21-4-Raws_Source -nnn-AAA (Time xxx).txt	Ex.: 21-4-Raws_Source -22-Na (Time 300).txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name Type
Structure	2 columns	
Channels	Channel(s): 2	ADC Channel Decimal float
		Counts Integer
Data Format	Decimal separator: point (.)	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.dat Structure (histogram Binary)

Typical Filename	21-4-Raws_Source -nnn-AAA (Time xxx).dat	Ex.: 21-4-Raws_Source -22-Na (Time 300).dat
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
File Properties recorder	Not applicable	
Organisation:		Name Type
Structure	Cluster of 2 elements	
Elements	Arrays	Not applicable Decimal float
		Not applicable Long Integer
Data Format	Not applicable	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.txt Structure (Dataset Conditions)

Filename	Dataset Conditions.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
File Properties recorder	Dataset Conditions: string record (*) without Header and Footer
Organisation:	No further structures are present
Data Format	Decimal separator: point (.)

Note: Conditions are replicated in the file names. Never change the names of the .bin files.

Activity 21.5: B1.5- Scintillator Crystals Comparison: Light Yield and Decay Time
.txt Structure

Typical Filename	21-5-Raws_Crystal -AAA (Time xxx).txt		Ex.: 21-5-Raws_Crystal-BGO (Time 300).txt	
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).			
File Properties	Dataset Conditions: string record (*)			
recorder				
Organisation:		Name	Type	
Structure	2 columns			
Channels	Channel(s): 2	ADC Channel	Decimal float	
		Counts	Integer	
Data Format	Decimal separator: point (.)			
Data Type Recorded	Array(s)			
Length	Depends on the number of bins present in the Histogram			

.dat Structure (histogram Binary)

Typical Filename	21-5-Raws_Crystal -AAA (Time xxx).dat		Ex.: 21-5-Raws_Crystal-BGO (Time 300).dat
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.		
File Properties recorder	Not applicable		
Organisation:		Name	Type
Structure	Cluster of 2 elements		
Elements	Arrays	Not applicable	Decimal float
		Not applicable	Long Integer
Data Format	Not applicable		
Data Type Recorded	Array(s)		
Length	Depends on the number of bins present in the Histogram		

.txt Structure (Dataset Conditions)

Filename	Dataset Conditions.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
File Properties recorder	Dataset Conditions: string record (*) without Header and Footer
Organisation:	No further structures are present
Data Format	Decimal separator: point (.)

Note: Conditions are replicated in the file names. Never change the names of the .bin files.

Activity 21.6: B1.6- γ -Radiation Absorption
.txt Structure

Typical Filename	21-6-Raws_ Tower[mm] xx Thick.yy (Time zzz).txt	Ex.: 21-6-Raws_ Tower[mm] 50 Thick.20 (Time 300).txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		
Structure	2 columns	
Channels	Channel(s): 2	ADC Channel Decimal float
		Counts Integer
Data Format	Decimal separator: point (.)	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.dat Structure (histogram Binary)

Typical Filename	21-6-Raws_ Tower[mm] xx Thick.yy (Time zzz).dat	Ex.: 21-6-Raws_ Tower[mm] 50 Thick.20 (Time 300).dat
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
File Properties recorder	Not applicable	
Organisation:		
Structure	Cluster of 2 elements	
Elements	Arrays	Not applicable Decimal float
		Not applicable Long Integer
Data Format	Not applicable	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.txt Structure (Dataset Conditions)

Filename	Dataset Conditions.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
File Properties recorder	Dataset Conditions: string record (*) without Header and Footer
Organisation:	No further structures are present
Data Format	Decimal separator: point (.)

Note: Conditions are replicated in the file names. Never change the names of the .bin files.

Activity 21.7: B1.7- Photonuclear cross-section/Compton Scattering cross-section
.txt Structure

Typical Filename	21-7-Raws_Source -xxx-AAA (Time zzz).txt	Ex.: 21-7-Raws_Source -22-Na (Time 300).txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		
Structure	2 columns	
Channels	Channel(s): 2	ADC Channel Decimal float
		Counts Integer
Data Format	Decimal separator: point (.)	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.dat Structure (histogram Binary)

Typical Filename	21-7-Raws_Source -xxx-AAA (Time zzz).dat	Ex.: 21-7-Raws_Source -22-Na (Time 300).dat
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
File Properties recorder	Not applicable	
Organisation:		
Structure	Cluster of 2 elements	
Elements	Arrays	Not applicable Decimal float
		Not applicable Long Integer
Data Format	Not applicable	
Data Type Recorded	Array(s)	
Length	Depends on the number of bins present in the Histogram	

.txt Structure (Dataset Conditions)

Filename	Dataset Conditions.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
File Properties recorder	Dataset Conditions: string record (*) without Header and Footer
Organisation:	No further structures are present
Data Format	Decimal separator: point (.)

Note: Conditions are replicated in the file names. Never change the names of the .bin files.

10 Instructions for Cleaning

The equipment may be cleaned with isopropyl alcohol or deionized water and air dried. Clean the exterior of the product only.

Do not apply cleaner directly to the items or allow liquids to enter or spill on the product.

10.1 Cleaning the Touchscreen

To clean the touchscreen (if present), wipe the screen with a towelette designed for cleaning monitors or with a clean cloth moistened with water.

Do not use sprays or aerosols directly on the screen; the liquid may seep into the housing and damage a component. Never use solvents or flammable liquids on the screen.

10.2 Cleaning the air vents

It is recommended to occasionally clean the air vents (if present) on all vented sides of the board. Lint, dust, and other foreign matter can block the vents and limit the airflow. Be sure to unplug the board before cleaning the air vents and follow the general cleaning safety precautions.

10.3 General cleaning safety precautions

CAEN recommends cleaning the device using the following precautions:

- 1) Never use solvents or flammable solutions to clean the board.
- 2) Never immerse any parts in water or cleaning solutions; apply any liquids to a clean cloth and then use the cloth on the component.
- 3) Always unplug the board when cleaning with liquids or damp cloths.
- 4) Always unplug the board before cleaning the air vents.
- 5) Wear safety glasses equipped with side shields when cleaning the board.

11 Device decommissioning

After its intended service, it is recommended to perform the following actions:

- Detach all the signal/input/output cable
- Wrap the device in its protective packaging
- Insert the device in its packaging (if present)



THE DEVICE SHALL BE STORED ONLY AT THE ENVIRONMENT CONDITIONS SPECIFIED IN THE MANUAL, OTHERWISE PERFORMANCES AND SAFETY WILL NOT BE GUARANTEED

12 Disposal

The disposal of the equipment must be managed in accordance with Directive 2012/19 / EU on waste electrical and electronic equipment (WEEE).



The crossed bin symbol indicates that the device shall not be disposed with regular residual waste.

13 Technical Support

To contact CAEN specialists for requests on the software, hardware, and board return and repair, it is necessary a MyCAEN+ account on www.caen.it:

<https://www.caen.it/support-services/getting-started-with-mycaen-portal/>

All the instructions for use the Support platform are in the document:



A paper copy of the document is delivered with CAEN boards.

The document is downloadable for free in PDF digital format at:

https://www.caen.it/wp-content/uploads/2022/11/Safety_information_Product_support_W.pdf



CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy
Phone +39 0584 388 398
Fax +39 0584 388 959
info@caen.it
www.caen.it



CAEN GmbH

Eckehardweg 10
42653 - Solingen
Germany
Phone +49 212 254 40 77
Fax +49 212 254 40 79
info@caen-de.com
www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA
Phone: +1 (718) 981-0401
Fax: +1 (718) 556-9185
info@caentechnologies.com
www.caentechnologies.com

CAENspa INDIA Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri (W)
Mumbai, Mumbai City,
Maharashtra, India, 400053
info@caen-india.in
www.caen-india.in

