


# PRELIMINARY



User Manual UM7047

## i-Spector Digital SDK

Software Development Kit

Rev. 3 - July 12th, 2024

## Purpose of this Guide

This Manual contains a brief description of the Software Development Kit for i-Spector Digital (S2570) and i-Spector PSD (S2590)

## Change Document Record

Date	Revision	Changes
July 4 <sup>th</sup> , 2019	00	Initial release
February 17 <sup>th</sup> , 2020	01	Added Chapter Rad Cloud
September 23 <sup>rd</sup> , 2020	02	Revised ordering options
July 12 <sup>th</sup> , 2024	03	Removed Chapter Rad Cloud

## Symbols, abbreviated terms and notation

ADC	Analog to Digital Converter
API	Application Programming Interface
FPGA	Field Programmable Gate Array
MCA	Multichannel Analyzer
OS	Operating system
ROI	Region of Interest
LoRa	Long Range Communication Protocol
IoT	Internet of Things

## Reference Document

[RD1] GD6745 – i-Spector family Quick Start Guide

---

CAEN S.p.A.  
Via Vetraria, 11 55049 Viareggio (LU) - ITALY  
Tel. +39.0584.388.398 Fax +39.0584.388.959  
info@caen.it  
www.caen.it

© CAEN SpA – 2024

### Disclaimer

No part of this manual may be reproduced in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of CAEN spa.

The information contained herein has been carefully checked and is believed to be accurate; however, no responsibility is assumed for inaccuracies. CAEN spa reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it).

**MADE IN ITALY:** We remark that all our boards have been designed and assembled in Italy. In a challenging environment where a competitive edge is often obtained at the cost of lower wages and declining working conditions, we proudly acknowledge that all those who participated in the production and distribution process of our devices were reasonably paid and worked in a safe environment (this is true for the boards marked "MADE IN ITALY", while we cannot guarantee for third-party manufactures).



# Index

Purpose of this Guide .....	2
Change Document Record .....	2
Symbols, abbreviated terms and notation .....	2
Reference Document .....	2
<b>Index .....</b>	<b>3</b>
<b>List of Tables .....</b>	<b>3</b>
<b>1 Introduction .....</b>	<b>4</b>
<b>2 API usage .....</b>	<b>5</b>
Set Configuration .....	6
High Voltage Configuration .....	6
MCA Configuration .....	6
Get Instrument Status .....	7
Get Spectrum .....	8
Get Waveform .....	8
Get MCA configuration .....	9
<b>3 Testing API .....</b>	<b>10</b>
<b>4 Python SDK .....</b>	<b>11</b>
Required Modules .....	11
Download the SDK .....	11
Library Usage .....	11
<b>5 C# SDK .....</b>	<b>14</b>
Required Modules .....	14
Download the SDK .....	14
Library Usage .....	14
<b>6 Technical Support .....</b>	<b>15</b>

## List of Tables

Table 1.1: table of available models and accessories .....	4
Table 2.1: table of the HTTP GET/POST endpoints available for i-Spector Digital. ....	5
Table 2.2: table of the High Voltage configurable parameters .....	6
Table 2.3: table of the MCA configurable parameters .....	6

# 1 Introduction

i-Spector Digital S2570 and S2590 models can be managed via remote control using a specific **Software Development Kit (SDK)** based on HTTP API, Python or C# coding. The SDK is available for free download.

Models compatible with SDK are listed below:

Unit	Description	Product Code
S2570E	i-Spector Digital 1" (24x24 mm <sup>2</sup> ) -OEM	WS2570EXOAAA
S2570F	i-Spector Digital 30x30 mm <sup>2</sup> -OEM	WS2570FXOAAA
S2570E	i-Spector Digital 1" (24x24 mm <sup>2</sup> ) - Csl ASSEMBLY	WS2570EXAAAA
S2570F	i-Spector Digital 30x30 mm <sup>2</sup> - Csl ASSEMBLY	WS2570FXAAAA
S2590C	i-Spector PSD 24x24 mm <sup>2</sup> – ASSEMBLY	WS2590CXAAAA
S25X0	Assembly kit and service for i-Spector OEM – new version	WS25X0ASSBXA

**Table 1.1:** table of available models and accessories.

## 2 API usage

The i-Spector is controlled through **JSON-based REST APIs**. HTTP API allows to perform all operations supported by the web-based GUI from a remote host.

i-Spector Digital uses GET/POST HTTP call to both **retrieve spectrum/wave data** and **configure the instrument**. HTTP protocol, the same used to transport web page, has the big advantage of not requiring any special configuration on the network.

Several endpoints are available, as listed in the table below.

Endpoint (WEB PAGE)	Operation	Description
/set_config.cgi	POST	Set configuration of HV/MCA/PSD
/status.cgi	GET	Read the status of the instruments and of all available processing channels (MCA/PSD)
/spectrum.cgi	GET	Get spectrum data
/wavedump.cgi	GET	Get waveform data
/psd.cgi	GET	Get PSD data
/get_mca_config.cgi	GET	Readback MCA/PSD configuration
/resetspectrum.cgi	GET	Reset spectrum
/mca_run.cgi	GET	Start MCA acquisition
/mca_stop.cgi	GET	Stop MCA acquisition
/fb_settings.cgi	GET	Retrieve fabric configuration
/get_sysx.cgi	GET	Retrieve firmware version and installed options

**Table 2.1:** table of the HTTP GET/POST endpoints available for i-Spector Digital.

In order to access to the endpoint, it is sufficient to perform HTTP GET/POST to the unit address, for example *<i-spector\_ip>/set\_config.cgi*.

Default i-Spector IP: **192.168.50.2**

i-Spector API support multiple user simultaneously connected to the instruments. Meanwhile you are developing, keep your browser open on the i-Spector Web Interface page, in order see in live the performed operation. Please remind that, if you change parameters from API, you need to refresh browser page in order to see parameters changing also in the GUI.

## Set Configuration

In order to configure the i-Spector, the `/set_config.cgi` endpoint is available. This endpoint is a POST service which must be used performing a RAW POST operation, sending in the body of the POST the configuration JSON. It is not possible to use form-data or x-www-form-urlencoded because they are not supported by i-Spector.

There are two different possible configurations to be sent to i-Spector:

- High Voltage configuration
- MCA configuration

### High Voltage Configuration

In order to configure the High Voltage, the following JSON format can be used in the body of the HTTP POST.

```
{"command": "SET_CHANNEL_CONFIG", "channel_config": [{"id": 0, "HV_STATUS": true, "HV_VOLTAGE": 41.5, "MaxV": 46, "MaxI": 5, "RAMP": 20, "TCoeff": -34, "HV_MODE": "temperature", "HV_PWRON": true}], "store_flash": false}
```

All parameters are case-sensitive.

PARAMETERS	VALID VALUES	FUNCTION
<b>HV_STATUS</b>	true/false	Enable/Disable HV
<b>HV_VOLTAGE</b>	22..80	HV voltage. Please pay attention to do not destroy the sensor setting too high voltage
<b>MaxV</b>	22..80	Max output voltage
<b>MaxI</b>	0..9	Trip Current in mA
<b>RAMP</b>	1..100	Ramp speed of HV
<b>TCoeff</b>	-1000 .. 1000	[mV/°C] Temperature compensation coefficient
<b>MODE</b>	"digital"/"temperature"	Enable / Disable temperature compensation on HV
<b>HV_PWRON</b>	true/false	Power on/off the HV on instrument boot

**Table 2.2:** table of the High Voltage configurable parameters.

### MCA Configuration

In order to configure the MCA, the following JSON format can be used in the body of the HTTP POST.

```
{"command": "SET_CHANNEL_CONFIG", "mca_config": [{"id": 0, "trigger_thrs": 28, "trigger_inib": 300, "int_pre": 300, "int_val": 10, "int_gain": 80, "pileup_inib": 30, "pileup_pen": 30, "baseline_inib": 24, "baseline_len": 256, "taget_run": 0, "taget_value": 0, "reset_on_apply": true}], "store_flash": false}
```

All parameters are case-sensitive.

PARAMETERS	VALID VALUES	FUNCTION
<b>trigger_thrs</b>	10..1000 [int]	(LSB) Trigger threshold
<b>trigger_inib</b>	10..1000 [int]	(ns) Trigger inhibit after a trigger events. (set in in order to avoid double triggers)
<b>int_pre</b>	0..1000 [int]	(ns) Charge integrator pre-trigger integration extension
<b>int_val</b>	0..100 [float]	(us) Charge integrator integration time
<b>int_gain</b>	0..1000 [int]	Charge integrator GAIN
<b>pileup_inib</b>	0..100 [float]	(us) Pileup inhibition after a trigger
<b>pileup_pen</b>	0..100 [float]	(us) Pileup penalty if a pileup event occurs
<b>baseline_inib</b>	0..100 [float]	(us) Baseline inhibition after a trigger
<b>baseline_len</b>	1024,512,256, 128,64,32,16	Length is samples of the moving average used to calculate the baseline
<b>taget_run</b>	0,1,2	Acquisition run mode 0 – FREE 1 – TIME CONTRAINED (ms) 2 – TOTAL COUNTS ON SPECTRUM
<b>taget_value</b>	[int]	Referring to taget_run parameters, this field specify the run limit. For example, to run for 10 seconds set taget_run=1 and taget_value=10000
<b>reset_on_apply</b>	true/false	Reset spectrum when one or more configuration parameters are changed

**Table 2.3:** table of the MCA configurable parameters.

## Get Instrument Status

In order to get the status of the i-Spector, it is possible to perform a HTTP GET to the `/status.cgi` endpoint. In the following, the structure returned by the endpoint.

```
{
  "command": "GET_SYSTEM_STATUS",
  "Result": "ok",
  "ErrorCode": 0,
  "Reason": "",
  "current_status": {
    "system_status": {
      "temperature": 0,
      "eth_status": 0,
      "eth_ip": "192.168.50.2",
      "last_user_interact": -1,
      "power": "wall",
      "battery": false,
      "battery_life": 0,
      "battery_charge": 0,
      "battery_in_charge": false,
      "remaining_time": 0,
      "battery_voltage": 0,
      "battery_current": 0,
      "battery_temperature": 0,
      "alarm": 0,
      "httpcloud": 0,
      "loracloud": 0
    },
    "channels": [
      {
        "id": 0,
        "HV_STATUS": true,
        "HV_VOLTAGE": 41.5,
        "HV_MODE": "temperature",
        "COMPL_V": false,
        "COMPL_I": false,
        "Vout": 42.38652,
        "Vref": 1.954313,
        "Iout": 0.3540874,
        "IoutRAW": 0.050250001,
        "Temp": 50.199402,
        "SetPoint": 42.356781,
        "ICR": 1294,
        "OCR": 1254,
        "runtime": 4542,
        "livetime": 4540,
        "sattime": 0,
        "incnt": 5856370,
        "outcnt": 5646006,
        "live": 0.969088,
        "dead": 0.030912,
        "mca_running": 1,
        "mca_status": 0
      }
    ]
  }
}
```

## Get Spectrum

In order to get the spectrum measured by the i-Spector, it is possible to perform a HTTP GET to the */spectrum.cgi* endpoint. In the following, the structure returned by the endpoint.

```
{
  "command": "GET_SPECTRUM",
  "Result": "ok",
  "ErrorCode": 0,
  "Reason": "",
  "data": [
    3,
    2,
    4,
    7,
    ..,
    1883
  ]
}
```

The data field in the JSON body contains the spectrum data in a 4096 bin array.

## Get Waveform

In order to get the waveform measured by the i-Spector, it is possible to perform a HTTP GET to the */wavedump.cgi* endpoint. In the following, the structure returned by the endpoint.

```
{
  "command": "GET_WAVEDUMP",
  "Result": "ok",
  "ErrorCode": 0,
  "Reason": "",
  "data": [
    [
      3421,
      0,
      0,
      0,
      1,
      0,
      0
    ],
    [
      3425,
      0,
      0,
      0,
      1,
      0,
      0
    ],
    [
      3425,
      0,
      0,
      0,
      1,
      0,
      0
    ]
  ]
}
```



```
[
  3423,
  0,
  0,
  0,
  1,
  0,
  0
],
....
]
```

The data field in the JSON body contains the waveform data.

Each data element is an array of 7 elements:

- Analog data
- Trigger pulse
- Charge Integration Window
- PSD Tail integration Window
- Baseline restorer status
- Pile up rejector discard
- Pile up inhibition

## Get MCA configuration

In order to get the status of the i-Spector MCA, it is possible to perform a HTTP GET to the `/get_mca_config.cgi` endpoint. In the following, the structure returned by the endpoint.

```
{
  "command": "GET_CHANNEL_CONFIGURATION",
  "Result": "ok",
  "ErrorCode": 0,
  "Reason": "",
  "mca_config": [
    {
      "id": 0,
      "trigger_thrs": 28,
      "trigger_inib": 300.000000,
      "int_pre": 300.000000,
      "int_val": 10.000000,
      "int_gain": 80.000000,
      "pileup_inib": 30.000000,
      "pileup_pen": 30.000000,
      "baseline_inib": 24.000000,
      "baseline_len": 256,
      "rebinning": 4096,
      "reset_on_apply": true,
      "taget_run": 0,
      "taget_value": 0,
      "psd_gain": 1.000000,
      "psd_delay": 0.500000,
      "psd_int": 0.800000,
      "scaleTimeWave": 0
    }
  ]
}
```

## 3 Testing API

It is suggested to use CURL to POST/GET data. Below, some example of API usage.

- EXAMPLE: Send configuration for MCA

```
curl -d '{"command" : "SET_CHANNEL_CONFIG", "mca_config" : [{"id" : 0, "trigger_thrs" : 100, "trigger_inib" : 100, "int_pre" : 300, "int_val" : 10, "int_gain" : 100, "pileup_inib" : 10, "pileup_pen" : 10, "baseline_inib" : 10, "rebinnig" : 4096, "baseline_len" : 512, "taget_run" : 0, "taget_value" : 0, "reset_on_apply" : true}], "store_flash" : false}]' -H "Content-Type: application/json" -X POST http://192.168.50.2/set_config.cgi
```

- EXAMPLE: Read system status

```
curl -X POST http://192.168.50.2/status.cgi
```

- EXAMPLE: Download spectrum

```
curl -X POST http://192.168.50.2/spectrum.cgi
```

## 4 Python SDK

A SDK for Python language is also available. The SDK uses HTTP API communication to interface with the i-Spector and requires Python >3.2. It works on x86/ia64/ARM processor. It could be used on a standard PC as well as on a Raspberry PI.

### Required Modules

The SDK requires the following modules:

- requests [pip install requests]
- enum [pip install enum]

The Example file requires the following modules:

- pprint [pip install pprint]
- numpy [pip install numpy]
- matplotlib [pip install matplotlib]

### Download the SDK

Python SDK files can be download from Nuclear Instruments Github:

<https://github.com/NuclearInstruments/InspectorSDK-Python>

or cloned with git:

```
git clone https://github.com/NuclearInstruments/InspectorSDK-Python.git
```

The SDK includes the library (inspector\_sdk.py) and an example file (test\_inspector\_sdk.py)

### Library Usage

In order to use the library, import inspector\_sdk and open a connection creating a new inspector\_sdk object:

```
from inspector_sdk import inspector_sdk
I1 = inspector_sdk("192.168.50.2")
```

The following functions are available:

- **set\_hv\_basic(self, hv\_on, hv\_voltage):**

Parameter	Type	Description
hv_on	bool	Enable/Disable HV
hv_voltage	float	HV voltage value (V)

Return: NONE

- **def set\_hv\_compensation(self, mode, temp\_coeff):**

Parameter	Type	Description
mode	HVCompensation	DISABLE_COMPENSATION: no active temperature compensation ENABLE_COMPENSATION: active temperature compensation
temp_coeff	int	SiPM temperature compensation in mV/°C

Return: NONE

- **set\_hv\_cfg(self, ramp, maxI, maxV, on\_startup):**

Parameter	Type	Description
<b>ramp</b>	int	[V/s] HV ramp speed
<b>maxI</b>	int	[mA] HV trip current
<b>maxV</b>	int	[V] Protection maximum voltage
<b>on_startup</b>	bool	Power on/off the HV on instrument boot

Return: NONE

- **configureMCA(self, trigger\_threshold, trigger\_inibit, pre\_int\_time, int\_time, int\_gain, pileup\_inib, pileup\_penalty, baseline\_inib, baseline\_len, target\_run, target\_value):**

Parameter	Type	Description
<b>trigger_threshold</b>	10..1000 [int]	(LSB) Trigger threshold
<b>trigger_inibit</b>	10..1000 [int]	(ns) Trigger inhibit after a trigger events. (set in in order to avoid double triggers)
<b>pre_int_time</b>	0..1000 [int]	(ns) Charge integrator pre-trigger integration extension
<b>int_time</b>	0..100 [float]	(us) Charge integrator integration time
<b>int_gain</b>	0..1000 [int]	Charge integrator GAIN
<b>pileup_inib</b>	0..100 [float]	(us) Pileup inhibition after a trigger
<b>pileup_pen</b>	0..100 [float]	(us) Pileup penalty if a pileup event occurs
<b>baseline_inib</b>	0..100 [float]	(us) Baseline inhibition after a trigger
<b>baseline_len</b>	[BaselineLength]	Length is samples of the moving average used to calculate the baseline
<b>target_run</b>	[RunMode]	Acquisition run mode 0 – FREE 1 – TIME CONTRAINED (ms) 2 – TOTAL COUNTS ON SPECTRUM
<b>target_value</b>	[int]	Referring to taget_run parameters, this field specify the run limit. For example to run for 10 seconds set taget_run=1 and target_value=10000

Return: NONE

- **getChannelStatus(self):**

Return: Dictionary with channel status information, as shown below

```

▼ ChStatus = (dict) <class 'dict': { 'id': 0, 'HV_STATUS'
  01 'id' (277509408) = (int) 0
  02 'HV_STATUS' (286962416) = (bool) True
  03 'HV_VOLTAGE' (286962656) = (float) 41.5
  04 'HV_MODE' (275378560) = (str) 'temperature'
  05 'COMPL_V' (286907616) = (bool) False
  06 'COMPL_I' (286907680) = (bool) False
  07 'Vout' (286907200) = (float) 42.339127
  08 'Vref' (286907424) = (float) 1.970438
  09 'Iout' (286907328) = (float) 0.31180954
  10 'IoutRAW' (286904896) = (float) 0.045125003
  11 'Temp' (286907520) = (float) 49.013733
  12 'SetPoint' (286962736) = (float) 42.316467
  13 'ICR' (286904960) = (int) 1160
  14 'OCR' (286904736) = (int) 1130
  15 'runtime' (286906304) = (int) 1897
  16 'livetime' (286962776) = (int) 1896
  17 'sattime' (286905280) = (int) 0
  18 'incnt' (286907776) = (int) 2198310
  19 'outcnt' (286907744) = (int) 2127708
  20 'live' (286906400) = (float) 0.974138
  21 'dead' (286907392) = (float) 0.025862
  22 'mca_running' (286962816) = (int) 1
  23 'mca_status' (286962856) = (int) 0

```

In order to read a particular value using the `getChannelStatus` function:

```
I1 = inspector_sdk("192.168.50.2")

ChStatus = I1.getChannelStatus()
print(ChStatus["ICR"])
```

- **`getSystemStatus(self):`**

Return: Dictionary with system status

- **`getWave(self):`**

Return list of array of array containing the waveform information.

In order to extract a column of the matrix, it is possible to use a numpy matrix and select one column (column 0 contains analog values)

```
WaveMatrix = I1.getWave()

A = np.array(WaveMatrix)
wave_track = A[:,0]
```

Each data element is an array of 7 columns:

- [0] Analog data
- [1] Trigger pulse
- [2] Charge Integration Window
- [3] PSD Tail integration Window
- [4] Baseline restorer status
- [5] Pile up rejector discard
- [6] Pile up inhibition

- **`getSpectrum(self):`**

Return: Array with 4096 spectrum bins

- **`resetSpectrum(self):`**

Return: NONE

- **`def runSpectrum(self):`**

Start spectrum acquisition

Return: NONE

- **`def stopSpectrum(self):`**

Stop spectrum acquisition

Return: NONE

## 5 C# SDK

A SDK for C# language is also available. The SDK uses HTTP API communication to interface with the i-Spector.

### Required Modules

The SDK requires the following module:

- Newtonsoft JSON

It will be automatically downloaded from NuGet at compiling time.

The Example file requires the following modules:

- pprint [pip install pprint]
- numpy [pip install numpy]
- matplotlib [pip install matplotlib]

### Download the SDK

C# SDK files can be download from Nuclear Instruments Github:

<https://github.com/NuclearInstruments/IsectorSDK-CSHARP>

or cloned with git:

```
git clone https://github.com/NuclearInstruments/IsectorSDK-CSHARP.git
```

The SDK includes a DLL library and a C# example. The DLL can be imported in any programming language supporting C# (.NET) DLL, including VB.NET, Labview, Matlab.

### Library Usage

The function in C# SDK are the same as Python SDK. Refer to **Python SDK** for usage guide.

## 6 Technical Support

CAEN makes available the technical support of its specialists for request concerning the software and the hardware.  
Use the support form available at the following link:

<https://www.caen.it/support-services/support-form/>





CAEN SpA is acknowledged as the only company in the world providing a complete range of High/Low Voltage Power Supply systems and Front-End/Data Acquisition modules which meet IEEE Standards for Nuclear and Particle Physics. Extensive Research and Development capabilities have allowed CAEN SpA to play an important, long term role in this field. Our activities have always been at the forefront of technology, thanks to years of intensive collaborations with the most important Research Centres of the world. Our products appeal to a wide range of customers including engineers, scientists and technical professionals who all trust them to help achieve their goals faster and more effectively.

**CAEN S.p.A.**

Via Vetràia, 11  
55049 Viareggio  
Italy  
Tel. +39.0584.388.398  
Fax +39.0584.388.959  
info@caen.it  
www.caen.it

**CAEN GmbH**

Klingenstraße 108  
D-42651 Solingen  
Germany  
Tel. +49 (0)212 254 4077  
Mobile +49 (0)151 16 548 484  
Fax +49 (0)212 25 44079  
info@caen-de.com  
www.caen-de.com  
CAEN GmbH

**CAEN Technologies, Inc.**

1140 Bay Street - Suite 2 C  
Staten Island, NY 10305  
USA  
Tel. +1.718.981.0401  
Fax +1.718.556.9185  
info@caentechnologies.com  
www.caentechnologies.com