

Register your device

Register your device to your **MyCAEN+** account and get access to our customer services, such as notification for new firmware or software upgrade, tracking service procedures or open a ticket for assistance. **MyCAEN+** accounts have a dedicated support service for their registered products. A set of basic information can be shared with the operator, speeding up the troubleshooting process and improving the efficiency of the support interactions.

MyCAEN+ dashboard is designed to offer you a direct access to all our after sales services. Registration is totally free, to create an account go to <https://www.caen.it/become-mycaenplus-user> and fill the registration form with your data.



<https://www.caen.it/become-mycaenplus-user/>

Index

| | |
|---|---|
| Purpose of this Manual..... | 3 |
| Change Document Record..... | 3 |
| Symbols, abbreviated terms and notation | 3 |
| Reference Documents | 3 |
| Introduction | 4 |
| Block diagram | 4 |
| Back panel components..... | 5 |
| Front panel components | 5 |
| +12V DC External Power Supply | 6 |
| Technical specifications | 6 |

Purpose of this Manual

This User's Manual contains the full description of the SP5600 Power Supply and Amplification Unit.

Change Document Record

| Date | Revision | Changes |
|-----------------|----------|---|
| 2 February 2012 | 00 | Initial release |
| 10 January 2018 | 01 | Completed revised layout and Removed Firmware Upgrade Section |

Symbols, abbreviated terms and notation

Not available

Reference Documents

Not available

CAEN S.p.A.
Via Vetraila, 11 55049 Viareggio (LU) - ITALY
Tel. +39.0584.388.398 Fax +39.0584.388.959
info@caen.it
www.caen.it

© CAEN SpA – 2018

Disclaimer

No part of this manual may be reproduced in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of CAEN SpA.

The information contained herein has been carefully checked and is believed to be accurate; however, no responsibility is assumed for inaccuracies. CAEN SpA reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it.


MADE IN ITALY: We stress the fact that all the boards are made in Italy because in this globalized world, where getting the lowest possible price for products sometimes translates into poor pay and working conditions for the people who make them, at least you know that who made your board was reasonably paid and worked in a safe environment. (this obviously applies only to the boards marked "MADE IN ITALY", we cannot attest to the manufacturing process of "third party" boards).



Introduction

The SP5600 Power Supply and Amplification Unit (PSAU) integrates two SiPM's (Silicon Photo Multiplier) mounted in single interchangeable heads, allowing easy mounting and replacement. The PSAU supplies the bias for the sensors, features a variable amplification factor up to 50 dB and integrates a feedback circuit to stabilize the sensor gain against temperature variations. Moreover, the PSAU includes one leading edge discriminator per channel and a coincidence circuit (logic unit) for flexible event trigger logic. All parameters can be programmed and monitored via USB port.

Block diagram




The diagram above shows the SP5600 functional blocks

- Microcontroller: the device is handled by a Microcontroller, PC interfaced via USB; all functional parameters are handled by the Microcontroller and therefore programmable via USB.

- **SIPM Head:** the SIPM Head (one per channel) hosts the Silicon photo multiplier; a Temperature probe inside the Head senses variations, routing them to the Microcontroller, thus allowing this one to compensate for gain instability (such function can be disabled via software).
- **Bias Generator:** the Bias Generator (one per channel), provides the bias for the Silicon photo multiplier, programmable in the $0 \div 120V$ range with 16 bit resolution (step = 1.8mV).
- **Wideband Amplifier:** the Silicon photo multiplier output is collected by the Wideband Amplifier (one per channel), whose gain is programmable ($0 \div 50$ dB); the amplifier analog output ($\pm 2V$ dynamics) is available for pulse processing on output connector (CHx), while an ultra fast comparator, as the amplifier output exceeds a programmable threshold ($\pm 800mV$, 16 bit resolution), provides a discriminator signal to the Coincidence circuit (logic unit).
- **Coincidence circuit:** the Coincidence circuit may provide on DOUTx output a logic level (NIM or TTL) in presence of comparator output as well as the coincidence (within a $10 \div 320$ ns programmable window) of the two comparators output; output level and width are programmable via USB. DOUTx signals are useful for triggering the PSAU output digitization. The unit hosts also two scalers (one per channel, 16 bit counting depth), that count the incoming pulses; their value can be readout as the programmable counting gate signal is sent. A typical application of the counters is the “Staircase”: the Dark Count Rate variation as a function of the discriminator threshold.

Back panel components



CHx OUT: Analog output, MCX connector

DOUTx: Coincidence circuit output (NIM or TTL); LEMO connector and Led: the Led is off as the analog signal is under discriminator threshold, lights up “bright” as the signal is above threshold and “weak” in the vicinities of the threshold.


Link: 3M-7610-5002 connector (reserved for future implementation)

12V: 2.1mm DC Power Jack (+12V DC Input)

ON / OFF: power switch

NIM / TTL LEDs: Green Leds, signal the DOUTx logic level

Front panel components



CHx IN: FC Connector (input)

USB Port: B type USB connector (USB 2.0 and USB 1.1 compliant)

+12V DC External Power Supply

The module is powered via the external AC/DC stabilized power supply (Mod. FRA045-S12-4, 12V DC Output, 3.75A).


Technical specifications

| Packaging | |
|---|-------------------------------------|
| Size | 154x71x43 mm ³ |
| Power Requirements | |
| External power supply | 11.5V ÷ 18V (typ. 12V) @ 0.5A |
| SiPM bias | |
| V bias | 0 ÷ 120 V range; 100 µA max current |
| V bias Setting resolution | 16bit |
| Temperature Feedback Resolution | 0.1 °C |
| Wide band amplifier and fast leading edge discriminator | |
| Gain | 0 ÷ 50 dB |
| Gain setting step | 1 dB |
| Bandwidth (-3dB) | 100 kHz ÷ 500 MHz |
| Output dynamic range | ±2V |
| Discriminator threshold | ±800mV; 25µV min. step |
| Coincidence unit and scaler | |
| Coincidence window | From 10 to 320 ns |
| Max frequency | 100 MHz |
| Delay from analog IN to DOUT | 20 ns |
| DOUT width | From 10 to 320 ns |
| DOUT level | NIM or TTL |
| Internal counters | 16 bit |
| SiPM detachable header | |
| Temperature sensor | 0.1 °C |
| Link | |
| USB interface | USB2.0 and USB1.1 compliant |



CAEN S.p.A.

Via Vetraria 11
55049 - Viareggio
Italy
Phone +39 0584 388 398
Fax +39 0584 388 959
info@caen.it
www.caen.it



CAEN GmbH

Brunnenweg 9
64331 Weiterstadt
Germany
Tel. +49 (0)212 254 4077
Mobile +49 (0)151 16 548 484
info@caen-de.com
www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA
Phone: +1 (718) 981-0401
Fax: +1 (718) 556-9185
info@caentechnologies.com
www.caentechnologies.com

CAENspa INDIA Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri (W)
Mumbai, Mumbai City,
Maharashtra, India, 400053
info@caen-india.in
www.caen-india.in

