

The **A1425** is a fast and low noise preamplifier with AC coupled input. The fast rise time makes it suited for amplifying the signals from detectors as fast silicon sensors and diamond detectors.

It has been designed for **spectroscopy applications**, and can be combined with sub-nanosecond measurements of particle time-of flight. It is optimized for high-speed single MIP particle detection with diamond detectors, where the signal integrated change is extremely small. It integrates a high voltage input for detector bias **up to 1 kV**.

- **Fast, inverting preamplifier**
- **Sensitivity of 3.6 mV/fC**
- **ENC of 0.16 fC (1000 e)**
- **Output impedance of 50 Ω**
- **Output range 0 to 1 V**
- **Up to 1 kV detector bias voltage**



## Specification

### Polarity

Negative input signals, positive inverting output

### Sensitivity

3.6 mV/fC

### Output Type

Shaped, 12 ns FWHM

### Output voltage range

0 ÷ 1 V (50 Ω termination)

### Noise

ENC of 0.16 fC (1000 e)  
measured with no input signal

### Packaging

Shielded Box

### Dimensions (WxHxD):

65 x 24 x 95 mm<sup>3</sup> (without connectors)

65 x 24 x 110 mm<sup>3</sup> (including connectors)

Weight: 134 g

### Inputs

#### IN<sup>(\*)</sup>

Detector input (AC coupled)  
SMA 142-0711-811 Johnson connector.

#### HV

HV BIAS input / Detector bias voltage  
Range: ±1000 V  
SMA 142-0711-811 Johnson connector.

#### 12 V

Power supply input connector  
Power supply voltage: +12 V DC  
2.0 mm Center Pin Diameter  
DC Power Jack Locking Type  
KLDX-0202-A-LT KYCON connector



### Outputs

#### OUT

Amplifier Out (AC decoupled)  
Dynamics: 0 ÷ 1 V (50 Ω termination)  
Output impedance: 50 Ω  
SMA 142-0711-811 Johnson connector

### Power Requirements

+12 V 30 mA (typical)

The module is powered by an external AC-DC stabilized power supply provided with the amplifier and included in the delivered kit.

#### Note:

The power jack is a 2mm type, Recommended to be used with KYCON locking plug: KLDX-PA-0202-A-LT  
A suitable cable is the RS 656-3816 type (or similar).

### (\*)Safety and Operation requirements

Care must be taken in the use of A1425 with high voltage detectors. Please remember to:

- Turn down gradually bias voltage prior to connect or disconnect preamp input
- Avoid fast changes in bias voltage
- Avoid Detector breakdown or discharge

## Ordering Option

Ordering code	Description
WA1425XAAAAA	A1425 - Fast Charge Sensitive Preamplifier





**CAEN S.p.A.**  
Via Vetraia 11  
55049 - Viareggio  
Italy  
Phone +39 0584 388 398  
Fax +39 0584 388 959  
[info@caen.it](mailto:info@caen.it)  
[www.caen.it](http://www.caen.it)



**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt  
Germany  
Tel. +49 (0)212 254 4077  
Mobile +49 (0)151 16 548 484  
[info@caen-de.com](mailto:info@caen-de.com)  
[www.caen-de.com](http://www.caen-de.com)

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA  
Phone: +1 (718) 981-0401  
Fax: +1 (718) 556-9185  
[info@caentechnologies.com](mailto:info@caentechnologies.com)  
[www.caentechnologies.com](http://www.caentechnologies.com)

**CAENspa INDIA** Private Limited

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri (W)  
Mumbai, Mumbai City,  
Maharashtra, India, 400053  
[info@caen-india.in](mailto:info@caen-india.in)  
[www.caen-india.in](http://www.caen-india.in)

