



- **Fast, low noise inverting preamplifier specifically designed for Scintillation Detectors**
- **Variable sensitivity from 0.8 to 10 mV/pC**
- **Fast output for timing measurements**
- **Test input for calibration**

The **A1424** is a preamplifier specifically designed for **Scintillation Detectors** widely used in Nuclear and High Energy Physics where low noise, fast response and high counting rates are required.

The **A1424** relies on an inverting Charge Sensitive Preamplifier which integrates both positive and negative input charge pulses coming from the Photodetector (e.g. PMT) coupled to the Scintillator.

It provides a voltage signal in the ± 4 V range on 50 Ω termination (± 8 V on 1 k Ω) with exponential decay ($\tau = 50$ μ s) as Energy output (ENERGY). The height of the resulting pulse is proportional to the integrated charge. The sensitivity of the Charge Sensitive Preamplifier can be set via a 10 position rotary switch ranging from **0.8 to 10 mV/pC**.

The **A1424** is provided with a non-inverting buffer (gain ~ 1) which reproduces the input signals coming from the detector as fast output being useful for timing measurements (FAST). Moreover, a test input accepts positive and negative signals for calibration purposes (TEST).

Specification

Input Sensitivities

0.8, 0.9, 1.1, 1.3, 1.5, 1.7, 2.5, 3, 5, 10 mV/pC selectable via Rotary Switch.

Rise Time

ENERGY Out < 60 ns
FAST Out < 2.3 ns

Fall Time

Energy Out 50 μ s

Integral Nonlinearity

Energy Out < $\pm 0.02\%$

Temperature Coefficient

$\pm 0.01\%$ /C°

Noise (Energy Out)

Max. sensitivity (10 mV/pC) < 3.2 fC
Min. sensitivity (0.8 mV/pC) < 8.5 fC

Counting Rate

Gain shift < 0.15% on 130 mV pulses on IN input (sensitivity: 3 mV/pC), adding on TEST input 60 mV @40 kHz random pulses.

Inputs

IN

Accepts positive and negative input charge pulses from Scintillation detectors. Maximum Input Voltage: ± 4 V, input impedance: 50 Ω , **LEMO-00** connector.

TEST

Positive or negative inputs for the energy calibration via Ctest = 100 pF; **LEMO-00** connector.

Power

Input Power through a 2.1 m cable with a **D-type 9 pin** male connector.




Fig. 2: Power Supply Connector pin out.

Outputs

ENERGY

± 8 V (on 1 k Ω), 50 Ω back termination. The output voltage is proportional to the amount of input charge, **LEMO-00** connector.

FAST

Reproduces the input signals (gain ~ 1) coming from the detector as fast output for timing measurements. 50 Ω back termination, **LEMO-00** connector.

Packaging

Shielded Box

Dimensions: 55 mm x 25 mm x 95 mm
Weight: 260 g (with Power Supply cable)

Power Requirements

+12 V 15 mA
-12 V 15 mA




Fig. 1: A1424 block diagram.


Ordering Option

Ordering code	Description
WA1424XAAAAA	A1424 - Scintillation Preamplifier



CAEN S.p.A.

Via Vetràia 11
55049 - Viareggio
Italy
Phone +39 0584 388 398
Fax +39 0584 388 959
info@caen.it
www.caen.it



CAEN GmbH

Brunnenweg 9
64331 Weiterstadt
Germany
Tel. +49 (0)212 254 4077
Mobile +49 (0)151 16 548 484
info@caen-de.com
www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA
Phone: +1 (718) 981-0401
Fax: +1 (718) 556-9185
info@caentechnologies.com
www.caentechnologies.com

CAENspa INDIA Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri (W)
Mumbai, Mumbai City,
Maharashtra, India, 400053
info@caen-india.in
www.caen-india.in

