




# Register your device

Register your device to your **MyCAEN+** account and get access to our customer services, such as notification for new firmware or software upgrade, tracking service procedures or open a ticket for assistance. **MyCAEN+** accounts have a dedicated support service for their registered products. A set of basic information can be shared with the operator, speeding up the troubleshooting process and improving the efficiency of the support interactions.

**MyCAEN+** dashboard is designed to offer you a direct access to all our after sales services. Registration is totally free, to create an account go to <https://www.caen.it/become-mycaenplus-user> and fill the registration form with your data.



<https://www.caen.it/become-mycaenplus-user/>

## Purpose of this Guide

This document is the HERA User Guide. It contains information about the compatibility with CAEN Educational Kits, the installation of the software, its main functions, and a brief description of the GUI.

## Change Document Record

Date	Revision	Changes
March 1 <sup>st</sup> , 2021	00	Initial release

## Symbols, Abbreviated Terms and Notation

DPP	Digital Pulse Processing
OS	Operating System
PSAU	Power Supply & Amplification Unit
DGTZ	Digitizer
GUI	Graphical User Interface
DC	Direct Current

## Reference Documents

[RD1] GD2783 - First Installation Guide to Desktop Digitizers & MCA

All CAEN documents can be downloaded at: [www.caen.it/support-services/documentation-area](http://www.caen.it/support-services/documentation-area)

---

CAEN S.p.A.  
Via Vetraia, 11 55049 Viareggio (LU) - ITALY  
Tel. +39.0584.388.398 Fax +39.0584.388.959  
[info@caen.it](mailto:info@caen.it)  
[www.caen.it](http://www.caen.it)

© CAEN SpA – 2021

### Disclaimer

No part of this manual may be reproduced in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of CAEN SpA.

The information contained herein has been carefully checked and is believed to be accurate; however, no responsibility is assumed for inaccuracies. CAEN SpA reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it).

**MADE IN ITALY:** We remark that all our boards have been designed and assembled in Italy. In a challenging environment where a competitive edge is often obtained at the cost of lower wages and declining working conditions, we proudly acknowledge that all those who participated in the production and distribution process of our devices were reasonably paid and worked in a safe environment (this is true for the boards marked "MADE IN ITALY", while we cannot guarantee for third-party manufactures).



# Index

Purpose of this Guide.....	2
Change Document Record .....	2
Symbols, Abbreviated Terms and Notation .....	2
Reference Documents .....	2
<b>Index .....</b>	<b>3</b>
<b>List of Figures.....</b>	<b>3</b>
<b>List of Tables.....</b>	<b>3</b>
<b>1. Introduction.....</b>	<b>4</b>
<b>2. Requirements &amp; Installation.....</b>	<b>6</b>
<b>3. Software Description.....</b>	<b>12</b>
Selection by Experiment .....	14
Selection by Kit .....	15
Hardware Management.....	16
Data Storage.....	16
<b>4. Technical Support.....</b>	<b>37</b>

## List of Figures

Fig. 2.1: Tracking the PSAU port assignment on a PC running Windows 10. ....	11
Fig. 3.1: Main GUI of the HERA software. ....	12
Fig. 3.2: About window. ....	13
Fig. 3.3: Bias Setting Window. ....	14
Fig. 3.4: Selection by Experiment. ....	14
Fig. 3.5: Example of experimental activity. ....	15
Fig. 3.6: HERA: Selection by kit. ....	15
Fig. 3.7: Hardware Management. ....	16

## List of Tables

Tab. 1.1: Table of the CAEN Educational Kits supported by HERA.....	5
Tab. 2.1: Host PC requirements. ....	6
Tab. 3.1: Waveforms saving scheme. ....	17
Tab. 3.2: Histograms saving scheme.....	18
Tab. 3.3: Two channels saving scheme. ....	19
Tab. 3.4: Charge vs Time saving scheme. ....	20
Tab. 3.5: Counting saving scheme. ....	20
Tab. 3.6: Physical structure of the .dat files. Each coloured box indicates a single byte (8 bits). ....	23
Tab. 3.7: Data saving scheme of the Experiments Section A1. ....	24
Tab. 3.8: Data saving scheme of the Experiments Section B1. ....	28
Tab. 3.9: Data saving scheme of the Experiments Section B2. ....	32
Tab. 3.10: Data saving scheme of the Experiments Section C2. ....	35

# 1. Introduction

CAEN brings the experience acquired in more than 40 years of collaboration with the High Energy & Nuclear Physics community into the University educational laboratories.

CAEN realized different modular Educational Kits, mainly based on Silicon Photomultipliers (SiPM) state-of-the-art sensor of light. Thanks to the most advanced instrumentation developed by CAEN for the major experiments Worldwide, together with the University teaching experience at the University of Insubria (IT), a series of experiments covering several applications has been carried out.

HERA (Handy Educational Radiation Application) is a user-friendly software developed by CAEN to make the educational product complete in every perspective. HERA allows the user to manage almost all CAEN Educational kits indeed (see **Tab. 1.1**).

The software represents a modern and flexible platform for teaching the fundamentals of Statistics, Particles Detection, and Nuclear Imaging thanks to the simple graphical interfaces and the embedded documentation.

The initial access multiplicity makes the software very flexible and suitable both for expert users as well as for beginner ones. With a simple selection, the user can decide how to execute the activity by choosing the direct access to the suggested experiments or access to devices management.

In the first case, a complete experiments suite with different difficulty levels, from basic experiments to more complex applications, is provided and the data analysis is supported by advanced tools implemented in the software itself. These tools can be enabled or disabled structuring the experiment depending on the user level. With this approach, the experiments proposed can be performed at high school level (grade 11, 12) science classes up to undergraduate physics laboratory, PhD courses and training courses. The main goal is to inspire users and guide them towards the analysis and comprehension of different physics phenomena with a series of experiments based on state-of-the-art technologies, instruments, and methods.

In the second case, the user can manage all the parameters of both the Power Supply & Amplification Unit and the Digitizer. In few easy steps, the user can control the SiPM bias voltage and the signal gain and modify the thresholds and the digital outputs. In a dedicated window, the digitized signals can be monitored for real-time fine-tuning of the setup. Energy Spectra, Charge vs time, Counting, Staircase plots, and Waveforms are also displayed together with the possibility to save the data to file for further analysis.

An interesting summary of the HERA features classified according to the type of user is the following:

- **Beginner Users**
  - A wide range of experiments covering Nuclear and Particle Physics fields: From the radioactive decays ( $\beta$  and  $\gamma$ ) to the cosmic rays, from the light quanta to the advanced statistics.
  - Frame and Pre-setting for Each Experiment.
  - Data Saving & Analysis Tools.
  - Data can be retrieved offline for additional analysis.
  - Embedded Step-by-step guides.
- **Intermediate Users**
  - Full access to hardware management.
  - Access to Experiments Frame covering Nuclear and Particle Physics fields.
  - Pre-setting for Each Experiment.
  - Data Saving.
  - Embedded Help.
- **Expert Users**
  - Full control of the Power Supply & Amplification Unit: Bias Voltage and active feedback in temperature, Gain, Thresholds, Digital Output Width, two channels coincidence.
  - Full control of the Digitizer settings: Gate, Baseline, Trigger.
  - 6 plot tabs available: Waveform, Histogram, Charge Vs. Time, Two Channels Histogram, Frequency Scan, Counting.
  - Data saving in two formats: ASCII, Binary.
  - Embedded Help.


Supported Educational Kits		
Code	Description	Picture
WK5600XEAAA	SP5600E - Educational Photon kit	
WK5600XDAAA	SP5600D - Educational Beta kit	
WSP5609XAAA	SP5600C - Educational Gamma kit	
WK5600XANAA	SP5600AN - Educational Premium kit	

Tab. 1.1: Table of the CAEN Educational Kits supported by HERA.

## 2. Requirements & Installation

The HERA software, LabVIEW™ based, is compliant with Windows 7, 8 and 10 OS, both 32 and 64 bit.


### Requirements

OS	Hardware	CAEN drivers required
 <p>Microsoft® <b>Windows</b> 10 (64-bit)</p>	2 available USB2.0 ports	DT5720 USB driver (32/64-bit) SP5600 USB driver (32/64-bit)


**Tab. 2.1:** Host PC requirements.

- Download the USB drivers for both DT5720 and SP5600 compliant to the Windows version (32 or 64-bit) on CAEN website: Educational kit webpage > “Download” > “Software” tab > Driver section (login is required before the download).
- **Install the DT5720 drivers** following the instruction of the setup wizard. The OS will automatically recognize the DT5720 when it is connected to the PC. If the automatic installation fails, perform it manually from the Device Manager by selecting the driver update and pointing to the driver folder you downloaded from CAEN website.


For example (Windows 10 – 64bit), once connected and powered on the digitizer, you can do it going to Control Panel -> System & Security -> System -> Device Manager. In the Device Manager window, find the unknown **CAEN DT5xxx USB 1.0** in the list **Other Devices**:



Right click on **CAEN DT5xxx USB 1.0** and select **Update Driver Software** option in the scroll menu.




Select **Browse My Computer** for driver software.




Click **[Browse]** to point to the Windows drivers' folder you have previously unpacked, click **[OK]** to include the path in the search and click **[Next]** to continue.





When the driver installation will be completed, click **Close** to close the window.




Refer to [RD1] for detailed installation OS-dependent.


- Connect to the PC and power ON the **SP5600**; the PC will recognize as a new peripheral by the OS. Perform the driver installation manually from the Device Manager by selecting the driver update and pointing to the driver folder you downloaded from CAEN website. Finally, a COMM port will be associated to SP5600.

For example (Windows 10 – 64bit), once connected and powered on the SP5600, you can follow the previous instructions going to Control Panel -> System & Security -> System -> Device Manager -> Controller USB [Ports (COM)] Manager.


Right click on **USB Serial Device** and select **Update Driver Software** option in the scroll menu.




Select **Browse My Computer** for driver software.




Click **[Browse]** to point to the Windows drivers' folder you have previously unpacked, click **[OK]** to include the path in the search and click **[Next]** to continue.



When the driver installation will be completed, click **Close** to close the window.



Finally, a COM port will be associated to SP5600; please check the port number as shown in **Fig. 2.1**.



**Fig. 2.1:** Tracking the PSAU port assignment on a PC running Windows 10.




**Important Note:** HERA rel. 1.0.0 Build: 1.5.21.0103 or higher:

- does not require LabVIEW™ Run-Time Engine. or LabVIEW™ version 2018 (or higher). The installation of LabVIEW™ Run-Time Engine 2018 is already implemented in the HERA.
- does not work with a digitizer USB Driver release < 3.4.7, if running in a 32-bit Windows environment.

## Software Installation

Download the standalone HERA Software revision **1.0.0 - Build: 1.5.21.0103** full installation package on CAEN website: Educational kit webpage > “Download” > “Software” tab > Application SW section (login is required before the download).

Unpack the installation package, login as administrator, launch the setup file, and complete the Installation wizard.




The setup automatically creates a link on the PC Desktop.

## 3. Software Description

When the installation procedure has been completed, the user can run the program by clicking the correspondent icon.

The opening window is shown in **Fig. 3.1**. Via this main GUI, it is possible to visualize the devices status, server messages and, last but not least, to access to data (log file, data stored, configuration files, etc.).



**Fig. 3.1:** Main GUI of the HERA software.

Before running the software, the user should wait the hardware connection. The software recognises the hardware automatically and start the connection. Two connection indicators, "Online Hardware", are present on the opening window:

- **Green light** means that the connection is ok.
- **Red light** means that there is no connection.  
If the PSAU is power on, but the light colour is red, the software can be forced to search for a new connection via the rescanning procedure from the Verbose Menu: Hardware-> PSAU Server Scanning -> Restart Scanning.
- **Yellow light** means that either the DGTZ is not a DT5720A/C, or its firmware is not compliant with Hera software, and another firmware type is probably running on the board.

The special firmware compatible with HERA Software is the Digital Pulse processing for Charge Integration for SiPM Kit (DPP-CI for SiPM) for DT5720A and the Digital Pulse Processing for Charge Integration and Pulse Shape Discrimination (DPP-PSD) for DT5720C [COMING SOON]. The firmware can be download from CAEN Website. Without any licenses, it will run in a 30-minute-per-power-cycle fully functional trial version.

To upload the firmware on the digitizer, use the CAENUpgrader Software (free download on CAEN Website): <https://www.caen.it/products/caenupgrader/>

The Main GUI clearly shows that several ways of operative openings are available:

- **"Selection by Experiment"**: access to experiments frame covering Nuclear and Particle Physics fields.
- **"Selection by Kit"**: access to operative options allowed by the educational kit in use.
- **"Hardware Management"**: direct access to the management of the device parameters and data readout.

The chosen option can be run by double clicking on the relative box or by selecting it and then by the press on the "Select" button.

The user can easily access to the GUI description via the “Help” button. Each window of the software is equipped with a dedicated “Help” button that must be closed before starting any activity.

Two points in particular deserve to be mentioned: Software Access Mode and Detector Safety.

### Software Access Mode

Three access levels to the software are implemented and are available through the “Limit Access...” selection in the drop-down menu. The first one, “Level 1”, just gives access to the Hardware Management. The second one, “Level 2”, allows the user to access the Hardware Management and the guided procedures to perform the experiments listed in the CAEN Educational Handbook. This access level does not include analysis tools. The third one, “Level 3”, gives full access to all software functionalities and all the analysis tools are included.

The initial option, “Selection by Kit”, is accessible by all the three access levels.

The user needs the Master Password to change the access mode. The Password is unique, not changeable, and not declared in the embedded Help to give this type of modification power to the tutors only.

This functionality allows to the tutors deciding what each user can access and therefore, structuring the courses depending on the course attendee levels.

The Master Password is the build of the HERA release in use. The build is displayed in the "About HERA" window via the about label in the GUI verbose menu.




Fig. 3.2: About window.

### Detector Safety

The Bias Voltage Limits can be modified from the main window of the GUI only, before selecting the experiment or hardware. The user can set the detector safety condition via the “Hardware” drop-down menu.

This functionality is very important for preventing action related to the possible detectors damaging.

The SP5600 module houses two detectors (SiPM). The module provides independent bias voltages (up to 130 V) to the sensors with gain stabilization. The user can apply a safety measure to prevent detector damage due to a wrong and too high bias voltage. Via “Safety and Serial Number Setting...” selection, the user can set the recommended operating voltage for each channel and, discretionary, the serial number to identify the detector itself. The software stabilizes the maximum value of bias voltage that can be applied to the sensor as a percentage (2,5%) of the operating one. To change the voltage limit is requested to modify the value in the “Bias Setting window” (see Fig. 3.3).




Fig. 3.3: Bias Setting Window.

## Selection by Experiment

This option allows the user to access the experiment menu listed in the CAEN Educational Handbook.

By selecting the Physics topic of interest, a series of experiments can be performed. The software programs a predefined settings of the devices and gives a detailed guide into the “Help” button.

The option “Selection by Experiment” can be run by double click on the relative icon or by selecting it and then by pressing the “Select” button.





Fig. 3.4: Selection by Experiment.



The “Help” button is present in all the windows in use and provides guides and advices about the experimental procedures.



**Fig. 3.5:** Example of experimental activity.


## Selection by Kit

This option allows the user to access the experiment menu listed in the CAEN Educational Handbook.

By selecting the Physics topic of interest, a series of experiments can be performed. The software programs a predefined settings of the devices and gives a detailed guide into the “Help” button.

The option “Selection by Experiment” can be run by double click on the relative icon or by selecting it and then by pressing the “Select” button.

The “Help” button is present in all the windows in use and provides guides and advices about the experimental procedures.



**Fig. 3.6:** HERA: Selection by kit.




## Hardware Management

The main units of the Educational kit, which are common among all the systems, are:

- Power Supply and Amplification Unit (PSAU) - SP5600
- Desktop Digitizer (DGTZ) - DT5720A

The “Hardware Management” section allows the user to manage all the parameters of both PSAU and DGTZ giving the highest flexibility in the operating modes.

With few easy steps, the setting of bias voltage, gain, thresholds, and digital outputs are possible. The digitized signals can be monitored for a real-time fine-tuning of the set-up. Energy spectra, trends of the charge as a function of the time, signal frequency versus threshold, and frequency counting are also displayed in the visualization tabs of the main GUI.



**Fig. 3.7:** Hardware Management.

Plots and data can be saved on file for the offline analysis processes.

The details about PSAU and DGTZ are present in the Help associated with their management window.

## Data Storage

The HERA system allows the user to save data in several ways:

- Save data during the run (streaming mode).
- Save the data for an offline run (snapshot mode).
- Save an image.
- Export data in Excel.

The generation of the files due to the experimental activity deserves a dedicated discussion and it is referred to the step-by-step guides of each experiment.

All directories and files are generated under the following path: *C:\ProgramData\HERA\UserName*.

Where *UserName* is the name of the logged user. If no specific username is chosen, the default name used is “Guest”.

The list of the directories created by the HERA system during Hardware Management usage is the following:

- PSAU Temperatures
- Waveform
- Histograms

- Charge-Time
- TDMS

HERA generates several file formats: ASCII (.txt), binary (.dat), and another special file format of Labview (.TDMS) format. National Instruments defined a new flexible technical data management (TDM) data model, which is accessible through LabVIEW, LabWindows™/CVI™, Measurement Studio, SignalExpress, and DIAdem.

The TDMS file format saves both the raw data and the metadata in the binary format in one file with the .TDMS extension. When creating or opening a .TDMS file, HERA automatically creates a .TDMS\_index file, used to speed up random access to the .TDMS file.

The .TDMS files can be open via a simple Add-In for Microsoft Excel (<https://www.ni.com/example/27944/en/>) or by using the "Convert data files (.TDMS)" in the File Menu.

### PSAU Temperature

Regardless of the type of acquisition and the tab in use, a new file (ASCII format) is stored at every change of date (Log file type).

.txt Structure			
Typical Filename General  File Properties recorder Organisation: Structure Separation character Channels   Data Format Data Type Recorded Length	PSAU_Temperature_date(mmm-dd-yy).txt		Ex.: PSAU Temperature_ Jan-15-21.txt
	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).		
	Dataset Conditions: absent		
		Name	Type
	5 columns		
	TAB		
	Channel(s): 5		Date
			Date (O.S. format)
			Time
			24 h format
	Board Temperature		Decimal float
	Ch.0 Temperature		Decimal float
	Ch.1 Temperature		Decimal float
	Decimal separator: point ( . )		
	Single points of measures		
	Depends on the running time of the Main Program		

### Waveform tab

In addition to the waveforms export in a Bitmap Image to the Clipboard or "Excel" numerical data via the "Export" button, the Waveforms can be saved in both modes, streaming, and snapshot. The generated files are in .TDMS format (see **Tab. 3.1**).

Folder	Streaming Mode	Snapshot Mode
Waveform	.TDMS streaming	--
TDMS	--	.TDMS

**Tab. 3.1:** Waveforms saving scheme.

- Streaming Mode

TDMS Structure		
Typical Filename	Activity Acronym_Wave(Time or events xx)_date(mm-dd-yy)-T-time(hhmm).TDMS	Ex.: HRDW_Wave(#Evn 1000)_01-08-21-T-1155.TDMS
General	TDMS Structure (NI standard), readable by Excel with “TDM Importer Plugin”.	
File Properties recorder	Dataset Conditions: string record (*)	
Organisation:		Name
Existing Groups	Group(s): 1	Analog Waveforms
Channels	Channel(s): 2	Trace Ch. 1
		Trace Ch. 0
Channel Range	0..4095	
Data Format	DT_Float (floating point double precision, 64 bits)	
Data Type	Array	
Length	Depends on the acquisition time or # of triggers	

- Snapshot Mode

<b>TDMS Structure</b>		
<b>Typical Filename</b>	<i>Activity Acronym _ Wave.TDMS</i>	Ex.: <i>SiPM_Wave.TDMS</i>
<b>General</b>	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation:</b>		<b>Name</b>
<b>Existing Groups</b>	Group(s): 2	Analog Waveforms
		Digital Waveforms
<b>Channels</b>	"Analog Waveforms" Channels: 2	Trace Ch. 1
		Trace Ch. 0
	"Digital Waveforms" Channels: 5	Virtual Probe
		Gate
		Hold Off
		Over Threshold
		Flat
<b>Channel Range</b>	0..4095	
<b>Data Format</b>	DT_Long (Long 32 bits integer)	
<b>Data Type</b>	Arrays	
<b>Length</b>	Depends on the x scale extension of the Waveform plots originating the file	

### Histogram tab

The histograms can be export in a Bitmap Image to the Clipboard and in "Excel" numerical data via the "Export" button. Moreover, as the waveforms saving, the histograms can be saved in streaming and snapshot mode. The generated file formats are summarized in **Tab. 3.2**.

Folder	Streaming Mode	Snapshot Mode
Histogram	.txt (ASCII)	.txt (ASCII) [Under request]
TDMS	--	.TDMS [Under request]

**Tab. 3.2:** Histograms saving scheme.

- Streaming Mode

<b>.txt Structure</b>		
<b>Typical Filename</b>	<i>Activity Acronym_Charge_Histo(Time or events xx)_date(mm-dd-yy)—Time(hhmm).txt</i>	Ex.: <i>HRDW_Charge_Histo(Time 10)_02-10-21 Time 1610.txt</i>
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation:</b>		<b>Name</b>
<b>Structure</b>	2 columns	
<b>Channels</b>	Channel(s): 2	ADC Channel
		Counts
		Decimal float
		Integer
<b>Data Format</b>	Decimal separator: point ( . )	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histogram	

- Snapshot Mode

<b>.txt Structure</b>		
<b>Typical Filename</b>	<i>Activity Acronym_Charge_Histo_date(mm-dd-yy)-Time-time(hhmm).txt</i>	Ex.: <i>HRDW_Charge_Histo_01-25-21 Time 1048.txt</i>
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation:</b>		<b>Name</b>
		<b>Type</b>

<b>Structure</b>	2 columns
<b>Channels</b>	Channel(s): 2
	ADC Channel
	Counts
	Decimal float
	Integer
<b>Data Format</b>	Decimal separator: point ( . )
<b>Data Type Recorded</b>	Array(s)
<b>Length</b>	Depends on the number of bins present in the Histogram

<b>TDMS Structure</b>		
<b>Typical Filename</b>	<i>Activity Acronym_Charge_Histo.TDMS</i>	<i>Ex.: HRDW_Charge_Histo.TDMS</i>
<b>General</b>	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> empty	
<b>Organisation:</b>		<b>Name</b>
<b>Existing Groups</b>	Group(s): 1	Charge Histogram
<b>Channels</b>	Channel(s): 2	X coord
		Histo(X)
<b>Data Format</b>	DT_Float (floating point double precision, 64 bits)	
	DT_Long (long 32 bits integer)	
<b>Data Type Recorded</b>	Array	
<b>Length</b>	Depends on the number of bins present in the Histograms	

### Two Channels tab

The storage of the histograms in the *Two Channels tab* can be occurred in snapshot mode only, as showed in **Tab. 3.3**.

Folder	Streaming Mode	Snapshot Mode
Histogram	--	.txt (ASCII), 2 separate files
TDMS	--	.TDMS (2 files) [Under request]

**Tab. 3.3:** Two channels saving scheme.

- Snapshot Mode

<b>.txt Structure</b>		
<b>Typical Filename</b>	Hardware Management (generic)Chx(channel number)_date(mm-dd-yy)-Time-time(hhmm).txt A file is generated for each channel.	<i>Ex.: Hardware Management (generic)Ch1_02-09-21 Time 1238.txt</i>
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation:</b>		<b>Name</b> <b>Type</b>
<b>Structure</b>	2 columns	
<b>Channels</b>	Channel(s): 2	Rate bin      Decimal float
		Counts/bin      Integer
<b>Data Format</b>	Decimal separator: point ( . )	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histograms	

<b>TDMS Structure</b>		
<b>Typical Filename</b>	<i>Activity Acronym_2Ch_Charge_Histo.TDMS</i>	<i>Ex.: HRDW_2Ch_Charge_Histo.TDMS</i>
<b>General</b>	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> empty	
<b>Organisation:</b>		<b>Name</b>
<b>Existing Groups</b>	Group(s): 2	Channel 0 Histogram
		Channel 1 Histogram
<b>Channels</b>	Channel(s): total 4, (2 per group)	X coord

	Histo(X)
<b>Data Format</b>	DT_Float (floating point double precision, 64 bits) DT_Long (long 32 bits integer)
<b>Data Type Recorded</b>	Array
<b>Length</b>	Depends on the number of bins present in the Histograms

### Charge vs Time tab

The Charge vs Time data can be saved in streaming and snapshot mode. The generated file formats are summarized in **Tab. 3.4**.

Folder	Streaming Mode	Snapshot Mode
Charge-Time	.TDMS streaming	--
TDMS	--	.TDMS (2 files) [Under request]

**Tab. 3.4:** Charge vs Time saving scheme.

- Streaming Mode

<b>TDMS Structure</b>		
<b>Typical Filename</b>	<i>Activity Acronym_ChargeVSTime(Time or events xx)_date(mm-dd-yy)-T-time(hhmm).TDMS</i>	Ex.: <i>HRDW_ChargeVSTime(Time 5)_02-09-21-T-1245.TDMS</i>
<b>General</b>	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation:</b>		<b>Name</b>
<b>Existing Groups</b>	Group(s): 1	ChargeVSTime
<b>Channels</b>	Channel(s): 2	Charge DGTZ-Ch.0 Charge DGTZ-Ch.1
<b>Data Format</b>	DT_Float (floating point double precision, 64 bits)	
<b>Data Type Recorded</b>	Array	
<b>Length</b>	Depends on the acquisition time or # of triggers	

- Snapshot Mode

<b>TDMS Structure</b>		
<b>Typical Filename</b>	<i>Activity Acronym_ChargeVSTime.TDMS</i>	Ex.: <i>HRDW_ChargeVSTime.TDMS</i>
<b>General</b>	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> empty	
<b>Organisation:</b>		<b>Name</b>
<b>Existing Groups</b>	Group(s): 1	ChargeVSTime
<b>Channels</b>	Channel(s): 2	Trace Ch. 0 Trace Ch. 1
<b>Channel Range</b>	Full range	
<b>Data Format</b>	DT_Float (floating point double precision, 64 bits)	
<b>Data Type</b>	Array	
<b>Length</b>	Depends on the acquisition time or # of triggers	

### Counting tab

The Counting tab data can be saved in snapshot mode only. The generated file formats are summarized in **Tab. 3.5**.

Folder	Streaming Mode	Snapshot Mode
Histogram	--	3 files .txt (ASCII)
TDMS	--	.TDMS

**Tab. 3.5:** Counting saving scheme.

- Snapshot Mode

#### .txt Structure

<b>Typical Filename</b>	<ul style="list-style-type: none"> <li>• <i>ActivityAcronym_Counts_HistoFrequency_date (mm-dd-yy) Time time(hhmm).txt</i></li> <li>• <i>ActivityAcronym_Counts_HistoLeft_date(mm-dd-yy) Time time(hhmm).txt</i></li> <li>• <i>ActivityAcronym_Counts_HistoRight_date(mm-dd-yy) Time time(hhmm).txt</i></li> </ul>	Ex.: <ul style="list-style-type: none"> <li>• <i>HRDW_Counts_HistoFrequency_01-23-21 Time 1219.txt</i></li> <li>• <i>HRDW_Counts_HistoLeft_01-23-21 Time 1219.txt</i></li> <li>• <i>HRDW_Counts_HistoRight_01-23-21 Time 1219.txt</i></li> </ul>
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation: Structure Channels</b>	<b>Name</b>	<b>Type</b>
	2 columns (all)	
	Channels: 2 (Histo Left and Right)	X coord Histo(X)
		Decimal float Integer
<b>Channels</b>	Channels: (HistoFrequency)	Freq. or #of pulses Event Index
		Decimal float Integer
<b>Data Format</b>	Decimal separator: point ( . )	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of events present in the Plot	

#### TDMS Structure

<b>Typical Filename</b>	<i>ActivityAcronym_Counts_Histo.TDMS</i>	Ex.: <i>HRDW_Counts_Histo.TDMS</i>
<b>General</b>	TDMS Structure (NI standard), readable by Excel with "TDM Importer Plugin".	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> empty	
<b>Organisation: Existing Groups</b>	<b>Name</b>	
	Group(s): 3	Count Histogram Left Count Histogram Right Frequency
<b>Channels</b>	Channel(s): total 6, (2 per group)	X coord Histo(X)
<b>Data Format</b>	DT_Float (floating point double precision, 64 bits) DT_Long (long 32 bits integer)	
<b>Data Type Recorded</b>	Array	
<b>Length</b>	Depends on the number of bins present in the Histograms	

#### (\*) Dataset Conditions description

##### STRING Structure

<b>type</b>	ASCII STRING, readable by any text editor. Variable. Terminated by \n (new line == 0x13 character). The string can be read abruptly until footer, and reproduced in any desirable context (comment, reports, screen fields,...)
<b>HEADER</b>	*** Start of Header ***
<b>FOOTER</b>	*** End of Header ***
<b>Number of Records</b>	56
<b>Structure of records</b>	<i>Description &lt;tab&gt; value &lt;nl&gt; .OR. Description &lt;nl&gt; Value &lt;nl&gt;</i>
<b>Data Type Recorded</b>	Single points of measures
<b>Length</b>	Fixed: 56 + 5 lines.

<b>Contents</b>	<p>Example of a typical Dataset Condition header follows (with sample values)</p> <p>*** Start of Header ***</p> <p>Signature: HERA Writer V 1.0  Separator TAB  Decimal Separator .  Date: 05/01/2021  Time: 10.51.45  Username:  Description:  Dataset Channels: 2  Dataset Samples: 200</p> <p>Data taking Conditions:</p> <p>Dataset metadata  Date / Time creation: 20210105-T105144  User:  Activity code: 25  Sub-Activity code: 0  Sub-Activity: D2 After-Pulses studies</p> <p>PSAU settings:  Channel in use: 0</p> <p>Channels setting follow:  Ch.: 0  SiPM Serial: *****  Bias Voltage [V]: 55,00  Gain [dB]: 32,00  Threshold [mV]: -15,00  Channel Temperature [°C]: 25,50  T Compensation: Off  dV/dT 50,00  Ch.: 1  SiPM Serial: *****  Bias Voltage [V]: 55,00  Gain [dB]: 32,00  Threshold [mV]: -15,00  Channel Temperature [°C]: 21,00  T Compensation: Off  dV/dT 20,00</p> <p>DGTZ settings:  Model: DT5720A  Model #: 9  Serial #: 812  Channel 0 Status: On  Channel 1 Status: Off  DC Offset 0: 0  DC Offset 1: 0  Trigger Mode: FALSE  Trigger Rise Time: 8  Trigger Mean: 8  Trigger 0 Thresh.: 0  Trigger 1 Thresh.: 10  Gate Mode: FALSE  Gate Width Pre Hold: 340   154   304  Baseline Mean Thresh. NoFlatTime: 1024   8   4008  Coincidence Status: FALSE  Coincidence on GPO Time: 0   0  *** End of Header ***</p>
-----------------	---

**Notes:** <nl> stays for “new line character”. <tab> stays for Tab character.

**Note:**

Legend of the *Activity Acronyms*

- *HRDW*: Hardware Management
- *SiPM*: SiPM Experiments
- *ADV-AfterP*: Advanced Statistics Experiment (After-Pulses)
- *BETA*: Beta Spectroscopy Experiments
- *GAMMA*: Gamma Spectroscopy Experiments
- *PHOTONS*: Photons Experiments
- *COSMICS*: Cosmic Rays Experiments



### Files generated during experimental activities.

The files generated during experimental activities are saved in different data formats. In addition to the previous data saving, each experimental activity generates a directory every time the activity is undergoing or has been completed. The file and directory names, the structure, and contents of these directories must not be changed. Moreover, no files can be added to those directories because it would affect data analysis procedures with different issues, including wrong results or the inability in performing the analysis.

The file formats are ASCII (.txt) and binary (.dat).

The .dat files contain direct binary copy of data in memory. Data represent a single histogram and is composed of two arrays of the same number of elements. This number depends on the number of bins included in the Histogram saved.

Arrays represent respectively the bin values sequence (float) and the counts per bin (integer).

No header or footer is present, so no data length information is present.

Arrays are aligned one after the other and data representation is, in the order: Double Precision Float (64 bits) and Long integer (32 bits). Therefore, the physical structure of the file is the following:

First array	Element 0	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	Element 1	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	Element 2	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	...	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
First array	Element n	DB_F 7	DB_F 6	DB_F 5	DB_F 4	DB_F 3	DB_F 2	DB_F 1	DB_F 0
Second array	Element 0	I32 3	I32 2	I32 1	I32 0				
Second array	Element 1	I32 3	I32 2	I32 1	I32 0				
Second array	Element 2	I32 3	I32 2	I32 1	I32 0				
Second array	...	I32 3	I32 2	I32 1	I32 0				
Second array	Element n	I32 3	I32 2	I32 1	I32 0				

**Tab. 3.6:** Physical structure of the .dat files. Each coloured box indicates a single byte (8 bits).

Since no array length is prepend, the only way to locate and separate the two blocks is to consider that the first one must occupy the 2/3 of the total number of bytes and the second one the remaining 1/3.

Data are in Little Endian coding (Windows). Double Float is in 64-bit IEEE double-precision format.

For example, a .dat file of 72,000 bytes, contains:

1. The DB\_Float array in the first 48,000 bytes
2. The I32 array in the following 24,000 bytes

And this principle must be used to locate and separate them.

- **Section A1: Silicon Photomultipliers**

The following table reports the organization of the data files generated by HERA during the experimental activities of Section A: Silicon Photomultipliers.

File generating experiment	Folder	Generation during the run	Description
<b>A1.1</b> <u>Histogram TAB</u>	<b>20-1-DateTime-</b> <b>Histo_aaaa..</b>	<ul style="list-style-type: none"> <li>• .dat (histogram binary)</li> <li>• .txt (ASCII)</li> <li>• Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>• Histogram of charge in binary</li> <li>• ASCII translation of the histogram</li> </ul>



			<ul style="list-style-type: none"> <li>Dataset Condition.txt logfile</li> </ul>
<b>A1.2</b> <u>Histogram TAB</u>	<b>20-2-DateTime-</b> <b>Histo_aaaa..</b>	<ul style="list-style-type: none"> <li><i>n</i>.dat (histogram binary) +</li> <li><i>n</i>.txt (ASCII) +</li> <li>Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>Histogram of charge in binary</li> <li>ASCII translation of histogram</li> <li>Dataset Condition.txt logfile.</li> </ul>
<b>A1.2</b> <u>Counting TAB</u>	<b>20-2-DateTime-</b> <b>Count_aaaa..</b>	<ul style="list-style-type: none"> <li><i>n</i>.txt (counts summary)</li> <li>Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>Summary of measured counts</li> <li>Dataset Conditions.txt logfile</li> </ul>
<b>A1.3</b> <u>Histogram TAB</u> (generates a series of directories. User can flag the directories with a mnemonic to find them during analysis)	<b>20-3-DateTime-</b> <b>Histo_aaaa..</b>	<ul style="list-style-type: none"> <li><i>n</i>.dat (histogram binary)</li> <li><i>n</i>.txt (ASCII)</li> <li>Dataset Conditions .txt</li> <li>Temperat. Monitor.wvf</li> <li>Temperat. Monitor.txt</li> </ul>	<ul style="list-style-type: none"> <li>Histogram of charge in binary</li> <li>ASCII translation of histogram</li> <li>Dataset Condition.txt logfile</li> <li>Temperature trends (binary waveform)</li> <li>Temperature trends (ASCII)</li> </ul>
<b>A1.3</b> <u>Counting TAB</u> (generates a series of directories. User can flag the directories with a mnemonic to find them during analysis)	<b>20-3-DateTime-</b> <b>Count_aaaa..</b>	<ul style="list-style-type: none"> <li><i>n</i>.txt (ASCII)</li> <li>Dataset Condition .txt</li> <li>Temperat. Monitor.wvf</li> <li>Temperat. Monitor.txt</li> </ul>	<ul style="list-style-type: none"> <li>Summary of measured counts</li> <li>Dataset Condition.txt logfile</li> <li>Temperature trends (binary waveform)</li> <li>Temperature trends (ASCII)</li> </ul>
in background for all experiments	PSAU Temperatures	.txt (ASCII) New file at every change of date (Log file type	

**Tab. 3.7:** Data saving scheme of the Experiments Section A1.

### Activity 20.1: A1.1- SiPM Characterization

#### .txt Structure

<b>Typical Filename</b>	20-1-Raws_SiPM Basic.txt		
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).		
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)		
<b>Organisation:</b>		<b>Name</b>	<b>Type</b>
<b>Structure</b>	2 columns		
<b>Channels</b>	Channel(s): 2	ADC Channel	Decimal float
		Counts	Integer
<b>Data Format</b>	Decimal separator: point ( . )		
<b>Data Type Recorded</b>	Array(s)		
<b>Length</b>	Depends on the number of bins present in the Histogram		

#### .dat Structure (histogram Binary)

<b>Typical Filename</b>	20-1-Raws_SiPM Basic.dat		
<b>General</b>	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.		
<b>File Properties recorder</b>	<b>Not applicable</b>		
<b>Organisation:</b>		<b>Name</b>	<b>Type</b>
<b>Structure</b>	Cluster of 2 elements		
<b>Elements</b>	Arrays	Not applicable	Decimal float
		Not applicable	Long Integer
<b>Data Format</b>	Not applicable		
<b>Data Type Recorded</b>	Array(s)		
<b>Length</b>	Depends on the number of bins present in the Histogram		

#### .txt Structure (Dataset Conditions)

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Activity 20.2: A1.2 - Dependence of the SiPM Properties on the Bias Voltage [Histogram TAB]**

.txt Structure				
<div>Typical Filename</div> <div>General</div> <div>File Properties recorder</div> <div>Organisation:</div> <div>Structure</div> <div>Channels</div> <div> </div> <div>Data Format</div> <div>Data Type Recorded</div> <div>Length</div>	20-2-Raws_Bias[V] xx,xx.txt		Ex.: 20-2-Raws_Bias[V] 55,20.txt	
	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).			
	Dataset Conditions: string record (*)			
			Name	Type
	2 columns			
	Channel(s): 2		ADC Channel	Decimal float
			Counts	Integer
	Decimal separator: point ( . )			
	Array(s)			
Depends on the number of bins present in the Histogram				

<b>.dat Structure (histogram Binary)</b>				
<div>Typical Filename</div> <div>General</div> <div>File Properties recorder</div> <div>Organisation:</div> <div>Structure</div> <div>Elements</div> <div> </div> <div>Data Format</div> <div>Data Type Recorded</div> <div>Length</div>	20-2-Raws_Bias[V] xx,xx.dat		Ex.: 20-2-Raws_Bias[V] 55,20.dat	
	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.			
	Not applicable			
			Name	Type
	Cluster of 2 elements			
	Arrays		Not applicable	Decimal float
			Not applicable	Long Integer
	Not applicable			
	Array(s)			
	Depends on the number of bins present in the Histogram			

.txt Structure (Dataset Conditions)	
Filename	Dataset Conditions.txt
General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
File Properties recorder	Dataset Conditions: string record (*) without Header and Footer
Organisation:	No further structures are present
Data Format	Decimal separator: point ( . )

Conditions are replicated in the file names. Never change the names of the .bin files.

**Activity 20.2: A1.2 - Dependence of the SiPM Properties on the Bias Voltage [Counting TAB]**

<b>.txt Structure (Dark Count vs Bias)</b>		<i>Summary of measured counts</i>	
Typical Filename	20-2 Dark_Bias[V]xx,yy.txt	Xx,yy indicates Bias conditions in	
	<u>Never change filename and contents of this file.</u>	Volts	
General	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).		
File Properties recorder	Dataset Conditions: absent		
Organisation:	Name and value per line	<i>Name</i>	<i>Type</i>
Structure	2 columns, 7 rows (fixed)		
Separation character	TAB		
Channels	Channel(s):	Rate 0.5	Decimal float
		Rate 0.5 error	Decimal float
		Rate 1.5	Decimal float
		Rate 1.5 error	Decimal float
		OCT	Decimal float
		OCT Error	Decimal float
Data Format	Decimal separator: O.S. dependent ( “ , ” or “ . ” )		

<b>Data Type Recorded</b>	Single points of measures during execution of the experiment
<b>Length</b>	Fixed

**.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Activity 20.3: A1.3 - Temperature Effects on SiPM Properties [Histogram TAB]**
**.txt Structure**

<b>Typical Filename</b>	20-3-Raws_Bias[V] xx,xx -Temp[°C] xx.txt	Ex.: 20-3-Raws_Bias[V] 55,00 -Temp[°C] 30.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation:</b>		
<b>Structure</b>	2 columns	<b>Name</b> <b>Type</b>
<b>Channels</b>	Channel(s): 2	ADC Channel Decimal float Counts Integer
<b>Data Format</b>	Decimal separator: point ( . )	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histogram	

**.dat Structure (histogram Binary)**

<b>Typical Filename</b>	20-3-Raws_Bias[V] xx,xx -Temp[°C] xx.dat	Ex.: 20-3-Raws_Bias[V] 55,00 -Temp[°C] 30.dat
<b>General</b>	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
<b>File Properties recorder</b>	<b>Not applicable</b>	
<b>Organisation:</b>		
<b>Structure</b>	Cluster of 2 elements	<b>Name</b> <b>Type</b>
<b>Elements</b>	Arrays	Not applicable Decimal float Not applicable Long Integer
<b>Data Format</b>	Not applicable	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histogram	

**.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

Conditions are replicated in the file names. Never change the names of the .bin files.

The file "Temperature Monitor.wvf" is used internally by the Analysis procedure. Not provided for custom analysis. Never change the name file of this file.

The "Temperature Monitor.txt" file = ASCII translation of the Waveform TM file.

**.txt Structure (Temperature Monitor)**

<b>Filename</b>	Temperature Monitor.txt
<b>General</b>	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).

<b>File Properties recorder</b>	<b>Dataset Conditions:</b> absent		
<b>Organisation:</b>		<b>Name</b>	<b>Type</b>
<b>Structure</b>	4 columns		
<b>Separation character</b>	TAB		
<b>Channels</b>	Channel(s): 4	Timestamp	Date Time (O.S. format)
		PSAU Board Y[0]	Decimal float
		PSAU Ch. 0 Y[1]	Decimal float
		PSAU Ch. 1 Y[2]	Decimal float
<b>Data Format</b>	Decimal separator: O.S. dependent ( " , " or " . " )		
<b>Data Type Recorded</b>	Single points of measures during execution of the experiment		
<b>Length</b>	Depends on the running time of the Experiment		

### Activity 20.3: A1.3 - Temperature Effects on SiPM Properties [Counting TAB]

<b>.txt Structure</b>		<i>Summary of measured counts</i>	
<b>Typical Filename</b>	20-3-Dark_Bias[V] xx,xx -Temp[°C] xx	Ex.: 20-3-Dark_Bias[V] 55,00 -Temp[°C] 28.txt	
<b>General</b>	<i>Never change filename and contents of this file.</i>		
	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).		
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> absent		
<b>Organisation:</b>	Name and value per line	<i>Name</i>	<i>Type</i>
<b>Structure</b>	2 columns, 7 rows (fixed)		
<b>Separation character</b>	TAB		
<b>Channels</b>	Channel(s):	Rate 0.5	Decimal float
		Rate 0.5 error	Decimal float
		Rate 1.5	Decimal float
		Rate 1.5 error	Decimal float
		OCT	Decimal float
		OCT Error	Decimal float
<b>Data Format</b>	Decimal separator: O.S. dependent ( “ , ” or “ . ” )		
<b>Data Type Recorded</b>	Single points of measures during execution of the experiment		
<b>Length</b>	Fixed		

<b>.txt Structure (Dataset Conditions)</b>	
<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

- Section B1: Gamma Spectroscopy**

The following table reports the organization of the data files generated by HERA during the experimental activities of Section B1: Gamma Spectroscopy.

File generating experiment	Folder	Generation during the run	Description
<b>B1.1</b> <u>Counting TAB</u>	21-1-DateTime-Count_aaaa..	<ul style="list-style-type: none"> <li>.txt (ASCII)</li> <li>Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>Summary of measured counts</li> <li>Dataset Condition.txt logfile</li> </ul>
<b>B1.2</b> <u>Counting TAB</u>	21-2-DateTime-Count_aaaa..	<i>Not yet implemented!</i>	
<b>B1.3</b> <u>Histogram TAB</u>	21-3-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> <li>.txt (ASCII)</li> <li>.dat (histogram binary)</li> <li>Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>Histogram of charge in binary</li> </ul>

			<ul style="list-style-type: none"> <li>• ASCII translation of the histogram</li> <li>• Dataset Conditions.txt logfile</li> </ul>
<b>B1.4</b> <u>Histogram TAB</u>	<b>21-4-DateTime-Histo_aaaa..</b>	<ul style="list-style-type: none"> <li>• <b>n</b>.dat (histogram binary)</li> <li>• <b>n</b>.txt (ASCII)</li> <li>• Dataset Conditions .txt</li> </ul>	<ul style="list-style-type: none"> <li>• Histogram of charge in binary</li> <li>• ASCII translation of histogram</li> <li>• Dataset Condition.txt logfile</li> </ul>
<b>B1.5</b> <u>Histogram TAB</u>	<b>21-5-DateTime-Histo_aaaa..</b>	<ul style="list-style-type: none"> <li>• <b>n</b>.dat (histogram binary)</li> <li>• <b>n</b>.txt (ASCII)</li> <li>• Dataset Conditions .txt</li> </ul>	<ul style="list-style-type: none"> <li>• Histogram of charge in binary</li> <li>• ASCII translation of histogram</li> <li>• Dataset Condition.txt logfile</li> </ul>
<b>B1.6</b> <u>Histogram TAB</u>	<b>21-6-DateTime-Histo_aaaa..</b>	<ul style="list-style-type: none"> <li>• <b>n</b>.dat (histogram binary)</li> <li>• <b>n</b>.txt (ASCII)</li> <li>• Dataset Conditions .txt</li> </ul>	<ul style="list-style-type: none"> <li>• Histogram of charge in binary</li> <li>• ASCII translation of histogram</li> <li>• Dataset Condition.txt logfile</li> </ul>
<b>B1.7</b> <u>Histogram TAB</u>	<b>21-7-DateTime-Histo_aaaa..</b>	<ul style="list-style-type: none"> <li>• <b>n</b>.dat (histogram binary)</li> <li>• <b>n</b>.txt (ASCII)</li> <li>• Dataset Conditions .txt</li> </ul>	<ul style="list-style-type: none"> <li>• Histogram of charge in binary</li> <li>• ASCII translation of histogram</li> <li>• Dataset Condition.txt logfile</li> </ul>
in background for all experiments	<b>PSAU Temperatures</b>	.txt (ASCII) New file at every change of date (Log file type)	

**Tab. 3.8:** Data saving scheme of the Experiments Section B1.

#### **Activity 21.1: B1.1- Detecting $\gamma$ -Radiation**

<b>.txt Structure</b>	<i>Summary of measured counts</i>		
<b>Typical Filename</b>	20-2 Dark_Gamma Spectroscopy.txt <i>Never change filename and contents of this file.</i>		
<b>General</b>	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).		
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> absent		
<b>Organisation:</b>	Name and value per line	<b>Name</b>	<b>Type</b>
<b>Structure</b>	2 columns, 7 rows (fixed)		
<b>Separation character</b>	TAB		
<b>Channels</b>	Channel(s):	Rate Src	Decimal float
		Rate Src error	Decimal float
		Rate NO Src	Decimal float
		Rate NO Src error	Decimal float
<b>Data Format</b>	Decimal separator: O.S. dependent ( " , " or " . " )		
<b>Data Type Recorded</b>	Single points of measures during execution of the experiment		
<b>Length</b>	Fixed		

#### **.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

### Activity 21.3: B1.3- $\gamma$ Spectrum and Energy Resolution

#### **.txt Structure**

Typical Filename		21-3-Raws_En.Spectrum –(Time or Events xxx).txt	Ex.: 21-3-Raws_En. Spectrum - (Time 300).txt	
General		ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).		
File Properties recorder		Dataset Conditions: string record (*)		
Organisation:		NameType		
Structure		2 columns		
Channels		Channel(s): 2	ADC Channel	Decimal float
			Counts	Integer
Data Format		Decimal separator: point ( . )		
Data Type Recorded		Array(s)		
Length		Depends on the number of bins present in the Histogram		

#### **.dat Structure (histogram Binary)**

Typical Filename	21-3-Raws_En.Spectrum –(Time or Events xxx).dat	Ex.: 21-3-Raws_En. Spectrum - (Time 300).dat	
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.		
File Properties recorder	Not applicable		
Organisation:		Name	Type
Structure	Cluster of 2 elements		
Elements	Arrays	Not applicable	Decimal float
		Not applicable	Long Integer
Data Format	Not applicable		
Data Type Recorded	Array(s)		
Length	Depends on the number of bins present in the Histogram		

#### **.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

### Activity 21.4: B1.4- System Calibration: Linearity and Resolution

#### **.txt Structure**

Typical Filename		21-4-Raws_Source -nnn-AAA (Time xxx).txt		Ex.: 21-4-Raws_Source -22-Na (Time 300).txt	
General		ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).			
File Properties recorder		Dataset Conditions: string record (*)			
Organisation:				Name	Type
Structure		2 columns			
Channels		Channel(s): 2		ADC Channel	Decimal float
				Counts	Integer
Data Format		Decimal separator: point ( . )			
Data Type Recorded		Array(s)			
Length		Depends on the number of bins present in the Histogram			

**.dat Structure (histogram Binary)**

Typical Filename	21-4-Raws_Source -nnn-AAA (Time xxx).dat		Ex.: 21-4-Raws_Source -22-Na (Time 300).dat
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.		
File Properties recorder	Not applicable		
Organisation:		Name	Type
Structure	Cluster of 2 elements		
Elements	Arrays	Not applicable	Decimal float
		Not applicable	Long Integer
Data Format	Not applicable		
Data Type Recorded	Array(s)		
Length	Depends on the number of bins present in the Histogram		

**.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Note:** Conditions are replicated in the file names. Never change the names of the .bin files.

**Activity 21.5: B1.5- Scintillator Crystals Comparison: Light Yield and Decay Time**
**.txt Structure**

File Properties recorder	Typical Filename	21-5-Raws_Crystal -AAA (Time xxx).txt		Ex.: 21-5-Raws_Crystal-BGO (Time 300).txt	
	General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).			
	Organisation:	Dataset Conditions: string record (*)			
	Structure				
	Channels	2 columns			
		Channel(s): 2		ADC Channel	Decimal float
				Counts	Integer
	Data Format	Decimal separator: point ( . )			
	Data Type Recorded	Array(s)			
Length	Depends on the number of bins present in the Histogram				

**.dat Structure (histogram Binary)**

Typical Filename	21-5-Raws_Crystal -AAA (Time xxx).dat		Ex.: 21-5-Raws_Crystal-BGO (Time 300).dat
General	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.		
File Properties recorder	Not applicable		
Organisation:		Name	Type
Structure	Cluster of 2 elements		
Elements	Arrays	Not applicable	Decimal float
		Not applicable	Long Integer
Data Format	Not applicable		
Data Type Recorded	Array(s)		
Length	Depends on the number of bins present in the Histogram		

### **.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Note:** Conditions are replicated in the file names. *Never change the names of the .bin files.*

### **Activity 21.6: B1.6- $\gamma$ -Radiation Absorption**

#### **.txt Structure**

<b>Typical Filename</b>	21-6-Raws_ Tower[mm] xx Thick.yy (Time zzz).txt	Ex.: 21-6-Raws_ Tower[mm] 50 Thick.20 (Time 300).txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	
<b>Organisation:</b>		
<b>Structure</b>	2 columns	
<b>Channels</b>	Channel(s): 2	ADC Channel    Decimal float Counts           Integer
<b>Data Format</b>	Decimal separator: point ( . )	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histogram	

#### **.dat Structure (histogram Binary)**

<b>Typical Filename</b>	21-6-Raws_ Tower[mm] xx Thick.yy (Time zzz).dat	Ex.: 21-6-Raws_ Tower[mm] 50 Thick.20 (Time 300).dat
<b>General</b>	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
<b>File Properties recorder</b>	<b>Not applicable</b>	
<b>Organisation:</b>		
<b>Structure</b>	Cluster of 2 elements	
<b>Elements</b>	Arrays	Not applicable    Decimal float Not applicable    Long Integer
<b>Data Format</b>	Not applicable	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histogram	

#### **.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Note:** Conditions are replicated in the file names. *Never change the names of the .bin files.*

### **Activity 21.7: B1.7- Photonuclear cross-section/Compton Scattering cross-section**

#### **.txt Structure**

<b>Typical Filename</b>	21-7-Raws_ Source -xxx-AAA (Time zzz).txt	Ex.: 21-7-Raws_ Source -22-Na (Time 300)txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)	



<b>Organisation:</b>		<b>Name</b>	<b>Type</b>
<b>Structure</b>	2 columns		
<b>Channels</b>	Channel(s): 2	ADC Channel	Decimal float
		Counts	Integer
<b>Data Format</b>	Decimal separator: point ( . )		
<b>Data Type Recorded</b>	Array(s)		
<b>Length</b>	Depends on the number of bins present in the Histogram		

#### **.dat Structure (histogram Binary)**

<b>Typical Filename</b>	21-7-Raws_Source-xxx-AAA (Time zzz).dat	Ex.: 21-7-Raws_Source-22-Na (Time 300).dat
<b>General</b>	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
<b>File Properties recorder</b>	<b>Not applicable</b>	
<b>Organisation:</b>		
<b>Structure</b>	Cluster of 2 elements	
<b>Elements</b>	Arrays	Not applicable
		Not applicable
<b>Data Format</b>	Not applicable	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histogram	

#### **.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Note:** Conditions are replicated in the file names. Never change the names of the .bin files.

- Section B2: Beta Spectroscopy**

The following table reports the organization of the data files generated by HERA during the experimental activities of Section B2: Beta Spectroscopy.

File generating experiment	Folder	Generation during the run	Description
<b>B2.1</b> <u>Counting TAB</u>	22-1-DateTime-Count_aaaa..	<ul style="list-style-type: none"> <li>.txt (ASCII)</li> <li>Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>Summary of measured counts</li> <li>Dataset Condition.txt logfile</li> </ul>
<b>B2.2</b> <u>Histogram TAB</u>	22-2-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> <li>.txt (ASCII)</li> <li>.dat (histogram binary)</li> <li>Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>Histogram of charge in binary</li> <li>ASCII translation of the histogram</li> <li>Dataset Conditions.txt logfile</li> </ul>
<b>B2.3</b> <u>Counting TAB</u>	22-3-DateTime-Count_aaaa..	<ul style="list-style-type: none"> <li><b>n</b> .txt (ASCII)</li> <li>Dataset Conditions .txt</li> </ul>	<ul style="list-style-type: none"> <li>Summary of measured counts</li> <li>Dataset Condition.txt logfile</li> </ul>
<b>B2.4</b> <u>Counting TAB</u>	22-4-DateTime-Count_aaaa..	<ul style="list-style-type: none"> <li><b>n</b> .txt (ASCII)</li> <li>Dataset Conditions .txt</li> </ul>	<ul style="list-style-type: none"> <li>Summary of measured counts</li> <li>Dataset Condition.txt logfile</li> </ul>
in background for all experiments	PSAU Temperatures	.txt (ASCII) New file at every change of date (Log file type)	

**Tab. 3.9:** Data saving scheme of the Experiments Section B2.

### Activity 22.1: B2.1- Response of a Plastic Scintillating Tile

<b>.txt Structure</b>		<i>Summary of measured counts</i>	
<b>Typical Filename</b>		20-2 Dark_Gamma Spectroscopy.txt <i>Never change filename and contents of this file.</i>	
<b>General</b>		ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>		<b>Dataset Conditions:</b> absent	
<b>Organisation:</b>		Name and value per line	<b>Name</b> <b>Type</b>
<b>.txt Structure</b>		<i>Summary of measured counts</i>	
<b>Typical Filename</b>		20-2 Dark_Posn[i].txt <i>Never change filename and contents of this file.</i>	<i>i indicates position of the radioactive source on the detector</i>
<b>General</b>		ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>		<b>Dataset Conditions:</b> absent	
<b>Organisation:</b>		Name and value per line	<b>Name</b> <b>Type</b>
<b>Structure</b>		2 columns, 7 rows (fixed)	
<b>Separation character</b>		TAB	
<b>Header</b>		Positional: Rates[Hz], Ratio[%]	
<b>Channels</b>		Channel(s):	Rate Src Rate Src error Rate NO Src Rate NO Src error Ratio Ratio error Decimal float Decimal float Decimal float Decimal float Decimal float Decimal float
<b>Data Format</b>		Decimal separator: O.S. dependent ( " , " or " . " )	
<b>Data Type Recorded</b>		Single points of measures during execution of the experiment	
<b>Length</b>		Fixed	

<b>.txt Structure (Dataset Conditions)</b>	
<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

### Activity 22.2: B2.2- $\beta$ Spectroscopy

<b>.txt Structure</b>	
<b>Typical Filename</b>	22-2-Raws_En.Spectrum -(Time or Events xxx).txt Ex.: 22-2-Raws_En. Spectrum - (Time 300).txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*)
<b>Organisation:</b>	<b>Name</b> <b>Type</b>
<b>Structure</b>	2 columns
<b>Channels</b>	Channel(s): 2 ADC Channel Counts Decimal float Integer
<b>Data Format</b>	Decimal separator: point ( . )
<b>Data Type Recorded</b>	Array(s)
<b>Length</b>	Depends on the number of bins present in the Histogram

**.dat Structure (histogram Binary)**

<b>Typical Filename</b>	22-2-Raws_En.Spectrum –(Time or Events xxx).dat	Ex.: 22-2-Raws_En. Spectrum - (Time 300).dat
<b>General</b>	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.	
<b>File Properties recorder</b>	<b>Not applicable</b>	
<b>Organisation:</b>		
<b>Structure</b>	Cluster of 2 elements	
<b>Elements</b>	Arrays	Not applicable
		Long Integer
<b>Data Format</b>	Not applicable	
<b>Data Type Recorded</b>	Array(s)	
<b>Length</b>	Depends on the number of bins present in the Histogram	

**.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Activity 22.3: B2.3-  $\beta$ -Radiation: Transmission through Matter**
**.txt Structure** *Summary of measured counts*

<b>Typical Filename</b>	22-3 Dark_NumS[n] x-Mpp MAI[mm]y.yy.txt	x indicates number of layers used; y,yy width in mm
	<u>Never change filename and contents of this file.</u>	Ex.: 22-3-Dark_NumS[n] 0-MPp[mm]1,00
<b>General</b>	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).	
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> absent	
<b>Organisation:</b>	Name and value per line	
<b>Structure</b>	2 columns, 7 rows (fixed)	
<b>Separation character</b>	TAB	
<b>Header</b>	Radiation Transmission: Rate[kHz], Dev[kHz]	
<b>Channels</b>	Channel(s):	Rate [kHz]
		Decimal float
<b>Data Format</b>	Decimal separator: O.S. dependent ( “ , “ or “ . “ )	Rate [kHz]
<b>Data Type Recorded</b>	Single points of measures during execution of the experiment	
<b>Length</b>	Fixed	

**.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Note:** Conditions are replicated in the file names. Never change the names of the .bin files.

**Activity 22.4: B2.4-  $\beta$ -Radiation as a Method to Measure Paper Sheet Grammage and Thin Layer Thickness**

<b>.txt Structure</b>		<i>Summary of measured counts</i>	
<b>Typical Filename</b>	22-3 Dark_NumS[n] x.txt	x indicates number of layers used. Ex.: 22-4-Dark_NumS[n] 0.txt	
	<u>Never change filename and contents of this file.</u>		
<b>General</b>	ASCII File, readable by any text editor. Fixed length records. Terminated by \n (new line == 0x13 character).		
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> absent		
<b>Organisation:</b>	Name and value per line	<b>Name</b>	<b>Type</b>
<b>Structure</b>	2 columns, 7 rows (fixed)		
<b>Separation character</b>	TAB		
<b>Header</b>	Radiation Transmission : Rate[kHz], Dev[kHz]		
<b>Channels</b>	Channel(s):	Rate [kHz]	Decimal float
		Rate [kHz]	Decimal float
<b>Data Format</b>	Decimal separator: O.S. dependent ( “ , “ or “ . “ )		
<b>Data Type Recorded</b>	Single points of measures during execution of the experiment		
<b>Length</b>	Fixed		

<b>.txt Structure (Dataset Conditions)</b>	
<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Note:** Conditions are replicated in the file names. Never change the names of the .bin files.

- Section C2: Photons**

The following table reports the organization of the data files generated by HERA during the experimental activities of Section C2: Photons.

File generating experiment	Folder	Generation during the run	Description
<b>C2.2</b> <u>Histogram TAB</u>	24-2-DateTime-Histo_aaaa..	<ul style="list-style-type: none"> <li>.txt (ASCII)</li> <li>.dat (histogram binary)</li> <li>Dataset Condition .txt</li> </ul>	<ul style="list-style-type: none"> <li>Histogram of charge in binary</li> <li>ASCII translation of the histogram</li> <li>Dataset Conditions.txt logfile</li> </ul>
in background for all experiments	PSAU Temperatures	.txt (ASCII) New file at every change of date (Log file type)	

**Tab. 3.10:** Data saving scheme of the Experiments Section C2.

**Activity 24.2: C2.2- Hands-on Photon Counting Statistics**

.txt Structure				
File Properties recorder	Typical Filename	24-18-Raws_LEDS[n] 0,00.txt		
	General	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).		
	Dataset Conditions:	string record (*)		
	Organisation:		Name	Type
	Structure	2 columns		
	Channels	Channel(s): 2	ADC Channel	Decimal float
			Counts	Integer
	Data Format	Decimal separator: point ( . )		
	Data Type Recorded	Array(s)		
Length	Depends on the number of bins present in the Histogram			

**.dat Structure (histogram Binary)**

<b>Typical Filename</b>	24-18-Raws_LEDS[n] 0,00.dat		
<b>General</b>	Direct Binary. Used by Analysis procedure only. Not recommended for custom analysis.		
<b>File Properties recorder</b>	<b>Not applicable</b>		
<b>Organisation:</b>		<b>Name</b>	<b>Type</b>
<b>Structure</b>	Cluster of 2 elements		
<b>Elements</b>	Arrays	Not applicable	Decimal float
		Not applicable	Long Integer
<b>Data Format</b>	Not applicable		
<b>Data Type Recorded</b>	Array(s)		
<b>Length</b>	Depends on the number of bins present in the Histogram		

**.txt Structure (Dataset Conditions)**

<b>Filename</b>	Dataset Conditions.txt
<b>General</b>	ASCII File, readable by any text editor. Variable length records. Terminated by \n (new line == 0x13 character).
<b>File Properties recorder</b>	<b>Dataset Conditions:</b> string record (*) without Header and Footer
<b>Organisation:</b>	No further structures are present
<b>Data Format</b>	Decimal separator: point ( . )

**Note:** Conditions are replicated in the file names. Never change the names of the .bin files.

## 4. Technical Support

CAEN makes available the technical support of its specialists for request concerning the software and the hardware. Use the support form available at the following link:

<https://www.caen.it/support-services/support-form/>



**CAEN S.p.A.**

Via Vetraria 11  
55049 - Viareggio  
Italy  
Phone +39 0584 388 398  
Fax +39 0584 388 959  
info@caen.it  
[www.caen.it](http://www.caen.it)

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt  
Germany  
Tel. +49 (0)212 254 4077  
Mobile +49 (0)151 16 548 484  
info@caen-de.com  
[www.caen-de.com](http://www.caen-de.com)

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA  
Phone: +1 (718) 981-0401  
Fax: +1 (718) 556-9185  
info@caentechnologies.com  
[www.caentechnologies.com](http://www.caentechnologies.com)

**CAENspa INDIA Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri (W)  
Mumbai, Mumbai City,  
Maharashtra, India, 400053  
info@caen-india.in  
[www.caen-india.in](http://www.caen-india.in)



Copyright © CAEN SpA. All rights reserved. Information in this publication supersedes all earlier versions. Specifications subject to change without notice.