



Register your device

Register your device to your **MyCAEN+** account and get access to our customer services, such as notification for new firmware or software upgrade, tracking service procedures or open a ticket for assistance. **MyCAEN+** accounts have a dedicated support service for their registered products. A set of basic information can be shared with the operator, speeding up the troubleshooting process and improving the efficiency of the support interactions.

MyCAEN+ dashboard is designed to offer you a direct access to all our after sales services. Registration is totally free, to create an account go to <https://www.caen.it/become-mycaenplus-user> and fill the registration form with your data.



<https://www.caen.it/become-mycaenplus-user/>

Index

Purpose of this Manual	3
Change Document Record	3
Symbols, abbreviated terms and notation	3
Reference Documents.....	3
Introduction	4
Available items	4
Back panel components.....	4
Front panel components	5
+12V DC External Power Supply.....	5
5mm Ultra Violet LED specifications	5
SP5605 UV-LED Driver Technical specifications.....	6
Characterization Measurements.....	6

Purpose of this Manual

This User's Manual contains the full description of the SP5605 LED driver.

Change Document Record

Date	Revision	Changes
23 October 2013	00	Initial release
17 February 2020	01	Modified "5mm Ultra Violet LED specifications" paragraph

Symbols, abbreviated terms and notation

Not available

Reference Documents

Not available

CAEN S.p.A.
Via Vetraila, 11 55049 Viareggio (LU) - ITALY
Tel. +39.0584.388.398 Fax +39.0584.388.959
info@caen.it
www.caen.it

© CAEN SpA – 2020

Disclaimer


No part of this manual may be reproduced in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of CAEN SpA.

CAEN will repair or replace any product within the guarantee period if the Guarantor declares that the product is defective due to workmanship or materials and has not been caused by mishandling, negligence on behalf of the User, accident or any abnormal conditions or operations.

CAEN declines all responsibility for damages or injuries caused by an improper use of the Modules due to negligence on behalf of the User. It is strongly recommended to read thoroughly the CAEN User's Manual before any kind of operation. *CAEN reserves the right to change partially or entirely the contents of this Manual at any time and without giving any notice.*

Disposal of the Product *The product must never be dumped in the Municipal Waste. Please check your local regulations for disposal of electronics products.*

MADE IN ITALY : We stress the fact that all the boards are made in Italy because in this globalized world, where getting the lowest possible price for products sometimes translates into poor pay and working conditions for the people who make them, at least you know that who made your board was reasonably paid and worked in a safe environment. (this obviously applies only to the boards marked "MADE IN ITALY", we cannot attest to the manufacturing process of "third party" boards).



Introduction

The SP5605 is a fast UV (ultra-violet) LED driver and this is the first available UV LED tool for SiPM tests and characterization. The SP5605 provides a triggered light burst of intensity from few photons to about 27000 photons, where the sensor saturates. The SP5605 features tuneable intensity and repetition rate. Indeed the UV LED driver can be triggered either via the internal pulse generator, or via an external source. The sensor guarantees constant wavelength signal from 0 to 1kHz. The optical signal is routed to the sensor through a fibre, FC interfaced.

Available items

Ordering Code	Item	Description
WSP5605XAAAA	SP5605	Led Driver for SiPM development kit
Thorlabs	M00279695	Optical Fiber, 600 μm thick, 1m long

Back panel components



POWER: power switch for device turning On / Off; red led, lights up as device is ON
 12V: PCB DC 2.1mm Power Socket (RAPC722X Switchcraft, +12V DC Input)

TRIGGER:

- IN: 00-type LEMO connector; External trigger input (TTL, High input impedance)
- OUT: 00-type LEMO connector; Local trigger output (TTL, drives 50 Ohm)
- INT/EXT (switch):
 - INT enables the internal trigger
 - EXT allows to use External trigger input (IN).

FREQ trimmer: frequency setting within the selected range

Front panel components



AMPLITUDE: Vishay Spectrol 10 Turn Dial with Lock lever for pulse light intensity setting

LIGHT: FC Connector (output)

+12V DC External Power Supply

The module is powered via the external AC/DC stabilized power supply (Mod. Meanwell GS40A12-P1J 40W, 12V DC Output, 3.34A).

5mm Ultra Violet LED specifications

The SP5605 LED driver is provided with a UVTOP250 TO39 UV-LED from Roithner LaserTechnik.

Absolute Maximum Rating (Ta=25 °C)			
Item	Symbol	Value	Unit
Power Dissipation, DC	P _D	180	mW
Forward Current, DC	I _F	30	mA
Operating Temperature	T _{opr}	-30 ~ +55	°C
Storage Temperature	T _{stg}	-30 ~ +100	°C
Lead Soldering Temperature	T _{sol}	190°C /5sec	-

Electrical - Optical Characteristics (Ta=25 °C)						
Item	Symbol	Condition	Min.	Typ.	Max.	Unit
DC Forward Voltage ¹	V _F	I _F =20mA	-	6.0	7.5	V
Peak. Wavelength ²	λ _p	I _F =20mA	250		260	nm
Half Width (FWHM)	Δλ	I _F =20mA	-	11	15	nm
50% Power Angle ³	2θ _{1/2}	I _F =20mA	-	7	-	deg

¹ Forward voltage tolerance ± 3%

² Wavelength tolerance ± 3nm

³ Output power tolerance ± 10%

SP5605 UV-LED Driver Technical specifications

Dimensions	Width 80mm, height 42mm, depth 80mm
Power requirements	140mA@+12V (12V DC External Power Supply)
Width of pulse	4ns typical @ full amplitude (FWHM)
Trigger Input	TTL ,High input impedance,1Hz-300kHz frequency (intensity has to be calibrated for frequencies greater than 1 kHz)
Trigger Output	TTL, drives 50 Ohm ; width: ~20ns; rise/fall time: <2ns
Output frequency (intrinsic mode)	100Hz – 1500Hz
Output light wavelength	248 ± 8 nm (measured)
Output light intensity	0 - ~27000 photons

Characterization Measurements

Complete test and characterization of the SP5605 UV-LED Driver pulses has been performed at the University of Insubria (Como, Italy).

Precise measurements of the light wavelength has been performed by the use of a spectrophotometer. The measured wavelength is 248 ± 8 nm, as shown in the following histogram. The measurement is compatible with the LED specifications of 255 ± 3 nm.




Figure 1: Emitted light spectrum wavelength. The peak is at 248 nm

The characterization of the UV-LED Driver has been performed by means of the R821 PMT from Hamamatsu Solar Blind, with known spectral and gain response.

A custom optical fiber from Thorlabs (Multimode, 600 μ m thick) has been used to drive the input light from the UV-LED to the PMT. The PMT can work in single photo-electron mode, making possible to calibrate the spectra in number of photons. The PMT quantum efficiency at 250 nm wavelength is equal to 10%.

Picture below shows the typical pulse response of UV-LED from PMT (yellow), with its synchronization pulse (green). The amplitude is proportional to the UV-LED intensity, and ranges from few tens of mV to hundreds of mV. The pulse width is defined by the PMT specifications.



Figure 2: Oscilloscope view of the UV-LED signal from PMT (yellow), and its synchronization signal (green).

The CAEN V792N QDC has been used to evaluate the number of emitted photons by the UV-LED. An external TTL trigger has been used to drive the UV-LED light emission.

In the following figure it is shown the result of the measurement of emitted light variations versus the input trigger frequency. The input trigger frequency varied from 50Hz to 1kHz. In the range of 1Hz to 1kHz the emitted wavelength is guaranteed to be stable. The resulting light intensity variations are within 15%, as can be seen from the plot. The measurements are normalized to the highest intensity value.




Figure 3: UV-LED light intensity variations in the 1Hz-1kHz range, where the wavelength is guaranteed to be stable.

Calibration of the light intensity has been performed as well, by varying the potentiometer knob. To evaluate the number of photons, the V792N QDC has been previously calibrated. Moreover, the PMT gain and quantum efficiency are known, and those values are used to convert the QDC integral value into the number of detected photo-electrons first, and finally into the number of emitted photons.

In the figure below it is shown for each value of the potentiometer knob the corresponding value of emitted photons. The measurement has been performed with the input trigger frequency of 500Hz. For frequencies in the range of 1Hz to 1kHz, the number of photons remains compatible within the differences already observed in the previous measurement.




Figure 4: Potentiometer calibration. For each value of the know the corresponding photon number is shown in the y axis.

**CAEN S.p.A.**

Via Vetraria 11
55049 - Viareggio
Italy
Phone +39 0584 388 398
Fax +39 0584 388 959
info@caen.it
www.caen.it

**CAEN GmbH**

Brunnenweg 9
64331 Weiterstadt
Germany
Tel. +49 (0)212 254 4077
Mobile +49 (0)151 16 548 484
info@caen-de.com
www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA
Phone: +1 (718) 981-0401
Fax: +1 (718) 556-9185
info@caentechnologies.com
www.caentechnologies.com

CAENspa INDIA Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri (W)
Mumbai, Mumbai City,
Maharashtra, India, 400053
info@caen-india.in
www.caen-india.in



Copyright © CAEN SpA. All rights reserved. Information in this publication supersedes all earlier versions. Specifications subject to change without notice.