


Register your device

Register your device to your **MyCAEN+** account and get access to our customer services, such as notification for new firmware or software upgrade, tracking service procedures or open a ticket for assistance. **MyCAEN+** accounts have a dedicated support service for their registered products. A set of basic information can be shared with the operator, speeding up the troubleshooting process and improving the efficiency of the support interactions.

MyCAEN+ dashboard is designed to offer you a direct access to all our after sales services. Registration is totally free, to create an account go to <https://www.caen.it/become-mycaenplus-user> and fill the registration form with your data.



<https://www.caen.it/become-mycaenplus-user/>

Purpose of this Manual

This document is the N605 - 4 Channel 200MHz Constant Fraction Discriminator User's Manual; it contains information about the installation, the configuration and the use of the unit.

Change Document Record

Date	Revision	Changes
22 June 2015	0	Preliminary
2 July 2015	1	Updated: Front and Back panel connections, Technical specifications table
21 July 2015	2	Updated: Front and Back panel connections, Technical specifications table

Symbols, abbreviated terms and notation

T.B.D.

Reference Documents

T.B.D.

Disclaimer


No part of this manual may be reproduced in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of CAEN SpA.

CAEN will repair or replace any product within the guarantee period if the Guarantor declares that the product is defective due to workmanship or materials and has not been caused by mishandling, negligence on behalf of the User, accident or any abnormal conditions or operations.

CAEN declines all responsibility for damages or injuries caused by an improper use of the Modules due to negligence on behalf of the User. It is strongly recommended to read thoroughly the CAEN User's Manual before any kind of operation. *CAEN reserves the right to change partially or entirely the contents of this Manual at any time and without giving any notice.*

Disposal of the Product *The product must never be dumped in the Municipal Waste. Please check your local regulations for disposal of electronics products.*

Made In Italy : We stress the fact that all the boards are made in Italy because in this globalized world, where getting the lowest possible price for products sometimes translates into poor pay and working conditions for the people who make them, at least you know that who made your board was reasonably paid and worked in a safe environment. (this obviously applies only to the boards marked "Made in Italy", we cannot attest to the manufacturing process of "third party" boards).



Index

1. General description	4
2. Technical specifications	5
Packaging	5
Power requirements	5
Front and Back panel	6
Front and Back panel connections.....	7
Internal components	8
Veto Disable	8
Gate Mode	9
Output Mode	9
Width Range.....	9
AC / DC Coupling	9
Technical specifications table	10
3. Operating Modes	11
Delay setting	11
Threshold setting	11
Walk adjustment.....	11
Output width adjustment	11

1. General description



The N605 houses four independent 20% Constant Fraction Discriminators (CFD) in 1-unit wide NIM. The module is intended for high resolution time measurements, where fast signals with rise time as low as 0.8 ns are involved. Those may involve the use of radiation detectors made of fast scintillators, PMT, microchannel plates and fast semiconductor detectors.

The module guarantees exceptional low walk jitter, lower than 50 ps, for signals spreading over a wide range of amplitude, from -5 V to -50 mV.

The unit is provided with independent threshold discriminators to reject signal baseline noise. Front panel test point and screwdriver permits the precise measurement of the threshold, in the range from -15 mV down to -1 V. Each channel can be adjusted for the zero crossing level of the discriminator to minimize the time walk. Moreover the N605 has a CFD signal monitor output connector on the front panel to facilitate the correct adjustment of the working parameters.

Each channel has a fan-out of three timing output signals provided with NIM standard. The outputs can be selected to operate both in blocking and updating mode: the former means that, while a logic output is produced, further incoming pulses are neglected; the latter means that if a new input pulse arrives while a logic output is produced, such output is extended for the corresponding amount of time, as long as the input is over discriminator threshold.

The board is provided with fast VETO input and individual GATE inputs to operate coincidence or anticoincidence logic conditions with the timing signal.

2. Technical specifications

Packaging



NIM-standard single width module. Weight: ~1 kg.

Power requirements


Power absorption requirements are the following:

+6V	300mA
-6V	1600mA
+12V	-
-12V	100mA
+24V	-
-24V	-


Front and Back panel



Front and Back panel connections



IN	LEMO 00 connector	Analog Input	Analog signal from detector for constant-fraction timing. Linear range from -20 mV to -5 V. 50 Ohm, DC coupled; protected with diode clamps at ± 10 V. Input reflections <10% for input rise times >2 ns
DLY	2 LEMO 00 connectors	Delay stub	Allows for selecting the CF shaping delay. 50-Ohm cable to be inserted between connectors. For triggering at 20% fraction, the shaping delay duration must be set equal to input pulse rise time plus 1ns.
MON	LEMO 00 connector	Monitor output	Allows monitoring of the CF shaped signal. 50 Ohm, AC coupled. The MON output is proportional to the input signal CF minus the delayed input; the zero crossing threshold is compared to the input signal for minimized walk.
OUT	3 LEMO 00 connectors	Discriminator output	Three bridged outputs, updating or blocking, fast negative NIM output signals, mark the CF zero-crossing time. Amplitude -800 mV on 50-Ohm load. Each connector has its own 50-Ohm resistor in series with the common output driver. Unused output connectors must be terminated on 50 Ohm.
T	20-turn trimmer; test point	Threshold setting	Sets the minimum pulse amplitude that will produce a timing output. Variable from -15 to -1000 mV. The relevant test point monitors the discriminator threshold setting. The test point voltage is 10X the actual threshold setting. $Z_{out} \leq 2$ kOhm.
Z	20-turn trimmer; test point	Walk adjustment	Adjustment for fine-tuning the zero crossing discriminator threshold to achieve minimum walk. Adjustable over a ± 15 mV range. The relevant test point monitors the actual setting of the zero-crossing discriminator. The test point voltage is 100X the actual zero crossing discriminator threshold setting. Z_{out} : 1 kOhm.
W	20-turn trimmer	Output width adjustment	Adjustment for each discriminator channel sets the width of the three output logic pulses. The range of width adjustment depends on the positions of jumpers, see p. 9
ACT	Red Led	Activity	Signals presence of output



VETO	LEMO 00 connector	Veto input	Accepts negative NIM level to inhibit the timing outputs on all the channels chosen with the internal VETO Disable switch. Input impedance, 50 Ohm. For proper FAST VETO operation, the leading edge of the VETO signal must precede the IN1 (IN2, IN3, or IN4) signal by 3 ns and have a width equal to the CF Shaping Delay plus 5 ns
GND	test point	Probe GND	Front-panel test point provides a ground connection for test probes




GATE	BNC	Gate IN	Accepts a negative NIM level to gate the respective constant-fraction timing output. Coincidence or anticoincidence gating is selected via jumper (see p.9). Input impedance, 50Ohm. GATE INPUT leading edge should precede the Analog Input by 1 ns and have a width equal to the CF Shaping Delay plus 5 ns.
-------------	-----	---------	--


Internal components

Several settings can be made via internal jumpers (whose position is indicated by coloured boxes in the figures below), accessible by removing the unit side covers.


Component side:



Soldering side:




Veto Disable




One jumper per channel (Ch0: J30; Ch1: J31; Ch2: J32; Ch3: J33) on PCB component side (orange box); short pin 1 and pin 2 in order to enable VETO (the channel output is disabled in presence of VETO signal); short pin 2 and pin 3 in order to disable it (channel output ignores the presence of VETO signal).

Gate Mode



Two jumpers per channel (Ch0: J37, J36; Ch1: J34, J35; Ch2: J41, J40; Ch3: J38, J39) on PCB component side (green box), allow to choose either Coincidence or anticoincidence gating. Insert GATE# C jumper to select “coincidence” and GATE# A jumper to select “anticoincidence”.

Output Mode




The outputs pulse can be either updating or blocking. One short circuitable socket per channel (Ch0: J13; Ch1: J14; Ch2: J15; Ch3: J16) on PCB component side (blue boxes), allows to choose the desired modality: short pin 1 and pin 2 in order to select Blocking mode; short pin 2 and pin 3 in order to select Updating mode.

Width Range



The output width range selection is achieved through the combination of a couple of short circuitable sockets per channel, on PCB component side: “Width range A”, one per channel (Ch0: J9; Ch1: J10; Ch2: J11; Ch3: J12) [red boxes] and “Width range B” (Ch0: J90; Ch1: J92; Ch2: J94; Ch3: J96) [purple boxes]; the combination result is explained at page 11.

AC / DC Coupling



One short circuitable socket per channel (Ch0: J46; Ch1: J48; Ch2: J50; Ch3: J52) on PCB soldering side (pink boxes), allows to choose the desired modality: short pin 1 and pin 2 in order to select DC coupling; contact open for AC coupling.

Technical specifications table

Packaging	1U-wide NIM unit
Input Channels	4 inputs (negative polarity, 50 Ohm impedance)
Input linear range	-20mV ÷ -5V
Max input voltage	-10V
Max input frequency	250 MHz
Threshold range	-15 mV ÷ -1 V
Output width	From 2ns to 500ns (jumper setting)
Output rise time	<600ps
Output fall time	<950ps
Walk	<±30 ps on -50mV ÷ -5V input with 2.5ns rise/fall time <±25 ps on -50mV ÷ -2V input with 1.6ns rise/fall time
Constant Fraction	20%
Double Pulse Resolution	<4 ns
Transmission Delay	Typically <7 ns with 3ns external delay.
Operating Temperature Range	0 to 50°C.
Threshold Temperature Sensitivity	t.b.d.
Transmission delay Temperature Sensitivity	<±10 ps/°C from 0 to 50°C.

3. Operating Modes

Delay setting

The Constant Fraction Discrimination technique is based on summing a delayed, full height input signal to an inverted and attenuated signal. The resulting signal is fed into a zero-crossing comparator, thus obtaining a precise timing information that eliminates any walk error induced by constant rise time and varying amplitude signals.

For correct operation, the maximum of the attenuated pulse has to cross the delayed pulse at the selected fraction. This condition leads to a delay setting equal to the input rise time plus 1ns; the correct delay value is obtained by inserting a cable segment of proper length between the front panel Delay connectors.

Threshold setting

The N605 produces an output signal as the analog input signal exceeds the threshold, i.e. the lowest energy of interest.

To adjust the threshold level of a channel, measure the dc voltage between the front-panel Threshold test point and ground for the relevant channel. The ground test point is located next to Veto Input connector. The nominal range of voltages at the test point is -150 mV to -10 V, corresponding to the actual threshold which is 10% of the test point voltage, i.e. a 5V reading leads to a 500mV threshold. Use a screwdriver to adjust the threshold level with the "T" trimmer.

Walk adjustment

To adjust the Walk characteristics, connect the signal source to be used to the Input connector in the active channel and connect the signal from the constant-fraction Monitor connector to the oscilloscope (bandwidth > 300 MHz) through a 50-Ohm cable. Select the CF Shaping Delay according to the instruction given above. The Bipolar signal can be observed on the oscilloscope, externally triggered by the corresponding discriminator output. Adjust the Walk (Z) control, which sets the zero-crossing reference, so that the bipolar constant-fraction signals for all input amplitudes cross the baseline at about the same time. The relevant test point can be used for fixing the zero-crossing reference; the nominal range of voltages at the test point is -1.5 V to 1.5 V, corresponding to the actual threshold which is 1% of the test point voltage. Use a screwdriver to adjust the "Z" trimmer.


Output width adjustment

In order to adjust the output width, provide an input pulse that exceeds the threshold at a rate lower than 500 KHz and measure the width of the relevant output pulse. Terminate the unused output connectors on 50 Ohm. Use a screwdriver to set the "W" trimmer for the output width in the active channel. The use of Output Mode and Width jumpers (see p.9), allows to adjust the width in the following ranges:

Width jumper A	Width jumper B	Min. Width updating	Max. Width updating
Open	Closed	2 ns	60 ns
Closed	Closed	10 ns	240 ns
Closed	Open	240 ns	500 ns

**CAEN S.p.A.**

Via Vetràia 11
55049 - Viareggio
Italy
Phone +39 0584 388 398
Fax +39 0584 388 959
info@caen.it
www.caen.it

**CAEN GmbH**

Brunnenweg 9
64331 Weiterstadt
Germany
Tel. +49 (0)212 254 4077
Mobile +49 (0)151 16 548 484
info@caen-de.com
www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA
Phone: +1 (718) 981-0401
Fax: +1 (718) 556-9185
info@caentechnologies.com
www.caentechnologies.com

CAENspa INDIA Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri (W)
Mumbai, Mumbai City,
Maharashtra, India, 400053
info@caen-india.in
www.caen-india.in



Copyright © CAEN SpA. All rights reserved. Information in this publication supersedes all earlier versions. Specifications subject to change without notice.