

# Register your device

Register your device to your **MyCAEN+** account and get access to our customer services, such as notification for new firmware or software upgrade, tracking service procedures or open a ticket for assistance. **MyCAEN+** accounts have a dedicated support service for their registered products. A set of basic information can be shared with the operator, speeding up the troubleshooting process and improving the efficiency of the support interactions.

**MyCAEN+** dashboard is designed to offer you a direct access to all our after sales services. Registration is totally free, to create an account go to <https://www.caen.it/become-mycaenplus-user> and fill the registration form with your data.



<https://www.caen.it/become-mycaenplus-user/>

# Technical Information Manual

Revision n. 3  
20 October 2010

**MOD. N914**

*8 FOLD CHARGE  
PREAMPLIFIER &  
DISCRIMINATOR*

**NPO:**  
**00122/06:N914x.MUTx/03**

CAEN will repair or replace any product within the guarantee period if the Guarantor declares that the product is defective due to workmanship or materials and has not been caused by mishandling, negligence on behalf of the User, accident or any abnormal conditions or operations.

**CAEN declines all responsibility for damages or injuries caused by an improper use of the Modules due to negligence on behalf of the User. It is strongly recommended to read thoroughly the CAEN User's Manual before any kind of operation.**



*CAEN reserves the right to change partially or entirely the contents of this Manual at any time and without giving any notice.*

### **Disposal of the Product**

*The product must never be dumped in the Municipal Waste. Please check your local regulations for disposal of electronics products.*



**MADE IN ITALY** : We stress the fact that all the boards are made in Italy because in this globalized world, where getting the lowest possible price for products sometimes translates into poor pay and working conditions for the people who make them, at least you know that who made your board was reasonably paid and worked in a safe environment. (this obviously applies only to the boards marked "MADE IN ITALY", we can not attest to the manufacturing process of "third party" boards).

---

## ***TABLE OF CONTENTS***

<b>1. GENERAL DESCRIPTION.....</b>	<b>4</b>
1.1 OVERVIEW .....	4
<b>2. TECHNICAL SPECIFICATIONS.....</b>	<b>6</b>
2.1 PACKAGING.....	6
2.2 POWER REQUIREMENTS .....	6
2.3 FRONT PANEL .....	7
2.4 INPUT/OUTPUT CONNECTIONS .....	8
2.5 FRONT PANEL TRIMMERS.....	9
2.6 INTERNAL JUMPERS .....	9
2.7 TECHNICAL SPECIFICATION TABLE .....	11

---

## ***LIST OF FIGURES***

FIG. 1.1: MOD. N914 BLOCK DIAGRAM .....	5
FIG. 2.1: MOD. N914 FRONT AND REAR PANEL.....	7
FIG. 2.2: INTERNAL JUMPERS LOCATION .....	10

---

## ***LIST OF TABLES***

TABLE 2.1: POWER REQUIREMENTS.....	6
TABLE 2.2: MOD. N914 TECHNICAL FEATURES .....	11

---

# 1. General description

---

## 1.1 Overview

The Model N914 is a 8 channel charge preamplifier and integrator, implemented in a single-width NIM module.


It is designed to process signals generated by Photomultipliers.

The module is provided with LEMO 00 connectors for both input and output signals, except the T# output which is provided through a 16 pin flat cable connector.

The Photomultiplier signal at the input is buffered and processed in several ways:

- The OL# output provides the signal processed by a gateless charge preamplifier and integrator
- The OH# output provides the signal processed by a gateless charge preamplifier and integrator, followed by a 10x gain stage
- The OLSUM output provides the analog sum of the OL# outputs
- The OHSUM output provides the analog sum of the OH# outputs
- The ASUM output provides the analog sum of the CH# inputs
- The MALU output provides a step function (Majority) with a height proportional to the number of channels simultaneously over a threshold settable individually or globally.
- The T\_OUT# output provides a discriminator output converted into a digital signal to be sent to a TDC for arrival times recording

Fig. 1.1 shows a block diagram of the module.



**Fig. 1.1: Mod. N914 block diagram**

---

## 2. Technical specifications

---

### 2.1 Packaging

The Model N914 is housed in a single width NIM module.

---

### 2.2 Power requirements

**Table 2.1: Power requirements**

+24 V	110 mA
+12 V	60 mA
+6 V	1.45 A.
-6 V	900 mA



## 2.3 Front panel




Fig. 2.1: Mod. N914 front and rear panel

## 2.4 Input/Output connections

<b>IN#:</b>	<u>Mechanical specifications:</u> LEMO 00 type connectors <u>Electrical specifications:</u> Accepts negative charge pulses from PMT detectors. Impedance 25 Ohm.
<b>OUTL#:</b>	<u>Mechanical specifications:</u> LEMO 00 type connectors <u>Electrical specifications:</u> Integrator output with -0.83 mV/pC sensitivity. Impedance: 50 Ohm, dynamics -2.1 V.
<b>OUTH#:</b>	<u>Mechanical specifications:</u> LEMO 00 type connectors <u>Electrical specifications:</u> Integrator output with -8.3 mV/pC sensitivity. Impedance: 50 Ohm, dynamics -1.7 V.
<b>ASUM:</b>	<u>Mechanical specifications:</u> LEMO 00 type connectors <u>Electrical specifications:</u> Provides the analog sum of the CH# inputs; Impedance: 50 Ohm, dynamics -1.7 V.
<b>OLSUM:</b>	<u>Mechanical specifications:</u> LEMO 00 type connectors <u>Electrical specifications:</u> Provides the analog sum of the OL# outputs; Impedance: 50 Ohm, dynamics -1.7 V.
<b>OHSUM:</b>	<u>Mechanical specifications:</u> LEMO 00 type connectors <u>Electrical specifications:</u> Provides the analog sum of the OL# outputs; Impedance: 50 Ohm, dynamics -1.7 V.
<b>MALU:</b>	<u>Mechanical specifications:</u> LEMO 00 type connectors <u>Electrical specifications:</u> PMT signal on each channel input is passed to a 10x gain stage and then a discriminator (with settable threshold), whose output is passed to a one-shot device, which converts it into a digital signal of tunable length. The resulting pulses are summed and then delivered by the MALU output, which therefore provides a step function (Majority) with a height proportional to the number of channels simultaneously over threshold; step amplitude: 40 mV / Ch. Impedance: 50 Ohm, dynamics +4 V.
<b>T_OUT#:</b>	<u>Mechanical specifications:</u> 8+8-pin, 3M 3408-5202 Header-type connector <u>Electrical specifications:</u> PMT signal on each channel input is passed to a 10x gain stage and then a discriminator (settable threshold), whose output is passed to a one-shot device, which converts it into a digital signal of tunable length. This signal will be sent to a TDC for arrival times detection. LVDS standard

---

## 2.5 Front panel trimmers

- LOC THR:** Each input channel is provided with a screwdriver trimmer for setting the LOCAL discriminator threshold for the production of MALU and T\_OUT# signals; such threshold can be set in a  $-1.5 \text{ mV} \div -60 \text{ mV}$  range.
- COM THR:** On the lower section of the front panel, a screwdriver trimmer and test point allows setting the COMMON discriminator threshold for the production of MALU and T\_OUT# signals; actual threshold is 1% of test point value. Internal jumpers (see § 2.6) allows to select between Local and Common threshold. The relevant test point provides a  $-150 \text{ mV} \div -3000 \text{ mV}$  output, corresponding to a  $-1.5 \text{ mV} \div -30 \text{ mV}$  set threshold.
- TDC WDT:** This trimmer allows to set the T\_OUT# width ( $250\text{ns} \div 5\mu\text{s}$  range)
- MAJ WDT:** This trimmer allows to set the duration of the majority steps ( $250\text{ns} \div 5\mu\text{s}$  range).

---

## 2.6 Internal jumpers

Eight internal jumpers (J8, J16, J24 and J32: they are four on each printed board) allow to select whether to use the Local (INT) or Common (EXT) threshold setting (via front panel trimmers). Each jumper handles one channel (Printed board left: J8: CH0, J16: CH2, J24: CH4, J32: CH6; Printed board right: J8: CH1, J16: CH3, J24: CH5, J32: CH7). Refer to Fig 2.2 for the jumpers' location.



**Fig. 2.2: Internal jumpers location<sup>1</sup>**

<sup>1</sup> N914 has two identical printed board inside

## 2.7 Technical specification table

**Table 2.2: Mod. N914 Technical Features**

<b>Packaging</b>	One unit wide NIM unit
<b>Input sensitivity</b>	OUTL#: 0.83 mV/pC OUTH#: 8.3 mV/pC
<b>T# width</b>	250ns÷5µs
<b>MALU width</b>	250ns÷5µs
<b>MALU step amplitude</b>	40 mV/channel
<b>Integrator time constant</b>	1.2 nF // 100 kOhm
<b>ASUM out bandwidth</b>	150 MHz
<b>Discriminator Sensitivity</b>	1.5 mV

**CAEN S.p.A.**

Via Vetràia 11  
55049 - Viareggio  
Italy  
Phone +39 0584 388 398  
Fax +39 0584 388 959  
info@caen.it  
[www.caen.it](http://www.caen.it)

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt  
Germany  
Tel. +49 (0)212 254 4077  
Mobile +49 (0)151 16 548 484  
info@caen-de.com  
[www.caen-de.com](http://www.caen-de.com)

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA  
Phone: +1 (718) 981-0401  
Fax: +1 (718) 556-9185  
info@caentechnologies.com  
[www.caentechnologies.com](http://www.caentechnologies.com)

**CAENspa INDIA Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri (W)  
Mumbai, Mumbai City,  
Maharashtra, India, 400053  
info@caen-india.in  
[www.caen-india.in](http://www.caen-india.in)



Copyright © CAEN SpA. All rights reserved. Information in this publication supersedes all earlier versions. Specifications subject to change without notice.