



**LOGIC GATES**  
AND/OR/Majority  
NIM/TTL translator  
Veto  
Coincidence



**COUNTERS and TIMERS**  
Scaler  
Rate Meter  
Counter  
Start/Stop



**TIMING**  
Events list  
ToF  
ToT




**DIGITAL PULSE GENERATOR**  
NIM/TTL generator  
Poisson statistics  
Pattern generator

### KEY FEATURES

- ◆ Double-width NIM unit
- ◆ Wide range of **user-selectable functionalities**: logic functions, counting, timing, digital pulse generator.
- ◆ 4 independent sections, each with:
  - 6 NIM/TTL input
  - 4 NIM/TTL output
- ◆ Input stage with **Gate&Delay** and **Voltage Threshold** for analog input signals
- ◆ **2.8" touch screen display** with user-friendly widgets to configure each section and monitor real-time data
- ◆ Ethernet (1 Gbps) and USB2.0 connectivity
- ◆ **Web-based Graphical User Interface** to remotely control the module and access advanced features



The N1081A is a **laboratory tool** that incorporates in a single NIM module the most common functionalities that you need to implement the **logic capabilities** of your experiment.

The module is organized in four sections, with 6 inputs and 4 outputs each (selectable impedance) accepting TTL/NIM signals. Each section is configurable independently according to one of the available **pre-programmed functions**, like scaler, counter, time stamping, digital pulse generator, etc. The section programming can be done using the **2.8" touch screen display** or via the **web-interface**. Each section integrates a programmable voltage threshold for analog input signals and an asynchronous Gate&Delay. This allows the user to trim at best the needed parameters and to perform accurate measurements.



On the touch screen interface, each function is associated to a **widget**, for quick configuration and monitor purposes, and an online help with cabling instructions is directly accessible on the display.

The **web-interface** allows the user to remotely configure the instrument, monitor the functions output, dump data on file or history chart and access to the most advanced features, like the logic analyser and the configuration of advanced functions. No software installation is required!



## TECHNICAL SPECIFICATIONS

<b>Form factor</b>	2U NIM
<b>Power consumption</b>	1.5 A @ 12V (Typ.)
<b>Logic Functions</b>	Logic Gates (AND/OR/Majority, ...) Fan IN-OUT NIM/TTL translator Gate&Delay Veto Coincidence Scaler, MCS Rate Meter Counter Timer Timing (timestamping, ToF, ToT) Digital Pulse Generator
<b>I/Os</b>	4 independent sections: ◇ 6 NIM/TTL input: Gate&Delay 5ns resolution, Programmable Voltage Threshold for analog signals ◇ 4 NIM/TTL output
<b>I/O delay</b>	20 ns
<b>Connectors</b>	LEMO
<b>Counter max. frequency</b>	30 MHz TTL 40 MHz NIM
<b>Pulse Generator</b>	NIM/TTL Max. frequency 10 MHz / 20 ns pulse width
<b>Integration gate</b>	1 $\mu$ s ... 3600 s
<b>Connectivity</b>	Ethernet ( Gbps), USB 2.0
<b>Touch Screen Display</b>	2.8" LCD, Transmissive type
<b>Software</b>	◇ Touch screen widgets ◇ Web-based Graphical User Interface




### Ordering Option

Ordering code	Description
WN1081AXAAAA	N1081A—Four-fold programmable Logic Unit

## Register your device

Register your device to your **MyCAEN+** account and get access to our customer services, such as notification for new firmware or software upgrade, tracking service procedures or open a ticket for assistance. **MyCAEN+** accounts have a dedicated support service for their registered products. A set of basic information can be shared with the operator, speeding up the troubleshooting process and improving the efficiency of the support interactions.

**MyCAEN+** dashboard is designed to offer you a direct access to all our after sales services. Registration is totally free, to create an account go to <https://www.caen.it/become-mycaenplus-user> and fill the registration form with your data.





**CAEN S.p.A.**

Via Vetraria 11  
55049 - Viareggio  
Italy  
Phone +39 0584 388 398  
Fax +39 0584 388 959  
info@caen.it  
**www.caen.it**



**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt  
Germany  
Tel. +49 (0)212 254 4077  
Mobile +49 (0)151 16 548 484  
info@caen-de.com  
**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA  
Phone: +1 (718) 981-0401  
Fax: +1 (718) 556-9185  
info@caentechnologies.com  
**www.caentechnologies.com**

**CAENspa INDIA Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri (W)  
Mumbai, Mumbai City,  
Maharashtra, India, 400053  
info@caen-india.in  
**www.caen-india.in**

